Cortisol

Enzyme immunoassay for the quantitative determination of Cortisol in human serum or plasma

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1. INTRODUCTION

Cortisol is a steroid hormone released from the adrenal cortex in response to a hormone called ACTH (produced by the pituitary gland), it is involved in the response to stress; it increases blood pressure, blood sugar levels, may cause infertility in women, and suppresses the immune system.

Cortisol acts through specific intracellular receptors and has effects in numerous physiologic systems, including immune function, glucose-counter regulation, vascular tone, substrate utilization and bone metabolism. Cortisol is excreted primarily in urine in an unbound (free) form.

Cortisol is bound with high affinity in plasma from corticosteroid-binding globulin (CBG, transcortin) and from albumin. Only free cortisol is available to most receptors.

The amount of cortisol present in the serum undergoes diurnal variation, with the highest levels present in the early morning, and lower levels in the evening, several hours after the onset of sleep. Highest levels are at about 6 – 8 a.m. and lowest levels are at about midnight. These normal endogenous functions are basis for the physiological consequences of chronic stress- prolonged cortisol secretion causes muscle wastage, hyperglycaemia, and suppresses immune/ inflammatory responses. The same consequences arise from long-term use of glucocorticoid drugs.

2. INTENDED USE

Competitive immunoenzymatic colorimetric method for quantitative determination of Cortisol in human serum or plasma.

3. PRINCIPLE OF THE ASSAY

Microtiter strip wells are precoated with anti-Cortisol antibodies (solid-phase). Cortisol in the sample competes with added horseradish peroxidase labelled Cortisol (enzyme-labelled antigen) for antibody binding. After incubation a bound/free separation is performed by solid-phase washing. The immune complex formed by enzyme-labelled antigen is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product. The intensity of this product is **inversely** proportional to the amount of Cortisol in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorption at 450 nm is read using an ELISA microwell plate reader.

4. MATERIALS

4.1. Reagents supplied

- **Anti-Cortisol IgG Coated Wells:** 12 breakapart 8-well snap-off strips coated with anti-Cortisol IgG; in resealable aluminium foil.
- **Stop Solution:** 1 bottle containing 15 ml sulphuric acid, 0.15 mol/l (avoid any skin contact).
- **Cortisol-HRP Conjugate.:** 1 bottle containing 21 ml of horseradish peroxidase labelled Cortisol.
- **TMB Substrate Solution:** 1 bottle containing 15 ml 3, 3', 5, 5'-tetramethylbenzidine (H₂O₂-TMB 0.26 g/l) (avoid any skin contact).
- **Wash solution 10 x conc.:** 1 bottle containing 50 ml of a 10 x concentrated solution of phosphate buffer 0.2 M.
- **Cortisol control:** 1 bottle containing 1 ml of a lot-specific control solution. The concentration is indicated on the label of the bottle.
- **Cortisol Standards:** 5 bottles, 1 ml each
 - Standard 0: 0 ng/ml
 - Standard 1: 10 ng/ml
 - Standard 2: 50 ng/ml
 - Standard 3: 150 ng/ml
 - Standard 4: 500 ng/ml

4.2. Materials supplied

- 1 Cover foil
- 1 instructions for use (IFU)

4.3. Materials and Equipment needed

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450 nm, 620-630 nm
- Manual or automatic equipment for rinsing wells
- Pipettes to deliver volumes between 10 and 1000 µl
- Vortex tube mixer
- Distilled water
- Disposable tubes
- Timer

5. STABILITY AND STORAGE

The reagents are stable up to the expiry date stated on the label when stored at 2...8 °C in the dark.

6. REAGENT PREPARATION

It is very important to bring all reagents, samples and standards to room temperature (22...28°C) for at least 30 minutes before starting the test run!

At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.

6.1. Coated snap-off Strips

The ready to use break apart snap-off strips are coated with anti-Cortisol IgG antibodies. Store at 2...8 °C. Open the bag only when it is at room temperature. *Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C; stability until expiry date.*

6.2. Cortisol-HRP Conjugate

The conjugate is ready to use. Mix gently for 5 minutes with rotating mixer. After first opening stable for another 6 months at 2...8°C.

6.3. Cortisol Standards

The standards are ready to use and have the following concentration of Cortisol:

Standard 0:	0 ng/ml
Standard 1:	10 ng/ml
Standard 2:	50 ng/ml
Standard 3:	150 ng/ml
Standard 4:	500 ng/ml

After first opening stable for another 6 months at 2...8°C.

6.4. TMB Substrate Solution

The bottle contains 15 ml of a tetramethylbenzidine/hydrogen peroxide system. The reagent is ready to use and has to be stored at 2...8°C in the dark. *The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.*

6.5. Stop Solution

The bottle contains 15 ml 0.15 mol/l sulphuric acid solution. This ready to use solution has to be stored at 2...8°C.

6.6. Wash Solution

Dilute the concentrated solution with distilled water to a final volume of 500 ml prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted solution is stable for 30 days at 2...8°C. In the concentrated solution it is possible to observe the presence of crystals, in this case mix at room temperature until complete dissolution of crystals. For greater accuracy dilute the whole bottle of concentrated wash solution to 500 ml and take care that all crystals are transferred by washing the bottle, then mix until crystals are completely dissolved.

6.7. Cortisol Control

The bottle contains 1 ml of a lot-specific control solution. The concentration is indicated on the label.

7. SPECIMEN COLLECTION AND PREPARATION

The determination of Cortisol can be performed in plasma as well as in serum.

Store the sample at -20°C if the determination is not performed on the same day as the sample collection.

Avoid repeated freezing and thawing.

Treatment of the patient with corticosteroids, natural or synthetic steroids can impair Cortisol determination.

8. ASSAY PROCEDURE

8.1. Test Preparation

Please read the instructions for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instructions for use as described. Prior to commencing the assay, the distribution and identification plan for all specimens and controls should be carefully established. Select the required number of microtiter strips or wells and insert them into the holder. Please allocate at least:

1 well	(e.g. A1)	for the substrate blank
2 wells	(e.g. B1+C1)	for standard 0
2 wells	(e.g. D1+E1)	for standard 1
2 wells	(e.g. F1+G1)	for standard 2
2 wells	(e.g. H1+A2)	for standard 3
2 wells	(e.g. B2+C2)	for standard 4
2 wells	(e.g. D2+E2)	for control

It is recommended to determine standards, control and patient samples in duplicate.

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each standard and each patient sample.

1. Dispense 20 µl standards, control and samples into their respective wells.
2. Add 200 µl Cortisol-HRP Conjugate to each well. Leave well A1 for substrate blank.
3. Cover wells with the foil supplied in the kit.
4. **Incubate for 1 hour at 37 °C.**

When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µl of diluted wash solution. Avoid overflows from the reaction wells. During each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

In case you use an automatic washer, it is advised to do 6 washing steps.

Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.

5. Dispense 100 µl TMB Substrate Solution into all wells.
6. **Incubate for exactly 15 min at room temperature (22...28°C) in the dark.**
7. Dispense 100 µl Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution. Shake the microplate gently. *Any blue colour developed during the incubation turns into yellow.*
8. Measure the absorbance (E) of the specimen at 450 nm against a reference wavelength of 620-630 nm or against blank within 5 minutes.

8.2. Measurement

Adjust the ELISA Microwell Plate Reader **to zero** using the **substrate blank in well A1**.

If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard and patient sample.

Where applicable calculate the **mean absorbance values** of all duplicates.

9. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of Cortisol for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the standard curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

10. RESULTS

10.1. Calculation of results

Calculate the mean absorbance for each point of the standard curve and each sample. Plot the mean value of absorbance of the standards against concentration. Draw the best-fit curve through the plotted points. (e. g.: Four Parameter Logistic). Interpolate the values of the samples on the standard curve to obtain the corresponding values of the concentrations expressed in ng/ml.

10.2. Reference values

Serum or plasma cortisol reference values are:

60 – 230 ng/ml between 8:00 -10:00 a.m.

30 – 150 ng/ml at 4:00 p.m.

Patient treated with ACTH 280 – 600 ng/ml

Patient treated with dexamethasone 0 – 50 ng/ml

11. SPECIFIC PERFORMANCE CHARACTERISTICS

11.1. Precision

Intra Assay Variation

Within run variation was determined by replicate determination (16x) of three different control sera in one assay. The within assay variability is $\leq 5,1\%$.

Inter Assay Variation

Between run variation was determined by replicate measurements (10x) of three different control sera in different lots. The between assay variability is $\leq 11.0\%$.

11.2. Cross Reactivity

The cross reaction of the antibody calculated at 50% according to Abraham is:

Cortisol	100 %
Prednisolone	46.2 %
11-Deoxycortisol	4 %
Cortisone	3.69 %
Prednisone	3.10 %
Corticosterone	1.34 %
11 α OH Progesterone	1 %
Progesterone	<0.01%
Aldosterone	<0.01%
Pregnenolone	<0.01%
17 beta Estradiol	<0.01%
Estrone 3-sulfate	<0.01%
Estriol	<0.01%
Testosterone	<0.01%
Spirolactone	<0.01%
DHEA	<0.01%
DHEA-S	<0.01%
Androstenedione	<0.01%
Androsterone	<0.01%
DHT	<0.01%
Danazol	<0.01%
Cholesterol	<0.01%
Dexamethasone	<0.01%

11.3. Specificity: interfering substances

Interference by Bilirubin, Hemoglobin and Triglycerides has been investigated in the Cortisol ELISA kit:

Substance	Assayed Conc.	Interference
Bilirubin	0.2 mg/ml	No
Hemoglobin	2 mg/ml	No
Triglycerides	6 mg/ml	No

No interference has been observed with the substances under investigation; anyway, according to good laboratory practices, avoid using highly lipemic or haemolysed samples.

11.4. Specificity: plasma and SST tube

Interference in plasma and SST (serum separation tube) samples has been evaluated. Serum obtained from the same patient has been used as reference.

Sample	Interference
SST (serum separation tube)	No
EDTA plasma	No
Lithium heparin plasma	No
Sodium-heparin plasma	No

No interference has been observed.

11.5. Analytic Sensitivity

The lowest detectable concentration of Cortisol that can be distinguished from the 0 standard is 2.42 ng/ml at the 95 % confidence limit.

11.6. Accuracy

The recovery of 15 - 35 - 70 ng/ml of Cortisol added to samples gave an average value (\pm SD) of 101.61% \pm 4.42% with reference to the original concentrations.

The dilution test performed on three sera diluted 2 - 4 and 8 times gave an average value (\pm SD) of 106.18% \pm 4.68%.

11.7. Correlation

The Cortisol ELISA kit has been compared to a commercially available chemiluminescence Cortisol kit. 61 serum samples have been analysed.

The linear regression curve has been calculated:

$$Y = 1.09 \cdot X - 19.23$$

$$r^2 = 0.904$$

12. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use. Not for internal or external use in humans or animals.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive. Nevertheless, all materials should still be regarded and handled as potentially infectious.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Maximum precision is required for dispensation of the reagents.
- This method allows the determination of Cortisol from 10 to 500 ng/ml
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Avoid the exposure of TMB substrate to direct sunlight, metal or oxidants. Do not freeze the solution.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Microbiologically contaminated samples should not be used in the assay. Highly lipemic or haemolysed specimens should also not be used.

- Plate readers measure vertically. Do not touch the bottom of the wells.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- Some reagents contain small amounts of Proclin 300® as preservative. Avoid contact with skin or mucosa.
- The ELISA is only designed for qualified personnel who are familiar with good laboratory practice.

12.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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


Cortisol

Test Preparation

Prepare reagents and samples as described.
Establish the distribution and identification plan for all specimens and controls.
Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate blank	Standard 0 - 4	Control	Sample
Standard 0 - 4	-	20 µl	-	-
Control	-	-	20 µl	-
Sample	-	-	-	20 µl
Conjugate	-	200 µl	200 µl	200 µl
Cover wells with foil supplied in the kit Incubate for 1 h at 37°C Wash each well three times with 300 µl diluted wash solution In case you use an automatic washer, it is advised to do 6 washing steps.				
TMB Substrate	100 µl	100 µl	100 µl	100 µl
Incubate for exactly 15 min at room temperature (22...28°C) in the dark				
Stop Solution	100 µl	100 µl	100 µl	100 µl
Shake the microplate gently Photometric measurement at 450 nm, 620-630 nm				

Gold Standard Diagnostics Frankfurt GmbH

Waldstrasse 23 A6

63128 Dietzenbach, Germany

Tel.: +49 6074 23698-0

Fax: +49 6074 23698-900

E-Mail: info.frankfurt@eu.goldstandarddiagnostics.com

Website: clinical.goldstandarddiagnostics.com

DNOV001_CORT_IFU_rev01_fromLot_6129AN