Androstenedione

CE

Enzyme immunoassay for the quantitative determination of Androstenedione in human serum or plasma

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ENGLISH

1. INTRODUCTION

Androstenedione (also known as 4-androstenedione) is a 19-carbon steroid hormone produced in the adrenal glands and the gonads as an intermediate step in the biochemical pathway that produces the androgen testosterone and the estrogens estrone and estradiol. It is the common precursor of male and female sex hormones. Some androstenedione is also secreted into the plasma, and may be converted in peripheral tissues to testosterone and estrogens.

Androstenedione has relatively weak androgenic activity, estimated at ~20% of testosterone. However, serum androstenedione levels often exceed testosterone in both normal and disease states. Secretion and production rates also exceed those of testosterone in women in whom significant extra-adrenal conversion of androstenedione to testosterone occurs.

In premenopausal women the adrenal glands and ovaries each produces about half of the total androstenedione (about 3 mg/day). After menopause androstenedione production is about halved, primarily due to the reduction of steroid secreted by the ovary. Nevertheless, androstenedione is the principal steroid produced by the postmenopausal ovary.

Measurement of serum androstenedione provides a useful marker of androgen biosynthesis. Elevated androstenedione levels have been demonstrated in virilizing congenital adrenal hyperplasia. Serum androstenedione levels are also increased in polycystic ovary syndrome, and in case of hirsutism in women. Elevated serum androstenedione levels may also occur in adrenal and ovarian virilizing tumors.

2. INTENDED USE

Competitive immunoenzymatic colorimetric method for quantitative determination of Androstenedione in human serum or plasma.

3. PRINCIPLE OF THE ASSAY

Microtiterplates are coated with anti-Androstenedione antibodies (solid-phase). Androstenedione in the sample competes with added horseradish peroxidase labelled Androstenedione (enzyme-labelled antigen) for antibody binding. After incubation a bound/free separation is performed by solid-phase washing. The immune complex formed by enzyme-labelled antigen is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product. The intensity of this product is **inversely** proportional to the amount of Androstenedione in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorption at 450 nm is read using an ELISA Microtiterplate reader. Androstenedione concentration in the sample is calculated through a calibration curve.

4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 breakapart 8-well snap-off strips coated with anti-Androstenedione in resealable aluminium foil.
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, < 5 % (w/v); (avoid any skin contact).
- **Conjugate:** 1 bottle containing 21 mL of horseradish peroxidase labelled Androstenedione.
- **TMB-Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (H_2O_2 -TMB 0.26 g/L) (avoid any skin contact).
- **Wash Solution 10x conc.:** 1 bottle containing 50 mL of a 10x concentrated solution (NaCl 160 g/L, Tween-20 10 g/L, 0.2 M phosphate buffer, pH 7.4)
- **Control:** 1 bottle containing 1 mL of a control solution, ready for use. The concentration is indicated on the label of the vial.
- **Standards:** 6 bottles, 1 mL each
 - Standard 0: 0 ng/mL
 - Standard 1: 0.1 ng/mL
 - Standard 2: 0.4 ng/mL
 - Standard 3: 1.2 ng/mL
 - Standard 4: 4.0 ng/mL
 - Standard 5: 10.0 ng/mL

4.2. Materials supplied

- 1 Cover foil
- 1 Instructions for use (IFU)

4.3. Materials and Equipment needed

- 37 °C Incubator
- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450 nm, 620 – 630 nm
- Manual or automatic equipment for rinsing Microtiterplates
- Pipettes
- rotation mixer
- Distilled water
- Disposable tubes
- Timer

5. STABILITY AND STORAGE

The reagents are stable up to the expiry date stated on the label when stored at 2...8 °C in the dark.

6. REAGENT PREPARATION

It is very important to bring all reagents, samples, control and standards to room temperature (22...28 °C) for at least 30 minutes before starting the test run! At the end of the assay, store immediately the reagents at 2...8°C; avoid long exposure to room temperature.

6.1. Microtiterplate

The ready to use break apart snap-off strips are coated with anti-Androstenedione antibodies. Open the bag only when it is at room temperature. *Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8°C; stability until expiry date.*

6.2. Conjugate

Androstenedione-HRP Conjugate is a ready to use solution.

6.3. Standards

The standards are ready to use and have the following concentration:

Standard 0:	0 ng/mL
Standard 1:	0.1 ng/mL
Standard 2:	0.4 ng/mL
Standard 3:	1.2 ng/mL
Standard 4:	4.0 ng/mL
Standard 5:	10.0 ng/mL

Before using, leave the standards on a rotating mixer for at least 5 min. After first opening stable for another 6 months at 2...8 °C.

6.4. Control

Before use mix 5 min. with a rotating mixer. The bottle contains 1 mL of a lot-specific control solution, ready to use. The concentration is indicated on the label.

6.5. TMB Substrate Solution

The bottle contains 15 mL of a tetramethylbenzidine/hydrogen peroxide system. The reagent is ready to use and has to be stored at 2...8 °C in the dark. *The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.*

6.6. Stop Solution

The bottle contains 15 mL sulphuric acid solution (H314). This ready to use solution has to be stored at 2...8 °C.

6.7. Wash Solution

Dilute the concentrated wash solution concentrate with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2...8°C. It is possible to observe the presence of crystals in the concentrated wash solution; in this case mix at room temperature until the crystals are completely dissolved; for greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care to transfer completely the crystals, then mix until crystals are completely dissolved.

7. SPECIMEN COLLECTION AND PREPARATION

The determination of Androstenedione can be performed in human plasma as well as in serum. Store the sample at -20 °C if the determination is not performed on the same day as sample collection. If samples are stored frozen, mix thawed samples well before testing. *Avoid repeated freezing and thawing.*
Dilute samples > 10 ng/mL 1:2 with standard 0.

7.1. Precaution

- Some reagents contain small amounts of Proclin 300® as preservative. Avoid the contact with skin or mucosa.
- The TMB substrate contains an irritant, which may be harmful if inhaled, ingested, or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Do not use heavily haemolysed or highly lipemic samples.
- Maximum precision is required for dispensation of the reagents.
- This method allows the determination of Androstenedione from 0.1 – 10.0 ng/mL.
- Avoid the exposure of TMB substrate to direct sunlight, metal or oxidants. Do not freeze the solution.
- Treatment of the patient with cortisone, natural or synthetic steroids can impair Androstenedione determination.

8. ASSAY PROCEDURE

8.1. Test Preparation

Please read the test instructions for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instructions for use as described. Prior to commencing the assay, the distribution and identification plan for all specimens and standards should be carefully established. Select the required number of microtiter strips or wells and insert them into the holder. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than one plate is used, it is recommended to repeat the dose response curve. Please allocate at least:

1 well	(e.g. A1)	for the substrate blank
2 wells	(e.g. B1+C1)	for standard 0
2 wells	(e.g. D1+E1)	for standard 1
2 wells	(e.g. F1+G1)	for standard 2
2 wells	(e.g. H1+A2)	for standard 3
2 wells	(e.g. B2+C2)	for standard 4
2 wells	(e.g. D2+E2)	for standard 5
2 wells	(e.g. F2+G2)	for control

It is recommended to determine standards, control and patient samples in duplicate.

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each standard and each patient sample.

1. Dispense 25 µL standards, control and samples into their respective wells.
2. Add 200 µL Conjugate to each well. Leave well A1 for substrate blank.
3. Cover wells with the foil supplied in the kit.
4. **Incubate for 1 hour at 37 °C.**
5. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL diluted wash solution. Avoid overflows from the reaction wells.

Important note: During each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: In case you use automatic equipment, wash the wells at least 5 times.

Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.

6. Dispense 100 µL TMB Substrate Solution into all wells.
7. **Incubate for exactly 15 min at room temperature (22...28°C) in the dark.**
8. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution. Shake the Microtiterplate gently.
Any blue colour developed during the incubation turns into yellow.
9. Measure the absorbance of the specimen at 450 nm against a reference wavelength of 620-630 nm or against blank within 5 minutes.

8.2. Measurement

Adjust the ELISA Microtiterplate Reader **to zero** using the **substrate blank in well A1**.

If - due to technical reasons - the ELISA Microtiterplate reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard and patient sample.

Bichromatic measurement using a reference wavelength of 620-630 nm is recommended.

Where applicable calculate the **mean absorbance values** of all duplicates.

9. RESULTS

9.1. Calculation of Results

Calculate the mean absorbance for each point of the standard curve and each sample. Plot the mean value of absorbance of the standards against concentration. Draw the best-fit curve through the plotted points (e.g., Four Parameter Logistic). Interpolate the values of the samples on the standard curve to obtain the corresponding values of the concentrations expressed in ng/mL.

9.2. Reference Values

The following values for serum or plasma Androstenedione should be considered as a guideline:

Men		0.60 – 2.7 ng/mL
Women	follicular phase	0.75 – 3.1 ng/mL
	luteinic phase	0.94 – 3.2 ng/mL

Please pay attention to the fact that the determination of a range of expected values for a “normal” population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore, each laboratory should consider the range given by the manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

10. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of Androstenedione for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the standard curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

11. SPECIFIC PERFORMANCE CHARACTERISTICS

11.1. Precision

Intra Assay Variation

Within run variation was determined by replicate determination (20x) of three different control sera in one assay. The within assay variability is $\leq 10.0\%$.

Inter Assay Variation

Between run variation was determined by replicate measurements (10x) of three different control sera in different lots. The between assay variability is $\leq 9.5\%$.

11.2. Specificity

The cross reaction of the antibody calculated at 50% according to Abraham is:

Androstenedione	100 %
5 α -dihydrotestosterone	0.05 %
DHEA	0.05 %
Epitestosterone	0.04 %
DHEA-S	0.027 %
Cortisol	0.008 %
Progesterone	0.007 %
Estrone	0.007 %
Testosterone	< 0.001 %
17B-Estradiol	< 0.001 %
Estriol	< 0.001 %
Aldosterone	< 0.001 %

11.3. Sensitivity

The lowest detectable concentration of Androstenedione that can be distinguished from the zero standard is 0.01 ng/mL at the 95 % confidence limit.

11.4. Accuracy

The recovery of 0.4 – 0.8 – 1.6 – 3.2 ng/mL Androstenedione added to a sample gave an average value ($\pm$ SD) of 100.91 % $\pm$ 5.61 % with reference to the original concentrations. The dilution test performed on three sera diluted 2 – 4 – 8 – 16 times gave an average value ($\pm$ SD) of 107.18 % $\pm$ 3.03.

11.5. Correlation with RIA

The Androstenedione ELISA was compared to a commercially available chemiluminescent method. 60 serum samples were analysed accordingly in both test systems.

The linear regression curve was calculated:

$$Y = 0,92 \cdot X - 0,02$$

$$r = 0,84$$

12. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the test kits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use by professional persons. Not for internal or external use in humans or animals.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive. Nevertheless, all materials should still be regarded and handled as potentially infectious.
- Do not interchange reagents or Microtiterplates of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- The incomplete or inaccurate removal of liquid from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- For reproducible results it is important that the time of reaction in each well is held constant. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than ten minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate.
- Addition of the TMB substrate solution initiates a kinetic reaction, which is terminated by the addition of the stop solution. Therefore, the TMB substrate and the stop solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Follow Good Laboratory Practice (GLP) for handling blood reagents.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- The ELISA is only designed for qualified personnel following the standards of good laboratory practice (GLP).

12.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

For information about the packaging materials refer to PACKAGING MATERIALS.

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Judd H. and Yen S. (1973) J.Clin.Endoc.& Metab 36, 475


Abraham G. (1974) J. Clin.Endoc 39, 340.

Hillier S.G. and de Zwart F. (1985) Intraovarian Steroids Metabolism and Puberty: roles of Androgens in Adolescence Female, 79th Year book Medical Publishers Inc: Chicago.

Venturoli S. et al (1987) Fertility and Sterility 48 (1), 78.

Venturoli S. et al (1986) Hormone Res .24,269.

PACKAGING MATERIALS/VERPACKUNGSMATERIALIEN/MATÉRIELS D'EMBALLAGE/MATERIALI D'IMBALLAGGIO/MATERIALES DE EMBALAJE

 PAP 21	 PAP 21	 PAP 22	MTP  ALU / LDPE 90
--	--	--	--

SYMBOLS KEY / SYMBOLSCHLÜSSEL / EXPLICATION DES SYMBOLES / LEGENDA / SÍMBOLOS

	Manufactured by / Hergestellt von / Fabriqué par / Prodotto da / Fabricado por
IVD	In Vitro Diagnostic Medical Device / In Vitro Diagnosticum / Dispositif médical de diagnostic in vitro / Diagnostico in vitro / Producto para diagnóstico In vitro
LOT	Lot Number / Chargenbezeichnung / Numéro de lot / Lotto / Número de lote
	Expiration Date / Verfallsdatum / Date de péremption / Scadenza / Fecha de caducidad
	Storage Temperature / Lagertemperatur / Température de conservation / Temperatura di conservazione / Temperatura de almacenamiento
	Keep away from sunlight / Vor direkter Sonneneinstrahlung schützen / Tenir à l'écart de la lumière du soleil / Tenere lontano dalla luce solare / Mantener alejado de la luz solar
CE	CE Mark / CE-Zeichen / Marquage CE / Marchio CE / Marca CE
UDI	Unique Device Identifier / Eindeutige Produktidentifizierung / identification unique des dispositifs / identificazione unica del dispositivo / identificación única del producto
REF	Catalogue Number / Katalog Nummer / Référence du catalogue / Numero di codice / Número de Catálogo
	Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter le mode d'emploi / Consultare le istruzioni per l'uso / Consultar las Instrucciones de Uso
MTP	Microtiterplate / Mikrotiterplatte / Plaque de Microtitrage / Piastra di Microtitolazione / Placa de Microtitulación
CONJ	Conjugate / Konjugat / Conjugué / Coniugato / Conjugado
CONTROL	Control / Kontrolle / Contrôle / Controllo / Control
CAL	Calibrator resp. Standard / Kalibrator bzw. Standard / Calibrateur resp Etalon / Calibratore ossia Standard / Calibrador o bien Estándar
SOLN STOP	Stop Solution / Stopplösung / Solution d'arrêt / Soluzione bloccante / Solución de Parada
SUB TMB	TMB Substrate Solution/ TMB-Substratlösung/ Substrat TMB/ soluzione substrato TMB/ solución substrato TMB
WASH BUF 10x	Wash Solution 10x conc. / Waschlösung 10x konz. / Solution de lavage concentré x 10 / soluzione di lavaggio concentrazione x10 / solución de lavado concentrado x10
	Contains sufficient for "n" tests / Ausreichend für "n" Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para "n" tests

SUMMARY OF TEST PROCEDURE / KURZANLEITUNG TESTDURCHFÜHRUNG /
RÉSUMÉ DE LA PROCEDURE DE TEST / SCHEMA DELLA PROCEDURA /
RESUMEN DE LA TÉCNICA

SCHEME OF THE ASSAY

Androstenedione

Test Preparation

Prepare reagents and samples as described.
Establish the distribution and identification plan for all specimens and.
Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate blank	Standard 0 - 5	Control	Sample
Standard 0 – 5	-	25 µL	-	-
Control	-	-	25 µL	-
Sample	-	-	-	25 µL
Conjugate	-	200 µl µL	200 µL	200 µL
Cover wells with foil supplied in the kit Incubate for 1 h 37 °C Wash each well three times with 300 µl diluted wash solution In case you use automatic equipment, wash the wells at least 5 times.				
TMB Substrate	100 µL	100 µL	100 µL	100 µL
Incubate for exactly 15 min at room temperature in the dark				
Stop Solution	100 µL	100 µL	100 µL	100 µL
Shake microplate gently Photometric measurement at 450 nm, 620 – 630 nm				

Gold Standard Diagnostics Frankfurt GmbH

Waldstrasse 23 A6

63128 Dietzenbach, Germany

Tel.: +49 6074 23698-0

Fax: +49 6074 23698-900

E-Mail: info.frankfurt@eu.goldstandarddiagnostics.com

Website: clinical.goldstandarddiagnostics.com

DNOV008_ASD_IFU_rev01_fromLot_5964BN