Urinary Cortisol

CE

Enzyme immunoassay for the quantitative determination of free Cortisol in human urine

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1. INTRODUCTION

Cortisol is a steroid hormone released from the adrenal cortex in response to an hormone called ACTH (produced by the pituitary gland), it is involved in the response to stress; it increases blood pressure, blood sugar levels, may cause infertility in women, and suppresses the immune system.

Cortisol acts through specific intracellular receptors and has effects in numerous physiologic systems, including immune function, glucose-counter regulation, vascular tone, substrate utilization and bone metabolism. Cortisol is excreted primarily in urine in an unbound (free) form.

Cortisol is bound, in plasma, from corticosteroid-binding globulin (CBG, transcortin), with high affinity, and from albumin. Only free cortisol is available to most receptors.

These normal endogenous functions are the basis for the physiological consequences of chronic stress - prolonged cortisol secretion causes muscle wastage, hyperglycaemia, and suppresses immune / inflammatory responses. The same consequences arise from long-term use of glucocorticoid drugs.

The free cortisol fraction represents the metabolically active cortisol. In normal conditions, less than 1% is excreted in urine. In pathological conditions (syndrome of Cushing) the levels of free urinary cortisol are elevated, because the CBG doesn't bound the plasma cortisol in excess and it is removed with urine.

During pregnancy or estrogen treatment an increase of plasma cortisol is caused by an increment of the production of the transport protein, but the levels of free urinary cortisol remain normal to indicate correct suprarenal functionality.

This test is very useful to estimate the real suprarenal function, because it measures the free cortisol, it is the metabolically active form. Moreover the measurement of free urinary cortisol is the better parameter for the diagnosis of the Cushing's syndrome.

2. INTENDED USE

Competitive immunoenzymatic colorimetric method for quantitative determination of "free" Cortisol in urine.

3. PRINCIPLE OF THE ASSAY

Microtiter strip wells are precoated with anti-Cortisol antibodies (solid-phase). Cortisol in the sample competes with added horseradish peroxidase labelled Cortisol (enzyme-labelled antigen) for antibody binding. After incubation a bound/free separation is performed by solid-phase washing. The immune complex formed by enzyme-labelled antigen is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product. The intensity is **inversely** proportional to the amount of Cortisol in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorption at 450 nm is read using an ELISA microwell plate reader.

4. MATERIALS

4.1. Reagents supplied

- **Anti-Cortisol IgG Coated Wells:** 12 breakapart 8-well snap-off strips coated with anti-Cortisol antibodies; in resealable aluminium foil.
- **Stop Solution:** 1 bottle containing 15 ml sulphuric acid, 0.15 mol/l (avoid any skin contact), ready to use.
- **Cortisol-HRP Conjugate:** 1 bottle containing 33 ml of horseradish peroxidase labelled Cortisol, ready to use.
- **TMB Substrate Solution:** 1 bottle containing 15 ml 3, 3', 5, 5'-tetramethylbenzidine (H₂O₂-TMB 0.26 g/l) (avoid any skin contact), ready to use.
- **Wash solution 10x conc:** 1 bottle containing 50 ml (phosphate buffer 0.2M, Proclin <0.0015%)
- **Low control:** 1 vial containing 1 ml, refer to label for exact value, ready to use.
- **High control:** 1 vial containing 1 ml, refer to label for exact value, ready to use.
- **Cortisol Standards:** 5 bottles, ready to use.

Standard 0:	(4 ml)	0 ng/ml
Standard 1:	(1 ml)	1 ng/ml
Standard 2:	(1 ml)	5 ng/ml
Standard 3:	(1 ml)	30 ng/ml
Standard 4:	(1 ml)	200 ng/ml

4.2. Materials supplied

- 1 Cover foil
- 1 Instructions for use (IFU)

4.3. Materials and Equipment needed

- 37° C Incubator
- ELISA microwell plate reader, equipped for the measurement of absorbance at 450 nm, 620-630 nm
- Automatic dispenser
- Distilled water
- Disposable tubes
- Timer
- Rotating mixer

5. STABILITY AND STORAGE

The reagents are stable up to the expiry date stated on the label when stored at +2...+8 °C in the dark. At the end of the assay, store immediately the reagents at 2...8 °C; avoid long exposure to room temperature.

6. REAGENT PREPARATION

It is very important to bring all reagents, samples and standards to room temperature (22...28 °C) for at least 30 minutes before starting the test run and mix well prior to use! At the end of the assay, store immediately the reagents at 2-8 °C; avoid long exposure to room temperature.

6.1. Coated snap-off Strips

The ready to use break-apart snap-off strips are coated with anti-Cortisol antibodies. Store at +2...+8 °C. Open the bag only when it is at room temperature. *Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C; stability until expiry date.*

6.2. Cortisol-HRP Conjugate

The conjugate is ready to use. Once opened, it is stable 6 months at 2-8 °C.

6.3. Cortisol Standards

Before use, leave 5 minutes on a rotating mixer. The standards are ready to use and have the following concentration of Cortisol:

Standard 0:	0 ng/ml
Standard 1:	1 ng/ml
Standard 2:	5 ng/ml
Standard 3:	30 ng/ml
Standard 4:	200 ng/ml

After first opening stable for another 6 months at 2...8 °C.

6.4. TMB Substrate Solution

The bottle contains 15 ml of a tetramethylbenzidine/hydrogen peroxide system. The reagent is ready to use and has to be stored at 2...8 °C in the dark. *The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.*

6.5. Stop Solution

The bottle contains 15 ml 0.15 mol/l sulphuric acid solution. This ready to use solution has to be stored at 2...8 °C.

6.6. Wash Solution

Dilute the concentrated solution with distilled water to a final volume of 500 ml prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted solution is stable for 30 days at 2...8 °C. In the concentrated solution it is possible to observe the presence of crystals, in this case mix at room temperature until complete dissolution of crystals. For greater accuracy dilute the whole bottle of concentrated wash solution to 500 ml and take care that all crystals are transferred by washing the bottle, then mix until crystals are completely dissolved.

6.7. Controls

The controls are ready to use.

7. SPECIMEN COLLECTION AND PREPARATION

The determination of Cortisol with this kit should be performed in urine samples

Important note: the kit has been designed to be used on untreated urine samples; acidification treatments of the urine that lead the pH to values below 5.0 could interfere with the assay and produce aberrant results.

It is not necessary to dilute the sample. The total volume of urine excreted during 24 hours should be collected and mixed in a single container.

Urine samples which are not to be assayed immediately should be stored at 2...8 °C or at – 20 °C for longer periods (maximum 6 months).

Samples with concentration greater than 200 ng/ml have not to be diluted; such samples have to be reported as > 200 ng/ml.

7.1. Precaution

- Some reagents contain a small amount of Proclin 300® as preservative. Avoid contact with skin or mucosa.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Avoid the exposure of TMB substrate to direct sunlight, metal or oxidants. Do not freeze the solution.
- This method allows the determination of Cortisol from 0,47 ng/ml (LOD) to 200 ng/ml.
- Treatment of the patient with corticosteroids, natural or synthetic steroids can impair Cortisol determination.

8. ASSAY PROCEDURE

8.1. Test Preparation

Please read the test instructions for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instructions for use as described. Prior to commencing the assay, the distribution and identification plan for all standards, specimens and controls should be carefully established. Select the required number of microtiter strips or wells and insert them into the holder. Please allocate at least:

1 well	(e.g. A1)	for the substrate blank
2 wells	(e.g. B1+C1)	for standard 0
2 wells	(e.g. D1+E1)	for standard 1
2 wells	(e.g. F1+G1)	for standard 2
2 wells	(e.g. H1+A2)	for standard 3
2 wells	(e.g. B2+C2)	for standard 4
2 wells	(e.g. D2+E2)	for low control
2 wells	(e.g. F2+G2)	for high control

It is recommended to determine standards and patient samples in duplicate.

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each standard and each patient sample.

1. Dispense 10 µl standards, controls and samples into their respective wells.
2. Add 300 µl Cortisol-HRP Conjugate to each well. Leave well A1 for substrate blank.
3. Cover wells with the foil supplied in the kit.
4. **Incubate for 1 hour at 37 °C.**
5. When incubation has been completed, remove the foil, aspirate the content of the wells. Wash the wells 3 times with 350 µl of diluted wash solution. Avoid overflows from the reaction wells. During each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: In case you use an automatic washer, it is advised to do 6 washing steps.

Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.

6. Dispense 100 µl TMB Substrate Solution into all wells.
7. **Incubate for exactly 15 min at room temperature (22 – 28 °C) in the dark.**
8. Dispense 100 µl Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution. Shake the microplate gently.
Any blue colour developed during the incubation turns into yellow.
9. Measure the absorbance of the specimen at 450 nm against reference wavelength of 620-630 nm or against Blank within 5 minutes.

8.2. Measurement

Adjust the ELISA Microwell Plate Reader **to zero** using the **substrate blank in well A1**.

If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard and patient sample.

Where applicable calculate the **mean absorbance values** of all duplicates.

9. RESULTS

9.1. Calculation of results

Calculate the mean absorbance for each point of the standard curve and each sample. Plot the mean value of absorbance of the standards against concentration. Draw the best-fit curve through the plotted points. (Four Parameter Logistic). Interpolate the values of the samples on the standard curve to obtain the corresponding values of the concentrations expressed in ng/ml.

To calculate the cortisol concentration in urine, calculate as above and correct for total volume of volume of urine collected in 24 hours:

$$\text{ng/ml} \times \text{Vol(ml) urine 24 h} / 1.000 = \mu\text{g Cortisol/24h}$$

9.2. Reference values

To determine the normal range for urine samples, 128 apparently healthy male and female adults were tested.

Result:

Normal range urine (24h)
1,5 - 63 $\mu\text{g}/24\text{h}$

10. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of Urinary Cortisol for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the standard curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

11. SPECIFIC PERFORMANCE CHARACTERISTICS

11.1. Analytical sensitivity

Analytical sensitivity was investigated through the LOB (limit of blank), the LOD (detection limit), the LOQ (quantification limit) and the anal sensitivity (A.S.).

The following table shows the criteria of the study and the results obtained.

	Criteria	Results (ng/ml)
LOB	60 replicates of Cal 0, used as "Blank", have been investigated in 5 different sessions over 3 days	0,28
LOD	6 urine samples with low cortisol concentration have been investigate over 10 assays in duplicate, performed in 5 days	0,47
LOQ	6 urine samples with low cortisol concentration have been investigate over 10 assays in duplicate, performed in 5 days	0,56
A.S.	20 replicates of Cal 0 and 5 replicates of Cal 1 have been assayed. A.S. has been calculated by linear regression	0,22

11.2. Precision and reproducibility (complex precision)

Precision and reproducibility have been assessed through 6 different urine samples with different concentration of Cortisol. The table below shows the Within Run and Total CV%.

Sample	n°	Mean (ng/ml)	Within Run CV%	Total CV%
PS2	20	112,141	6,6%	12%
PS4	20	64,563	8,1%	12%
CT High	20	50,577	7,3%	11%
PS5	20	25,878	7,6%	10%
PS6	20	9,269	7,6%	11%
CT Low	20	3,438	7,0%	9%

11.3. Analytical specificity

11.3.1. Interfering substances

Interference for albumin, acetylsalicylic acid, ibuprofen and ascorbic acid were studied by adding the interfering substance to the urine sample with a low and high Cortisol concentration, and by comparing its concentration to the unspiked sample.

The interference has been evaluated as "significant" if it causes a concentration bias >10% between spiked and unspiked sample.

The following table shows the results obtained:

Substance	Concentration	Interference
Albumin	5 mg/dl	No
Acetylsalicylic acid	3,62 mmol/l	No
Ibuprofen	2,42 mmol/l	No
Ascorbic Acid	5 mg/l	No

Conclusion: no interference has been found for albumin, acetylsalicylic acid, ibuprofen and ascorbic acid.

11.3.2. Cross Reactivity

The cross reaction of the antibody calculated at 50% according to Abraham is:

Cortisol	100 %
Prednisolone	46.2 %
11-Deoxycortisol	4 %
Cortisone	3.69 %
Prednisone	3.10 %
11 α OH Progesterone	1 %
Progesterone	< 0.1 %
Aldosterone	< 0.1 %
Pregnenolone	< 0.1 %
17 beta Estradiol	< 0.1 %
Estrone 3-sulfate	< 0.1 %
Estriol	< 0.1 %
Testosterone	< 0.1 %
Spironolactone	< 0.1 %
DHEA	< 0.1 %
DHEA-S	< 0.1 %
Androstenedione	< 0.1 %
Androsterone	< 0.1 %
DHT	< 0.1 %
Danazol	< 0.1 %
Cholesterol	< 0.1 %
Dexamethasone	< 0.1 %

11.4. Correlation

137 urine samples were tested with the Urinary Cortisol Elisa kit and with a LC-MS method (reference)
The linear regression curve is:

$$Y = 1,008X - 0,5019$$

$$r^2 = 0,83$$

12. LIMITATIONS OF THE PROCEDURE

Sample(s), which are contaminated microbiologically, should not be used in the assay. It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than one plate is used, it is recommended to repeat the dose response curve. Addition of the substrate solution initiates a kinetic reaction, which is terminated by the addition of the stop solution. Therefore, the addition of the substrate and the stopping solution should be added in the same sequence to eliminate any time deviation during reaction. Plate readers measure vertically. Do not touch the bottom of the wells. The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.

13. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive. Nevertheless, all materials should still be regarded and handled as potentially infectious.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE nevertheless these materials should be handled as potentially infectious.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Some reagents contain a small amount of Proclin300® as preservative. Avoid contact with skin or mucosa.
- The ELISA is only designed for qualified personnel who are familiar with good laboratory practice.

WARNING:	Sulphuric acid irritates eyes and skin. Keep out of the reach of children. Upon contact with the eyes, rinse thoroughly with water and consult a doctor!
-----------------	--

13.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

For information about the packaging materials refer to PACKAGING MATERIALS.

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Urinary Cortisol

Test Preparation

Prepare reagents and samples as described.
Establish the distribution and identification plan for all specimens and controls.
Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate blank	Standard 0	Standard 1	Standard 2	Standard 3	Standard 4	Sample/ Controls
Standard 0	-	10 µl	-	-	-	-	-
Standard 1	-	-	10 µl	-	-	-	-
Standard 2	-	-	-	10 µl	-	-	-
Standard 3	-	-	-	-	10 µl	-	-
Standard 4	-	-	-	-	-	10 µl	-
Sample/ Controls	-	-	-	-	-	-	10 µl
Conjugate	-	300 µl	300 µl	300 µl	300 µl	300 µl	300 µl
<p>Cover wells with foil supplied in the kit</p> <p>Incubate for 1 h at 37 °C</p> <p>Wash each well three times with 350 µl diluted wash solution</p> <p>Important note: Gently shake the plate for 5 sec at each washing step to ensure proper cleaning. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!</p> <p>Automatic washer: In case you use automatic equipment, wash the wells at least 6 times.</p>							
TMB Substrate	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl
<p>Incubate for exactly 15 min at room temperature (22-28 °C) in the dark</p>							
Stop Solution	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl
<p>Shake the microplate gently.</p> <p>Photometric measurement at 450 nm, 620 – 630 nm</p>							

Gold Standard Diagnostics Frankfurt GmbH

Waldstrasse 23 A6

63128 Dietzenbach, Germany

Tel.: +49 6074 23698-0

Fax: +49 6074 23698-900

E-Mail: info.frankfurt@eu.goldstandarddiagnostics.com

Website: clinical.goldstandarddiagnostics.com

DNOV010_uCORT_IFU_rev01_fromLot_6130AN