# Aldosterone

**Only for in-vitro diagnostic use**

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### 1. INTRODUCTION

---

Aldosterone is a steroid hormone produced by the adrenal cortex in the adrenal gland, is the most potent mineralocorticoid in humans, it regulate sodium and potassium balance in the blood.

Aldosterone secretion appears to be stimulated primarily through the renin-angiotensin system.

Acting on mineralocorticoid receptors (MR) on principal cells in the collecting ducts of the kidneys, it increases the permeability of their apical (luminal) membrane to potassium and sodium and activates their basolateral Na<sup>+</sup>/K<sup>+</sup> pumps, stimulating ATP hydrolysis, reabsorbing sodium (Na<sup>+</sup>) ions and water into the blood, and excreting potassium (K<sup>+</sup>) ions into the urine. Aldosterone regulates plasma bicarbonate (HCO<sub>3</sub><sup>-</sup>) levels and its acid/base balance.

Aldosterone is responsible for the reabsorption of about 2% of filtered sodium in the kidneys.

Plasma aldosterone levels normally vary with body position (upright>supine) and salt intake. Overall plasma aldosterone levels show a circadian rhythm which is similar to but less marked than cortisol, with peak levels in the early morning; about 75% of the daily production is secreted between 04:00 am and 10:00 am each day. Age-related levels tend to decline from fetal through adult life.

Abnormally high plasma aldosterone concentrations can occur in adenomas, glucocorticoid-responsive hyperaldosteronism, idiopathic.

Abnormally low aldosterone secretion occurs in a number of conditions including salt-wasting forms of congenital adrenal hyperplasia, nephropathy, and renal tubular acidosis.

### 2. INTENDED USE

---

Competitive immunoenzymatic colorimetric method for quantitative determination of Aldosterone in human serum or plasma.

### 3. PRINCIPLE OF THE ASSAY

---

Microtiter strip wells are precoated with anti-Aldosterone antibodies (solid-phase). Aldosterone in the sample competes with Aldosterone-HRP conjugate for antibody binding. After incubation a bound/free separation is performed by solid-phase washing. The formed immune complex is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product. The intensity of this product is **inversely** proportional to the amount of Aldosterone in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorption at 450 nm is read using an ELISA microwell plate reader.

### 4. MATERIALS

---

#### 4.1 Reagents supplied

- **Microtiterplate:** 12 breakapart 8-well snap-off strips coated with anti- Aldosterone, in resealable aluminium foil.
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.15 mol/L (avoid any skin contact).
- **Aldosterone conjugate:** 1 bottle containing 15 mL of horseradish peroxidase labelled aldosterone.
- **TMB Substrate Solution:** 1 bottle containing 15 mL H<sub>2</sub>O<sub>2</sub>-TMB 0.26g/L (avoid any skin contact).
- **Wash solution 10x conc.:** 1 bottle containing 50 mL 0.2 M Phosphate buffer pH 7.4
- **Aldosterone control A:** 1 bottle containing 1 mL of a lot-specific control solution. The concentration is indicated on the label of the bottle. Contains < 0.01 % (v/v) CMIT/MIT (3:1).
- **Aldosterone control B:** 1 bottle containing 1 mL of a lot-specific control solution. The concentration is indicated on the label of the bottle. Contains < 0.01 % (v/v) CMIT/MIT (3:1).
- **Aldosterone Standards:** 6 bottles, 1 mL each; contain < 0.01 % (v/v) CMIT/MIT (3:1).

Standard 0:	0 pg/mL
Standard 1:	20 pg/mL
Standard 2:	80 pg/mL
Standard 3:	300 pg/mL
Standard 4:	800 pg/mL
Standard 5:	2000 pg/mL

#### 4.2 Materials supplied

- 1 Cover foil
- 1 Instructions for use (IFU)

#### 4.3 Materials and Equipment needed

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450 nm, 620-630 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing wells
- Pipettes to deliver volumes between 10 and 1000 µL
- rotating mixer
- Distilled water
- Disposable tubes
- Timer

### 5. STABILITY AND STORAGE

---

The reagents are stable up to the expiry date stated on the label when stored at 2...8 °C in the dark.

## 6. REAGENT PREPARATION

---

*It is very important to bring all reagents, samples and standards to room temperature (22...28 °C) for at least 30 minutes before starting the test run! At the end of the assay, store immediately the reagents at 2-8 °C; avoid long exposure to room temperature.*

### 6.1 Microtiterplate

The ready to use break apart snap-off strips are coated with anti-Aldosterone antibodies. Store at 2...8 °C. Open the bag only when it is at room temperature. *Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C; stability until expiry date.*

### 6.2 Aldosterone Conjugate

The conjugate is ready to use.

### 6.3 Aldosterone Standards

Put the standards on a rotating mixer for at least 5 min. before use. The calibrators are ready to use. After first use the standards are still stable for another 6 months if stored at 2...8 °C.

### 6.4 Aldosterone controls

The bottles contain 1 mL of a lot-specific control solution. The concentration is indicated on the label.

### 6.5 TMB Substrate Solution

The bottle contains 15 mL of a tetramethylbenzidine/hydrogen peroxide system. The reagent is ready to use and has to be stored at 2...8 °C in the dark. *The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.*

### 6.6 Stop Solution

The bottle contains 15 mL 0.15 mol/l sulphuric acid solution. This ready to use solution has to be stored at 2...8 °C.

### 6.7 Wash Solution

Dilute the content of the vial "10X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2...8 °C. In concentrated wash solution it is possible to observe the presence of crystals, in this case mix at room temperature until complete dissolution of crystals; for greater accuracy dilute the whole bottle of concentrated wash solution to 500 mL on taking care also transfer the crystals, then mix until crystals are completely dissolved.

## 7. SPECIMEN COLLECTION AND PREPARATION

---

The determination of Aldosterone can be performed in human plasma as well as in human serum.

- *serum samples:*

SST tubes ("Serum Separation Tube") can be used to obtain serum without any interference with the assay.

- *plasma samples:*

Plasma samples can be obtained with EDTA, Lithium heparin and Sodium heparin without any interference with the assay.

If the assay is performed on the day of sample collection, the specimen should be kept at 2...8 °C; otherwise, it should be aliquoted and stored deep-frozen (-20 °C). If samples are stored frozen, mix thawed samples well before testing. *Avoid repeated freezing and thawing.*

Samples with concentration greater than 2000 pg/mL have not to be diluted; such samples have to be reported as ">2000 pg/mL".

## 8. ASSAY PROCEDURE

---

### 8.1 Test Preparation

Please read the test instructions for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instructions for use as described. Prior to commencing the assay, the distribution and identification plan for all specimens and standards should be carefully established. Select the required number of microtiter strips or wells and insert them into the holder. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than one plate is used, it is recommended to repeat the dose response curve. Please allocate at least:

1 well	(e.g. A1)	for the substrate blank
2 wells	(e.g. B1+C1)	for standard 0
2 wells	(e.g. D1+E1)	for standard 1
2 wells	(e.g. F1+G1)	for standard 2
2 wells	(e.g. H1+A2)	for standard 3
2 wells	(e.g. B2+C2)	for standard 4
2 wells	(e.g. D2+E2)	for standard 5
2 wells	(e.g. F2+G2)	for control A
2 wells	(e.g. H2+A3)	for control B

*It is recommended to determine standards, controls and patient samples in duplicate.*

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each standard and each patient sample.

Adjust the incubator to 37 °C.

1. Dispense 50 µL standards, controls and samples into their respective wells. Add 100 µL conjugate to each well. Leave well A1 for substrate blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour at 37 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL diluted wash solution. Avoid overflows from the reaction wells.

**Important note:** During each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

**Automatic washer:** In case you use automatic equipment, wash the wells at least 5 times.

*Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.*

5. Dispense 100 µL TMB Substrate Solution into all wells.
6. **Incubate for exactly 20 min at room temperature (+22 - +28 °C) in the dark.**
7. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution. Shake the microplate gently.  
*Any blue colour developed during the incubation turns into yellow.*
8. Measure the absorbance (E) of the specimen against blank at 450 nm against a reference wavelength of 620-630 nm or against blank within 5 minutes.

## 8.2 Measurement

Adjust the ELISA Microwell Plate Reader **to zero** using the **substrate blank in well A1**.

*If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!*

**Measure the absorbance** of all wells at **450 nm** and record the absorbance values for each standard and sample.

Where applicable calculate the **mean absorbance values** of all duplicates.

## 9. RESULTS

---

### 9.1 Calculation of Results

Calculate the mean absorbance for each point of the standard curve and each sample. Plot the mean value of absorbance of the standards against concentration. Draw the best-fit curve through the plotted points. (e. g.: Four Parameter Logistic). Interpolate the values of the samples on the standard curve to obtain the corresponding values of the concentrations expressed in pg/mL.

### 9.2 Reference Values

To determine the normal range for serum / plasma samples, EDTA plasma samples from 150 apparently healthy male and female adults have been tested. The posture of the patients, whether they were upright or supine before the collection of the specimen, is unknown.

Result:

Serum / plasma normal range
< 14.6 - 174 pg/mL

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore, each laboratory should consider the range given by the manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

## 10. QUALITY CONTROL

---

Each laboratory should assay controls at normal, high and low levels range of Aldosterone for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the standard curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

## 11. PERFORMANCE AND CHARACTERISTICS

### 11.1 Analytical Sensitivity

Analytical sensitivity has been investigated through the LOB (Limit of Blank), LOD (Limit of Detection), LOQ (limit of Quantitation) and Analytical Sensitivity (A.S.).

The table below shows the criteria of the study and the results obtained.

	Criteria	Results (pg/mL)
LOB	60 replicates of Cal 0, used as "Blank", have been investigated in 9 different sessions over 5 days	15.3
LOD	6 low serum samples have been analyzed in 10 different sessions over 5 days	41.5
LOQ	6 low serum samples have been analyzed in 10 different sessions over 5 days	79.1
A.S.	20 replicates of Cal 0 and 5 replicates of Cal 2 have been assayed. A.S. has been calculated by linear regression.	14.6

### 11.2 Precision and Reproducibility (Complex Precision)

To determine precision and reproducibility 5 different serum samples have been used.

The table below shows the Within Run and Total CV%.

Sample	n°	Mean (pg/mL)	Within Run CV%	Total CV%
PS1	20	1410.102	5.1%	7.2%
PS2	20	824.239	5.2%	8.4%
PS3	20	483.272	5.7%	8.4%
PS4	20	257.171	6.8%	12.3%
PS5	20	127.102	10.4%	13.9%

### 11.3 Sample type study

Sample type study has been evaluated by performing the assay in parallel with serum and plasma belonging to the same subject. 20 different subjects have been used.

Serum obtained with SST tubes, and plasma obtained with EDTA, lithium-heparin and sodium-heparin have been investigated.

The interference has been evaluated as "significant" if it causes a concentration bias >10% between the sample and the reference (EDTA plasma).

The table below shows the results obtained:

Sample	Bias	Interference
EDTA plasma	reference	/
Serum	9.1%	No
Serum (SST pipes)	6.4%	No
Lithium heparin plasma	5.1%	No
Sodium heparin plasma	7.6%	No

Conclusion: following the study, serum and plasma samples can be used indifferently with this kit.

### 11.4 Analytical Specificity

#### 11.4.1 Interfering Substances

Interference by conjugated Bilirubin, Hemoglobin and Triglycerides has been investigated by adding the interfering substance to the serum sample and by comparing its concentration to the unspiked sample.

Samples with low and high concentration of Aldosterone have been analyzed.

The interference has been evaluated as "significant" if it causes a concentration (bias) >10% between spiked and unspiked sample.

The table below shows the results obtained:

Substance	Conc. assayed	Interference
Triglycerides	600 mg/dL	No
Bilirubin, conjugated	33.1 mg/dL	No
Hemoglobin	6.6 mg/dL	No

Conclusion: following the study, there is no significant interference by conjugated Bilirubin, Hemoglobin and Triglycerides at the concentration assayed.

## 11.4.2 Cross-Reactivity

The following cross-reactants have been spiked in serum samples with low and high Aldosterone concentration. The respective low and high unspiked samples were used as reference.

The interference has been evaluated as "significant" if it causes a concentration bias >10 % between spiked and unspiked sample. The table below shows the results obtained:

Cross-reactant	Conc. assayed	Interference
Corticosterone	0.1 µg/mL	No
Androsterone	5 µg/mL	No
Cortisone	0.1 µg/mL	No
Diidrotestosterone (DHT)	0.5 µg/mL	No
Estradiol	10 µg/mL	No
Estrone	10 µg/mL	No
Estriol	10 µg/mL	No
Testosterone	0.2 µg/mL	No
DHEA-S	10 µg/mL	No

## 11.5 Accuracy

### 11.5.1 Linearity

Three tests have been performed using 3 pair of serum (with low and high Aldosterone concentration); 11 dilutions spanning the range low to high have been analyzed.

The assay is linear over the analyzed range (approximately 85.8 – 1520.9 pg/mL).

### 11.5.2 Recovery

Three tests have been performed using 3 serum samples with low aldosterone concentration, each one spiked with the Standard 5 in different percentages.

The mean recovery in the three samples has been 109 %, 102 % and 111 %.

## 11.6 Correlation

126 EDTA plasma samples have been tested with the Aldosterone ELISA kit and with the IDS - iSYS Aldosterone assay (reference method).

The linear regression curve is:

$$Y = 1.106 \cdot X - 36.06$$

$$r^2 = 0.94$$

## 12. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use by professional persons. Not for internal or external use in humans or animals.
- All components of human origin used for the production of these reagents have been tested for anti-HIV 1+2 antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive. Nevertheless, all materials should still be regarded and handled as potentially infectious.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- Do not use heavily haemolysed or highly lipemic samples.
- Maximum precision is required for dispensation of the reagents.
- Avoid the exposure of TMB substrate to direct sunlight, metal or oxidants.
- This method allows the determination of Aldosterone from 41,5 pg/mL (LOD) to 2000 pg/mL.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.

- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Microbiologically contaminated samples should not be used in the assay. Highly lipemic or haemolysed specimens should also not be used.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results.
- Addition of the substrate solution initiates a kinetic reaction, which is terminated by the addition of the stop solution. Therefore, the addition of the substrate and the stopping solution should be added in the same sequence to eliminate any time deviation during reaction
- The ELISA is only designed for qualified personnel who are familiar with good laboratory practice.
- The treatment with natural or synthetic steroids can affect the blood levels of Aldosterone

### 12.1 Safety note for reagents containing hazardous substances

Reagents may contain CMIT/MIT (3:1).

Therefore, the following hazard and precautionary statements apply.

**Warning**

H317	May cause an allergic skin reaction.
P261	Avoid breathing spray.
P280	Wear protective gloves/ protective clothing.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P362+P364	Take off contaminated and Wash it before reuse.

Further information can be found in the safety data sheet

### 12.2 Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

For information about the packaging materials refer to PACKAGING MATERIALS.

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
## BIBLIOGRAPHY / LITERATUR / BIBLIOGRAPHIE / BIBLIOGRAFÍA / BIBLIOGRAFIA

- Himathongkam, T. et al., Potassium-Aldosterone-Renin Interrelationships. *J. Clin. Endocrinol. Metab.* 41/1 :153-159, 1975.
- Lun, S. et al., A Direct Radioimmunoassay for Aldosterone in Plasma. *Clin. Chem.* 29/2:268-271, 1983.
- Carledge, S. and Lowson, N., Aldosterone and Renin Measurements. *Ann. Clin. Biochem.* 37:262-278, 2000.
- Sequeira, S. J. et al., Evaluation of an Aldosterone Radioimmunoassay: The Renin-Angiotensin-Aldosterone Axis as a function of Sex and Age. *Ann. Clin. Chem.* 37/11:1987-1989, 1991.
- Miller, M.A. et al., Extraction Method and Nonextracted Kit Comparison for measuring Plasma Aldosterone. *Clin. Chem.* 43/10: 1995-1997, 1997.
- Stabler and Siegel, A.L., Chemoluminescence Immunoassay of Aldosterone in serum. *Clin. Chem.* 37/11: 1987-1989, 1991.
- Vallotton M. B., Primary Aldosteronism. Part 1. Diagnosis of Primary Hyperaldosteronism. *Clin. Endocrinol.* 45: 47-52, 1996.
- Oelkers, W. et al., Diagnosis, Therapy Surveillance in Addison's Disease: Rapid Adrenocorticotropin (ACTH) Test and Measurement of Plasma ACTH, Renin Activity, and Aldosterone . *J. Clin. Endocrinol Metab.* 75: 259-264, 1992.
- Ad Dujaili, E.A.S. and Edwards, C. R. W., Optimization of a direct Radioimmunoassay for Plasma Aldosterone. *J. Steroid Biochem.* 14: 481-487, 1981.
- Corry, D. B. and Tuck, M. L., Secondary Aldosteronism. *Endocrinol. Metab. Clin. North Am.* 24: 511-528, 1995.
- Raizman E. et al, *Clinical Chemistry* 61:8, 1022-1027 (2015)
- Funder JW et al, *J. Clin. Endocrinol Metab.*, 93(9):3266-3281 (2008)

## PACKAGING MATERIALS / VERPACKUNGSMATERIALIEN / MATERIELS D'EMBALLAGE / MATERIALES DE EMBALAJE / MATERIAIS DE EMBALAGEM

 <p>PAP 21</p>	 <p>PAP 21</p>	 <p>PAP 22</p>	<div style="border: 1px solid black; padding: 2px; text-align: center; width: fit-content; margin: 0 auto;">MTP</div>  <p>ALU / LDPE 90</p>
---	---	--	--

**SYMBOLS KEY / SYMBOLSCHLÜSSEL / EXPLICATION DES SYMBOLES / SÍMBOLOS / TABELA DE SÍMBOLOS**

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Aldosterone

### Test Preparation

Prepare reagents and samples as described.  
 Establish the distribution and identification plan for all specimens and controls.  
 Select the required number of microtiter strips or wells and insert them into the holder.

### Assay Procedure

	Substrate blank	Standard 0 - 5	Control A + B	Sample
Standard 0 - 5	-	50 µL	-	-
Control A + B	-	-	50 µL	-
Sample	-	-	-	50 µL
Conjugate	-	100 µL	100 µL	100 µL
Cover wells with foil supplied in the kit <b>Incubate for 1 hour at 37 °C</b> Wash each well three times with <b>300 µL</b> diluted Wash Solution <b>In case you use automatic equipment, wash the wells at least 5 times.</b>				
TMB Substrate	100 µL	100 µL	100 µL	100 µL
<b>Incubate for exactly 20 min at room temperature in the dark</b>				
Stop Solution	100 µL	100 µL	100 µL	100 µL
Shake the microplate gently Photometric measurement at 450 nm, 620-630 nm				

**Gold Standard Diagnostics Frankfurt GmbH**

Waldstrasse 23 A6

63128 Dietzenbach, Germany

Tel.: +49 6074 23698-0

Fax: +49 6074 23698-900

E-Mail: [info.frankfurt@eu.goldstandarddiagnostics.com](mailto:info.frankfurt@eu.goldstandarddiagnostics.com)

Website: [clinical.goldstandarddiagnostics.com](http://clinical.goldstandarddiagnostics.com)

Date of issue: 2024-04-24  
 DNOV012\_ALD\_IFU\_rev01\_fromLot\_6086AN