LH

CE

Only for in-vitro diagnostic use

Instructions for use / Notice d'utilisation / Instrucciones de uso

English	2
Français	7
Español	13
Bibliography / Bibliographie / Bibliografía	19
Abbreviations / Abréviations / Abreviaciones	19
Packaging materials / Matériels d'emballage / Materiales de embalaje	19
Symbols Key / Explication des Symboles / Símbolos.....	20
Summary of Test Procedure / Résumé de la procédure de test / Resumen de la técnica	21

REF

DNOV030 (96 Determinations)

ENGLISH

1. INTRODUCTION

Luteinizing hormone (LH) is a glycoprotein consisting of two subunits with a molecular mass of 30,000 daltons. The α -subunit is similar to other pituitary hormones [follicle stimulating hormone (FSH), thyroid stimulating hormone (TSH) and chorionic gonadotropin (HCG)] while the β -subunit is unique. The β -subunit confers the biological activity to the molecule.


The α -subunit consists of 89 amino acid residues while the β -subunit contains 129 amino acids. The carbohydrate content is between 15% and 30%. The clinical usefulness of the measurement of luteinizing hormone (LH) in ascertaining the homeostasis of fertility regulation via the hypothalamic-pituitary-gonadal axis has been well established^{1,2}. In addition, the advent of in vitro fertilization (IVF) technology to overcome infertility associated problems has provided the impetus for rapid improvement in LH assay methodology from the technically demanding bioassay³ to the procedurally simple and rapid immunoenzymometric assays.

2. INTENDED USE

Immunoenzymatic colorimetric method (ELISA) for quantitative determination of LH in human serum or plasma from an adult population.

3. PRINCIPLE OF THE ASSAY

In this method, LH standards, patient specimens and/or controls (containing the native antigen) are first added to streptavidin coated wells. Biotinylated monoclonal and enzyme labelled antibodies are added and the reactants mixed: these antibodies have high affinity and specificity and are directed against distinct and different epitopes of LH. Reaction between the various LH antibodies and native LH occurs in the microwells without competition or steric hindrance forming a soluble sandwich complex. The interaction is illustrated by the following equation:

BtnAb(m) = Biotinylated Monoclonal Antibody (Excess Quantity)

AgLH = Native Antigen (Variable Quantity)

EnzAb(p) = Enzyme labelled Antibody (Excess Quantity)

EnzAb(p)-AgLH-BtnAb(m) = Antigen-Antibodies Sandwich Complex

ka = Rate Constant of Association

k-a = Rate Constant of Dissociation

Simultaneously, the complex is deposited to the well through the high affinity reaction of streptavidin and biotinylated antibody. This interaction illustrated below:

Streptavidin C.W. = Streptavidin immobilized on well.

Immobilized complex = Antibodies-Antigen sandwich bound.

After equilibrium is attained, the antibody-bound fraction is separated from unbound antigen by decantation or aspiration. The enzyme activity in the antibody-bound fraction is directly proportional to the native antigen concentration. The activity of the enzyme present on the surface of the well quantitated by reaction with a suitable substrate to produce colour. By utilizing several different calibrators of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 breakapart 8-well snap-off strips coated with Streptavidin; in resealable aluminium foil.
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.15 mol/L (avoid any skin contact), ready to use.
- **Conjugate:** 1 bottle containing 12 mL of horseradish peroxidase labelled anti-LH antibodies and Biotin labelled anti-LH antibodies, contains < 0.06 % (v/v) CMIT/MIT (3:1), BSA 0.1 %.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (H₂O₂-TMB 0.26g/L) (avoid any skin contact); contains < 0.0015 % (v/v) CMIT/MIT (3:1).
- **Wash solution 50x conc:** 1 bottle containing 20 mL NaCl 45 g/L, Tween-20 55 g/L; contains < 0.06 % (v/v) CMIT/MIT (3:1).
- **Control:** 1 bottle containing 1 mL of a lot-specific control solution. The concentration is mentioned on the label; contains < 0.06 % (v/v) CMIT/MIT (3:1), BSA 1 %.
- **Standards:** 6 bottles, 1 mL each; contain < 0.06 % (v/v) CMIT/MIT (3:1), BSA 1 % and have approx. the following concentrations:
 - Standard 0: 0 mIU/mL
 - Standard 1: 5 mIU/mL
 - Standard 2: 25 mIU/mL
 - Standard 3: 50 mIU/mL
 - Standard 4: 100 mIU/mL
 - Standard 5: 200 mIU/mL

For hazard and precautionary statements see 13.1.

4.2. Materials supplied

- 1 Cover foil
- 1 Instructions for use (IFU)

4.3. Materials and Equipment needed

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450, 620-630 nm
- Pipettes to deliver 50 µL with a precision better than 1.5%
- Dispenser(s) for repetitive deliveries of 0.100 mL and 0.300 mL volumes with a precision of better than 1.5%.
- Microplate automatic washer or a squeeze bottle (optional).
- Distilled or deionized water and a clean one litre cylinder for the dilution of the wash solution
- Absorbent Paper for blotting the microplate wells.
- Vacuum aspirator (optional) for wash steps.
- Timer
- Quality control materials.

5. STABILITY AND STORAGE

The originally closed reagents are stable up to the expiry date stated on the label when stored at +2...+8 °C in the dark.

6. REAGENT PREPARATION

It is very important to bring all reagents, samples and standards to room temperature (+22...+28°C) for at least 30 minutes before starting the test run! At the end of the assay, store immediately the reagents at 2-8°C; avoid long exposure to room temperature.

6.1. Microtiterplate

The ready to use break apart snap-off strips are coated with Streptavidin. *Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at +2...+8 °.*

6.2. Conjugate

The conjugate is ready to use.

6.3. Control

The bottle contains 1 mL of a lot-specific control solution. The concentration is indicated on the label. Before use, mix for 5 minutes with a rotating mixer. The control is ready to use.

6.4. Standards

Before use, mix for 5 minutes with a rotating mixer. The standards are ready to use. Once opened the standards are stable for 6 months at 2...8° C.

6.5. TMB Substrate Solution

The bottle contains 15 mL of a tetramethylbenzidine/hydrogen peroxide system. The reagent is ready to use and has to be stored at 2...8°C in the dark. *The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.*

6.6. Wash Solution

Dilute contents of wash buffer concentrate 50x to 1000 mL with distilled or deionised water in a suitable storage container. For smaller volumes respect the 1:50 ratio. The diluted wash solution is stable for 30 days at 2...8°C.

7. SPECIMEN COLLECTION AND PREPARATION

- The specimens shall be blood, serum in type and the usual precautions in the collection of venipuncture samples should be followed.
- For accurate comparison to established normal values, a fasting morning serum sample should be obtained.
- The blood should be collected in a venipuncture tube without additives or anti-coagulants. Allow the blood to clot. Centrifuge the specimen to separate the serum from the cells.
- Samples may be refrigerated at +2...+8°C for a maximum period of 5 days. If the specimen(s) cannot be assayed within this time, the sample(s) may be stored at -20°C for up to 30 days.
- Avoid repetitive freezing and thawing.
- When assayed in duplicate, 40 µL of the specimen is required.

7.1. Precaution

- Maximum precision is required for dispensation of the reagents.
- This method allows the determination of LH from 5 mIU/mL to 200 mIU/mL.
- Samples with concentration of LH higher than 200 mIU/mL have to be diluted with standard 0.

8. ASSAY PROCEDURE

8.1. Test Preparation

Please read the test instructions for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instructions for use as described Prior to commencing the assay, the distribution and identification plan for all specimens and standards should be carefully established. Select the required number of microtiter strips or wells and insert them into the holder.

Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than one plate is used, it is recommended to repeat the dose response curve. Please allocate at least:

1 well	(e.g. A1)	for the substrate blank
2 wells	(e.g. B1+C1)	for standard 0
2 wells	(e.g. D1+E1)	for standard 1
2 wells	(e.g. F1+G1)	for standard 2
2 wells	(e.g. H1+A2)	for standard 3
2 wells	(e.g. B2+C2)	for standard 4
2 wells	(e.g. D2+E2)	for standard 5
2 wells	(e.g. F2+G2)	for LH control

It is recommended to determine standards and patient samples in duplicate.

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each standard and each patient sample.

1. Dispense 20 µL standards, control and samples into their respective wells. Add 100 µL conjugate to each well. Leave well A1 for substrate blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour at room temperature (22...28°C).**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL diluted wash solution. Avoid overflows from the reaction wells.

Important note: During each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: In case you use automatic equipment, wash the wells at least 5 times.

Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.

5. Dispense 100 µL TMB Substrate Solution into all wells.
6. **Incubate for exactly 15 min at room temperature (22...28°C) in the dark.**
7. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution. Shake the microplate gently.
Any blue colour developed during the incubation turns into yellow.
8. Measure the absorbance of the specimen at 450 nm against a reference wavelength of 620-630 nm or against blank within 5 minutes.

9. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The control provided in the kit should be tested as unknown and are intended to assist in assessing the validity of results obtained with each assay plate.

We recommend the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the Certificate of Analysis (CoA) should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories⁸.

10. RESULTS

10.1. Calculation of results

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) curve fit, **including Standard 0 is required**. A smoothed spline fit including Standard 0 can be used. Other curve fitting algorithms are not recommended.

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Standard 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

In order for the assay results to be considered valid the kit standards and control must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid, and samples must be retested.

10.2. Reference values

The serum LH values are comprised in the following intervals:

Men	0.7 – 7.4 mIU/mL
Women	
follicular phase	0.5 – 10.5 mIU/mL
Ovulation phase	18.4 – 61.2 mIU/mL
luteal phase	0.5 – 10.5 mIU/mL
menopause	8.2 – 40.8 mIU/mL

It is important to keep in mind that establishment of a range of values which can be expected to be found by a given method for a population of "normal"-persons is dependent upon a multiplicity of factors: the specificity of the method, the population tested and the precision of the method in the hands of the analyst.

For these reasons each laboratory should depend upon the range of expected values established by the Manufacturer only until an in-house range can be determined by the analysts using the method with a population indigenous to the area in which the laboratory is located.

11. SPECIFIC PERFORMANCE CHARACTERISTICS

11.1. Precision

Intra Assay Variation

Within-run precision was determined by replicate determinations (20x) of three different control sera in one assay. The within-assay variability is $\leq 9.21\%$

Inter Assay Variation

Between-run precision was determined by replicate measurements of three different control sera in 15 different assays. The between-assay variability is $\leq 7.91\%$

11.2. Accuracy

The recovery test performed on three different samples, enriched with 5.63 - 11.25 - 22.5 - 45 - 90 mIU/mL of LH, gave a average value (\pm SD) of $97.17\% \pm 4.00\%$.

In the dilution test three different samples were diluted 2, 4, 8 and 16 times with Standard 0; the average value (\pm SD) obtained is $99.13\% \pm 7.37\%$.

11.3. Specificity

The cross-reactivity of the LH kit to selected substances was evaluated by adding the interfering substance to a serum matrix at various concentrations. The cross-reactivity was calculated by deriving a ratio between dose of interfering substance to dose of Luteinizing Hormone needed to produce the same absorbance.

Substance	Cross Reactivity
LH	100%
HCG	0.007%
Prolactin	None detected
FSH	None detected
TSH	None detected

11.4. Sensitivity

The minimal detectable concentration of LH by this assay is estimated to be 0.22 mIU/mL.

11.5. Correlation

The LH kit was compared to a commercially available LH kit. 36 serum samples were tested.

The regression curve is:

$$\text{GSD} = 0.91 * (\text{commercial kit}) + 0.05$$

$$R^2 = 0.971$$

The new LH kit was compared to the old LH kit. 36 serum samples were tested.

The regression curve is:

$$(\text{new}) = 1.08 * (\text{old}) - 1.22$$

$$R^2 = 0.981$$

11.6. Hook Effect

The LH assay shows no Hook Effect up to 400 mIU/mL.

12. LIMITATIONS OF THE PROCEDURE

- Bacterial contamination or repeated freeze-thaw cycles of the specimen may affect the absorbance values.
- LH is suppressed by estrogen but in woman taking oral contraceptives the level may be low or normal. Excessive dieting and weight loss may lead to low gonadotropin concentrations. Luteinising hormone is dependent upon diverse factors other than pituitary homeostasis; thus, the determination alone is not sufficient to assess clinical status.
- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.

- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays⁹. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

13. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use by professional persons. Not for internal or external use in Humans or Animals.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive. Nevertheless, all materials should still be regarded and handled as potentially infectious.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- Some reagents (standards, control, conjugate, and wash solution) contain small amounts of CMIT/MIT as preservative. Avoid contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Avoid the exposure of TMB substrate to direct sunlight, metal or oxidants. Do not freeze the solution.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- The ELISA is only designed for qualified personnel who are familiar with good laboratory practice.

13.1. Safety note for reagents containing hazardous substances

Reagents may contain CMIT/MIT (3:1) (refer to 4.1)

Therefore, the following hazard and precautionary statements apply.

Warning 	H317	May cause an allergic skin reaction.
	P261	Avoid breathing gas/mist/vapours/spray.
	P280	Wear protective gloves/protective clothing/eye protection/face protection.
	P302+P352	IF ON SKIN: Wash with plenty of water.
	P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
	P362+P364	Take off contaminated clothing and wash it before reuse.

Further information can be found in the safety data sheet.

13.2. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

For information about the packaging materials refer to PACKAGING MATERIALS.

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(96 Determinaciones)


BIBLIOGRAPHY / BIBLIOGRAPHIE / BIBLIOGRAFIA

1. Kosasa T.S., "Measurement of Human Luteinizing Hormone." *Journal of Reproductive Medicine*, 26 (1981) pg. 201-6.
2. Danzer H., Braunstein G.D., et al., "Maternal Serum Human Chorionic Gonadotropic Concentrations and Fetal Sex Predictions." *Fertility and Sterility*, 34 (1980) pg. 336-40.
3. Braunstein G.D., et al., "Serum Human Luteinizing Hormone Levels through Normal Pregnancy ", *American Journal of Obstetrics and Gynecology*, 126 (1976) pg. 678-81.
4. Goldstein D.P., and Kosasa T.S., "The Subunit Radioimmunoassay for LH Clinical Application." *Gynecology*, 6 (1975) pg. 145-84.
5. Batzer F., "Hormonal Evaluation of Early Pregnancy", *Fertility and Sterility*, 34 (1980) pg. 1-12.
6. Braunstein, G.D., et al., "First-Trimester Luteinizing Hormone Measurements as an Aid to the Diagnosis of Early Pregnancy Disorders", *American Journal of Obstetrics and Gynecology*, 131 (1978) pg. 25-32.
7. Lenton E., Neal L. and Sulaiman R., "Plasma Concentrations of Human Gonadotropin from the time of Implantation until the Second Week of Pregnancy ", *Fertility and Sterility*, 37 (1982) pg. 773-78.
8. Basic QC Practices On-line Course; <http://www.Westgard.com>.
9. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27–33

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
LH

Test Preparation

Prepare reagents and samples as described.
 Establish the distribution and identification plan for all specimens and controls.
 Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate blank	Standard 0 - 5	Control	Sample
Standard 0 - 5	-	20 µL	-	-
Control	-	-	20 µL	-
Sample	-	-	-	20 µL
Conjugate	-	100 µL	100 µL	100 µL
Incubate for 1 h at room temperature (+22 - +28 °C) Wash each well three times with 300 µL diluted wash solution In case you use automatic equipment, wash the wells at least 5 times				
TMB Substrate	100 µL	100 µL	100 µL	100 µL
Incubate for exactly 15 min at room temperature in the dark				
Stop Solution	100 µL	100 µL	100 µL	100 µL
Shake the microplate gently Photometric measurement at 450 nm, 620-630 nm				

Gold Standard Diagnostics Frankfurt GmbH

Waldstrasse 23 A6
 63128 Dietzenbach, Germany
 Tel.: +49 6074 23698-0
 Fax: +49 6074 23698-900
 E-Mail: info.frankfurt@eu.goldstandarddiagnostics.com
 Website: clinical.goldstandarddiagnostics.com

Date of issue: 2025-07-29
 DNOV030_LH_IFU_rev02_fromLot_6262D