beta-HCG

Enzyme immunoassay for the quantitative determination of beta-HCG in human serum or plasma

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ENGLISH

1. INTRODUCTION

Human chorionic gonadotropin (hCG) is a glycoprotein hormone secreted in pregnancy, that is made by the embryo soon after conception and later by the syncytiotrophoblast (part of the placenta). Its role is to prevent the disintegration of the corpus luteum of the ovary and thereby maintain progesterone production that is critical for a pregnancy in humans. hCG may have additional functions, for instance it is thought that it affects the immune tolerance of the pregnancy.

Pregnancy tests measure the levels of hCG in the blood or urine to indicate the presence or absence of an implanted embryo. In particular, pregnancy tests employ an antibody that is specific to the β -subunit of hCG (β hCG). This is important so that tests do not make false positives by confusing hCG with LH and FSH.

β hCG is also secreted by some cancers including teratomas, choriocarcinomas. But, elevated levels cannot prove the presence of a tumor, and low levels do not rule it out.

2. INTENDED USE

Immunoenzymatic colorimetric method (ELISA) for quantitative determination of beta-HCG in human serum or plasma.

3. PRINCIPLE OF THE ASSAY

The beta-HCG assay is based on the simultaneous capture of HCG by a monoclonal antibody immobilized on the microplate and directed against the β -HCG fraction, and another monoclonal antibody conjugated with peroxidase horseradish (HRP) and directed against the α -HCG fraction. After the incubation, the bound/free separation is performed by a simple solid-phase washing. The enzyme HRP in the bound-fraction reacts with the substrate (H_2O_2) and the TMB substrate and develops a blue color that changes into yellow when the stop solution (H_2SO_4) is added. The color intensity is proportional to the β -HCG concentration in the samples. The β -HCG concentration in the sample is calculated based on a calibration curve.

4. MATERIALS

4.1. Reagents supplied

- **Anti- β -HCG IgG Coated Wells:** 12 breakapart 8-well snap-off strips coated monoclonal antibody direct against β -HCG subunit, in resealable aluminium foil.
- **Stop Solution:** 1 bottle containing 15 ml sulphuric acid, 0.15 mol/l (avoid any skin contact).
- **Conjugate (conc.):** 1 bottle containing 1 ml of monoclonal antibody direct against α -HCG subunit conjugated with horseradish peroxidase (HRP)
- **Incubation buffer:** 1 bottle containing 50 ml phosphate buffer 50 mM, pH 7.4; BSA 1 g/l.
- **TMB Substrate Solution:** 1 bottle containing 15 ml 3, 3', 5, 5'-tetramethylbenzidine (H_2O_2 -TMB 0.26g/l) (avoid any skin contact).
- **Wash solution 50x conc.:** 1 bottle containing 20 ml (NaCl 45 g/l, Tween-20 55 g/l)
- **β -HCG Control:** 1 bottle containing 1 ml of a lot-specific, ready to use control. The concentration is mentioned on the label.
- **β -HCG Standards:** 6 bottles, 1 ml each. The standards are calibrated against the WHO 1st IRP 75/537 and have the following concentrations:

Standard 0:	0 mIU/ml
Standard 1:	1 mIU/ml
Standard 2:	5 mIU/ml
Standard 3:	20 mIU/ml
Standard 4:	100 mIU/ml
Standard 5:	400 mIU/ml

4.2. Materials supplied

- 1 Cover foil
- 1 Instructions for use (IFU)

4.3. Materials and Equipment needed

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450 nm, 620-630 nm
- Manual or automatic equipment for rinsing wells
- Pipettes to deliver volumes between 10 and 1000 μ l
- Rotating mixer
- Distilled water
- Disposable tubes
- Timer

5. STABILITY AND STORAGE

The closed reagents are stable up to the expiry date stated on the label when stored at 2...8 °C in the dark.

6. REAGENT PREPARATION

It is very important to bring all reagents, samples and standards to room temperature (22...28°C) for at least 30 minutes before starting the test run! At the end of the assay, store the reagents immediately at 2-8°C; avoid long exposure at room temperature.

6.1. Coated snap-off Strips

The ready to use break apart snap-off strips are coated with monoclonal antibody direct against β -HCG subunit. Store at 2...8 °C. *Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C; stability until expiry date.*

6.2. Conjugate

The bottle contains 1 ml of a concentrated solution with monoclonal antibody direct against α -HCG subunit conjugated with horseradish peroxidase (HRP)

Dilute immediately before use:

Add 10 μ l concentrated conjugate to 1.0 ml of Incubation Buffer. The quantity of diluted conjugate is proportional to the number of tests.

Mix gently with a rotating mixer for 10 min. Stable for 3 hours at room temperature (22-28° C).

6.3. Standards

Each of the 6 vials contains 1 ml standard solution of the concentration mentioned in 4.1. The standards are ready to use.

After first use the standard solutions are still stable for another 6 months if stored at 2...8 °C.

6.4. β -HCG control

The bottle contains 1 ml of a lot-specific control solution. The concentration is indicated on the label. The control is ready to use.

6.5. TMB Substrate Solution

The bottle contains 15 ml of a tetramethylbenzidine/hydrogen peroxide system. The reagent is ready to use and has to be stored at 2...8°C in the dark. *The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away. After first use the TMB substrate solution is still stable for another 6 months if stored at 2...8 °C.*

6.6. Stop Solution

The bottle contains 15 ml 0.15 M sulphuric acid solution (H314). This ready to use solution has to be stored at 2...8°C.

After first use stable until expiry date.

6.7 Wash Solution

Dilute content of concentrated wash solution to 1000 ml with distilled water in a suitable storage container. For smaller volumes respect the 1:50 ratio. The diluted wash solution is stable for 30 days at 2...8 °C.

7. SPECIMEN COLLECTION AND PREPARATION

Use human serum or plasma samples with this assay. If the assay is performed within 48 hours after sample collection, the specimens should be kept at 2...8°C; otherwise, they should be aliquoted and stored deep-frozen (-20 to -70°C). If samples are stored frozen, mix thawed samples well before testing. *Avoid repeated freezing and thawing.*

Samples with concentration of β -HCG over 400 mIU/ml have to be diluted with Incubation buffer.

7.1. Precaution

- Avoid the exposure of TMB substrate to direct sunlight, metal or oxidants. Do not freeze the solution.
- Maximum precision is required for dispensation of the reagents.
- This method allows the determination of beta-HCG from 1 to 400 mIU/ml.

8. ASSAY PROCEDURE

8.1. Test Preparation

Please read the test instructions for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instructions for use as described. Prior to commencing the assay, the distribution and identification plan for all specimens and standards should be carefully established. Select the required number of microtiter strips or wells and insert them into the holder. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than one plate is used, it is recommended to repeat the dose response curve.

Please allocate at least:

1 well	(e.g. A1)	for the substrate blank
2 wells	(e.g. B1+C1)	for standard 0
2 wells	(e.g. D1+E1)	for standard 1
2 wells	(e.g. F1+G1)	for standard 2
2 wells	(e.g. H1+A2)	for standard 3
2 wells	(e.g. B2+C2)	for standard 4
2 wells	(e.g. D2+E2)	for standard 5
2 wells	(e.g. F2+G2)	for control

It is recommended to determine standards, control and patient samples in duplicate.

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each standard and each patient sample.

1. Dispense 25 µl standards, control and samples into their respective wells. Add 100 µl diluted conjugate to each well. Leave well A1 for substrate blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour at room temperature.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µl diluted wash solution. Avoid overflows from the reaction wells. The soak time between each wash cycle should be >5sec. During each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbing paper towel. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step! (if you use automated equipment, wash the wells at least 5 times).

Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.

5. Dispense 100 µl TMB Substrate Solution into all wells.
6. **Incubate for exactly 15 min at room temperature in the dark.**
7. Dispense 100 µl Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution. Shake the microplate gently.
Any blue colour developed during the incubation turns into yellow.
8. Read the absorbance at 450 nm against a reference wavelength of 620 -630 nm or against blank within 5 minutes.

8.2. Measurement

Adjust the ELISA Microwell Plate Reader **to zero** using the **substrate blank in well A1**.

If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard and patient sample.

Where applicable calculate the **mean absorbance values** of all duplicates.

9. RESULTS

9.1. Calculation of results

Calculate the mean absorbance for each point of the standard curve and each sample. Plot the mean value of absorbance of the standards against concentration. Draw the best-fit curve through the plotted points. (Four Parameter Logistic or Sigmoide).

Interpolate the values of the samples on the standard curve to obtain the corresponding values of the concentrations expressed in mIU/ml.

9.2. Reference Values

Each laboratory must establish its own normal ranges based on patient population.

The serum or plasma β -HCG values are comprised in the following intervals:

Normal women:		< 8 mIU/ml
Pregnancy:	1. week	3 - 100 mIU/ml
	2. week	10 - 1000 mIU/ml
	3. week	100 - 10000 mIU/ml
	4. week	1000 - 100000 mIU/ml
	2. month	15000 - 200000 mIU/ml
	3. month	10000 - 100000 mIU/ml

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

10. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of β -HCG for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

11. SPECIFIC PERFORMANCE CHARACTERISTICS

11.1. Precision

Intra Assay Variation

Within run variation was determined by replicate determination (16x) of three different control sera in one assay. The within assay variability is $\leq 7.6\%$.

Inter Assay Variation

Between run variation was determined by replicate measurements of three different control sera in different lots. The between assay variability is $\leq 8.8\%$.

11.2. Specificity

The cross reactivity values of the β -HCG ELISA have been calculated on a weight/ weight basis.

β -HCG	100.0 %
hFSH	3.0 %
HCG	4.0 %
hTSH	0.02%

11.3. Analytic Sensitivity

The lowest detectable concentration of human β -HCG that can be distinguished from standard 0 is 0.09 mIU/ml at the 95% confidence limit.

11.4. Accuracy

The recovery of 6.25 – 12.5 – 25 – 50 mIU/ml of β -HCG added to sample gave an average value (\pm SD) of $99.2\% \pm 4.1\%$ with reference to the original concentrations.

11.5. Correlation with RIA

The NovaTec beta-HCG ELISA was compared to another commercially available beta-HCG assay. Serum samples of 49 females were analysed according in both test systems.

The linear regression curve was calculated:

$$y = 0.94x - 0.02$$

$$r = 0.96 \text{ (} r^2 = 0.92 \text{)}$$

11.6. Hook Effect

The beta-HCG ELISA, a competitive enzyme immunoassay, shows no Hook Effect up to 250.000 mIU/ml.

12. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the specimen may affect the absorbance values.

13. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use by professional persons. Not for internal or external use in Humans or Animals.
- Follow Good Laboratory Practice (GLP) for handling blood products.

- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive. Nevertheless, all materials should still be regarded and handled as potentially infectious.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- When using automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested. To improve the performance of the kit on ELISA automatic systems, it is recommended to increase the number of washes.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- The NovaTec ELISA is only designed for qualified personnel who are familiar with good laboratory practice.

WARNING: Some reagents contain small amounts of Proclin 300® as preservative. Avoid the contact with skin and mucous membranes!

13.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

For information about the packaging materials refer to PACKAGING MATERIALS.

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
Acosta, A.A.M.D. and Wright, G.L., Journal of Clinical Immunoassays (1983) 6, 41

Cohen, K.L., Metabolism (1977) 26, 1165 – 1177





Lindstedt (1982) Scand. J. Clin. Lab. Invest 42, 201

Chien (1987) J. Clin. Endocr. Metab. 64, 313

PACKAGING MATERIALS / MATÉRIELS D'EMBALLAGE / MATERIALES DE EMBALAJE

 PAP 21	 PAP 21	 PAP 22	<div style="border: 1px solid black; padding: 2px; text-align: center;">MTP</div>  ALU / LDPE 90
--	--	--	--

SYMBOLS KEY / EXPLICATION DES SYMBOLES / SIMBOLOS

	Manufactured by / Fabriqué par/ Fabricado por
	In Vitro Diagnostic Medical Device/ Dispositif médical de diagnostic <i>in vitro</i> / Producto para diagnóstico In vitro
	Lot Number / Numéro de lot / Número de lote
	Expiration Date / Date de péremption / Fecha de caducidad
	Storage Temperature / Température de conservation / Temperatura de almacenamiento
	Keep away from direct sun / Protéger de rayonnement solaire / Proteger de la luz solar.
	CE Mark/ Marquage CE / Marca CE
	Unique Device Identifier / identification unique des dispositifs / identificación única del producto
	Catalogue Number/ Référence du catalogue / Número de Catálogo
	Consult Instructions for Use / Consulter la notice d'utilisation / Consulte las Instrucciones de Uso
	Microplate / Microplaque/ Microplaca
	Conjugate conc. / Conjugué conc./ Conjugado conc.
	Control / contrôle / control
	Incubation buffer/ Tampon d'incubation / Buffer de incubación
	Calibrator resp. Standard / Callibrateur resp Etalon / Calibrador o bien Estándar
	Stop solution / Solution d'arrêt / Solución de parada
	TMB Substrate solution / Solution de TMB / Solución substrato TMB
	Washing solution 50x concentrated/ Solution de lavage concentré 50 x / solución de lavado concentrado x50
	Contains sufficient for "n" tests / Contenu suffisant pour "n" tests / Contenido suficiente para "n" tests

SCHEME OF THE ASSAY

beta-HCG

Test Preparation

Prepare reagents and samples as described.
Establish the distribution and identification plan for all specimens and controls.
Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate blank	Standard 0 - 5	Control	Sample
Standard 0 - 5	-	25 µl	-	-
Control	-	-	25 µl	-
Sample	-	-	-	25 µl
Diluted Conjugate	-	100 µl	100 µl	100 µl
Cover wells with foil supplied in the kit Incubate for 1 hour at room temperature (22...28°C) Wash each well three times with 300 µl diluted Wash Solution In case you use automated equipment, wash the wells at least 5 times.				
TMB Substrate	100 µl	100 µl	100 µl	100 µl
Incubate for exactly 15 min at room temperature in the dark				
Stop Solution	100 µl	100 µl	100 µl	100 µl
Shake the microplate gently Photometric measurement at 450 nm				

Gold Standard Diagnostics Frankfurt GmbH

Waldstrasse 23 A6

63128 Dietzenbach, Germany

Tel.: +49 6074 23698-0

Fax: +49 6074 23698-900

E-Mail: info.frankfurt@eu.goldstandarddiagnostics.com

Website: clinical.goldstandarddiagnostics.com

DNOV034-bHCG_IFU_rev01_fromLot_6183AN