Free T3

Only for in-vitro diagnostic use

Instructions for use / Notice d'utilisation / Instrucciones de uso

English	2
Français	7
Español	12
Bibliography / Bibliographie / Bibliografía	18
Abbreviations / Abréviations / Abreviaciones	18
Packaging materials/ Matériels d'emballage/ Materiales de embalaje	18
Symbols Key / Explication des symboles / Símbolos	19
Summary of Test Procedure / Résumé de la procédure de test / Resumen de la técnica	20

REF

DNOV051 (96 Determinations)

ENGLISH

1. INTRODUCTION

The thyroid hormone, triiodothyronine (T3), is produced by the thyroid gland. An important component in the synthesis is iodine. Thyroxine is converted to the active T3 (three to four times more potent than T4) within cells by deiodinases (5'-iodinase).

Thyroxine-binding globulin (TGB) is the major carrier protein for circulating thyroid hormone.

Only a very small fraction of the circulating hormone is free (unbound) 0.3%; this fraction is biologically active. Thus, measurements of free triiodothyronine concentrations correlate more reliably with clinical status than total triiodothyronine levels. For example, the increase in total triiodothyronine levels associated with pregnancy, oral contraceptives and oestrogen therapy result in higher total T3 levels while the free T3 (FT3) concentration remains basically unchanged.

The concentrations of the carrier proteins are altered in many clinical conditions, such as pregnancy. In normal thyroid function as the concentrations of the carrier proteins alters, the total triiodothyronine level changes so that the free triiodothyronine concentration remains constant.

The binding of T3 plays a key role in the feedback control of the thyroid, with FT3 acting on the pituitary to inhibit thyroid hormone secretion.

The thyronines act on the body to increase the basal metabolic rate, affect protein synthesis and increase the body's sensitivity to catecholamine (such as adrenaline) by permissiveness. The thyroid hormones are essential to proper development and differentiation of all cells of the human body. These hormones also regulate protein, fat, and carbohydrate metabolism, affecting how human cells use energetic compounds. Numerous physiological and pathological stimuli influence thyroid hormone synthesis. Thyrotoxicosis or hyperthyroidism is the clinical syndrome caused by an excess of circulating free thyroxine, free triiodothyronine, or both.

Both T3 and T4 are used to treat thyroid hormone deficiency (hypothyroidism).

Since conditions such as pregnancy, oestrogen therapy and other non-thyroid factors alter TBG concentrations, assessment of thyroid function through total T3 measurement may result in an erroneous diagnosis, because FT3 levels, are unaffected by binding protein changes.

2. INTENDED USE

Competitive immunoenzymatic colorimetric method for quantitative determination of free T3 (Triiodothyronine) in human serum or plasma from an adult population.

3. PRINCIPLE OF THE ASSAY

Microtiter strip wells are precoated with anti-T3 antibodies (solid phase). T3 in the sample competes with added horseradish peroxidase labelled T3 (enzyme-labelled antigen) for antibody binding. After incubation a bound/free separation is performed by solid-phase washing. The immune complex formed by enzyme-labelled antigen is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product. The intensity of this product is **inversely** proportional to the amount of T3 in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorption at 450 nm is read using an ELISA microwell plate reader.

FT3 concentration in the sample is calculated through a calibration curve.

4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 breakapart 8-well snap-off strips coated with anti-T3; in resealable aluminium foil.
- **Stop Solution:** 1 bottle containing 15 ml sulphuric acid, 0.15 mol/l (avoid any skin contact).
- **Conjugate:** 1 bottle containing 12 ml of horseradish peroxidase labelled T3.
- **Wash Solution 50x conc.:** 1 bottle containing 20 ml (NaCl 45 g/l, Tween-20 55 g/l); contains < 0.06 % (v/v) CMIT/MIT (3:1).
- **TMB Substrate Solution:** 1 bottle containing 15 ml 3,3',5,5'-tetramethylbenzidine (H₂O₂-TMB 0.26 g/l) (avoid any skin contact); contains ≤ 0.0015 % (v/v) CMIT/MIT (3:1).
- **Standards:** 6 bottles, 1 ml each
 - Standard 0: 0.0 pg/ml
 - Standard 1: 1.5 pg/ml
 - Standard 2: 3.1 pg/ml
 - Standard 3: 5.4 pg/ml
 - Standard 4: 8.5 pg/ml
 - Standard 5: 19.0 pg/ml

Exact concentrations are given on the labels on a lot specific basis.

For SI units: 1pg/ml x 1.536 = pmol/l.

For hazard and precautionary statements see 13.1.

4.2. Materials supplied

- 1 Cover foils
- 1 Instructions for use (IFU)

4.3. Materials and Equipment needed

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450 nm, 620-630 nm
- Manual or automatic equipment for rinsing wells
- Vortex tube mixer

- Distilled water
- Disposable tubes
- Timer

5. STABILITY AND STORAGE

The reagents are stable up to the expiry date stated on the label when stored in the original vial at 2...8 °C in the dark.

6. REAGENT PREPARATION

It is very important to bring all reagents, samples and standards to room temperature (22...28°C) for at least 30 minutes before starting the test run! At the end of the assay, store all reagents immediately at 2-8° C; avoid long exposure to room temperature.

6.1. Microtiterplate

The ready to use break apart snap-off strips are coated with anti-T3 IgG antibodies. Store at 2...8 °C. Open the bag only when it is at room temperature. *Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C; stability until expiry date.*

6.2. Conjugate

The T3-HRP conjugate is ready to use.

6.3. Wash Solution

Dilute concentrated Wash Solution to 1000 ml with distilled or deionised water in a suitable storage container. For smaller volumes respect the 1:50 ratio. The diluted wash solution is stable for 30 days at 2...8°C.

6.4. Standards

The standards are ready to use. Once opened the standards are stable for another 6 months at 2...8°C.

6.5. TMB Substrate Solution

The bottle contains 15 ml of a tetramethylbenzidine/hydrogen peroxide system. The reagent is ready to use and has to be stored at 2...8°C in the dark. *The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.*

6.6. Stop Solution

The bottle contains 15 ml 0.15 M sulphuric acid solution. This ready to use solution has to be stored at 2...8°C.

7. SPECIMEN COLLECTION AND PREPARATION

Specimen(s) may be refrigerated at 2...8°C for a maximum period of 48 hours. If the specimen(s) cannot be assayed within 48 hours, the sample(s) may be stored at temperatures of -20°C for up to 30 days. Avoid repetitive freezing and thawing of the samples.

When assayed in duplicate, 0.100 ml of the specimen is required.

7.1. Precaution

- Avoid the exposure of TMB substrate to direct sunlight, metal or oxidants. Do not freeze the solution.
- Maximum precision is required for dispensation of the reagents.

8. ASSAY PROCEDURE

8.1. Test Preparation

Please read the test instructions for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instructions for use as described. Prior to commencing the assay, the distribution and identification plan for all specimens and standards should be carefully established. Select the required number of microtiter strips or wells and insert them into the holder. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than one plate is used, it is recommended to repeat the dose response curve. Please allocate at least:

2 wells (e.g. A1+B1)	for standard 0
2 wells (e.g. C1+D1)	for standard 1
2 wells (e.g. E1+F1)	for standard 2
2 wells (e.g. G1+H1)	for standard 3
2 wells (e.g. A2+B2)	for standard 4
2 wells (e.g. C2+D2)	for standard 5

It is necessary to determine standards and patient samples in duplicate.

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each standard and each patient sample.

1. Dispense 50 µl standards and samples into their respective wells.
2. Add 100 µl T3-HRP Conjugate to each well.
3. Swirl the microplate gently for 20 – 30 sec to mix and cover wells with the foil supplied in the kit.
4. **Incubate for 1 hour at room temperature (22...28°C).**
5. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µl diluted Wash Solution. Avoid overflows from the reaction wells.

Important Note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: In case you use automatic equipment, wash the wells at least 5 times.

Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.

6. Dispense 100 µl TMB Substrate Solution into all wells.
7. **Incubate for exactly 15 min at room temperature (22...28°C) in the dark.**
8. Dispense 100 µl Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution. Mix gently for 15 – 20 sec.
Any blue colour developed during the incubation turns into yellow.
9. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against blank within 5 minutes.

8.2. Measurement

Adjust the ELISA Microwell Plate Reader **to zero** using the **standard 0**.

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard and patient sample.

Where applicable calculate the **mean absorbance values** of all duplicates.

9. RESULTS

9.1. Calculation of results

Calculate the mean absorbance for each point of the standard curve and each sample. Plot the mean value of absorbance of the standards against concentration. Draw the best-fit curve through the plotted points. (e.g. Four Parameter Logistic).

Interpolate the values of the samples on the standard curve to obtain the corresponding values of the concentrations expressed in pg/ml.

9.2. Reference values

A study of euthyroid adult population was undertaken to determine expected values for the Free T3 EIA Test System.

	Mean (pg/ml)	SD	Range (pg/ml)
Adult	2.8	0.7	1.4 – 4.2
Pregnancy	3.0	0.6	1.8 – 4.2

Please pay attention to the fact that the determination of a range of expected values for a “normal” population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore, each laboratory should consider the range given by the manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

10. QUALITY CONTROL

Each laboratory should assay controls at levels in the hypothyroid, euthyroid and hyperthyroid range for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the standard curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

In order for the assay results to be considered valid the kit calibrators and control must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid, and samples must be retested.

11. SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

11.1. Precision

Intra Assay Variation

Within run variation was determined by replicate determination (24x) of three different control sera in one assay. The within assay variability is ≤ 4.94 %.

Inter Assay Variation

Between run variation was determined by replicate measurements (12x) of three different control sera in different lots. The between assay variability is ≤ 13.19 %.

11.2. Cross Reactivity

The cross reactivity of the Triiodothyronine antibody to selected substances was evaluated by adding the interfering substance to a serum matrix at various concentrations. The cross-reactivity was calculated by deriving a ratio between doses of interfering substance to dose of Triiodothyronine needed to displace the same amount of tracer.

Substance	Concentration	Cross Reactivity
I-Triiodo-thyronine	-	1.0000
I-Thyroxine	10 µg/ml	< 0.0002
d-Thyroxine	10 µg/ml	< 0.0001
Iodo-thyrosine	10 µg/ml	< 0.0001
Diodo-thyrosine	10 µg/ml	< 0.0001
Triiodothyroacetic Acid	10 µg/ml	< 0.0001
Phenylbutazone	10 µg/ml	< 0.0001
Sodium Salicylate	10 µg/ml	N/D
Phenytoin	10 µg/ml	N/D
Oleic Acid	10 nmol/l	N/D
Albumin	50 mg/ml	N/D
Hemoglobin	10 µl/ml of packed red cells added to the serum	N/D

11.3. Analytic Sensitivity

The lowest detectable concentration of free T3 that can be distinguished from the standard 0 is 0.05 pg/ml at the 95 % confidence limit.

11.4. Interferences

- Several drugs are known to influence the binding of Triiodothyronine to the thyroid hormone carrier proteins or the metabolism to T3 and complicate the interpretation of free T3 assay results.
- Circulating autoantibodies to T3 and hormone-binding inhibitors may interfere.
- Heparin has been reported to have in vivo and in vitro effects on free T3 concentration. Therefore, do not use samples in which this anti-coagulant has been used.
- In severe nonthyroidal illness (NTI), the assessment of thyroid status becomes very difficult. TSH measurements are recommended to identify thyroid dysfunction.
- Familial dysalbuminemic conditions may yield erroneous results on direct T3 assays.
- **Not intended for newborn screening.**

11.5. Method comparison

The Free T3 ELISA was compared to another commercially available free T3 assay. 151 serum samples were analysed according in both test systems. The linear regression curve was calculated:

$$(\text{Free T3 ELISA}) = 0.923 * (\text{FT3 RIA}) + 0.350$$

$$r^2 = 0.903$$

12. LIMITATIONS OF THE PROCEDURE

"NOT INTENDED FOR NEWBORN SCREENING"

13. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use by professional persons. Not for internal or external use in humans or animals.
- When using automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested. To improve the performance of the kit on ELISA automatic systems, it is recommended to increase the number of washes.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive. Nevertheless, all materials should still be regarded and handled as potentially infectious.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.

- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Some reagents contain small amounts of CMIT/MIT as preservative. Avoid contact with skin or mucosa.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- The ELISA is only designed for qualified personnel who are familiar with good laboratory practice.

13.1. Safety note for reagents containing hazardous substances

Reagents may contain CMIT/MIT (3:1) (refer to 4.1)

Therefore, the following hazard and precautionary statements apply.

	Warning	H317	May cause an allergic skin reaction.
		P261	Avoid breathing gas/mist/vapours/spray.
		P280	Wear protective gloves/protective clothing/eye protection/face protection.
		P302+P352	IF ON SKIN: Wash with plenty of water.
		P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
		P362+P364	Take off contaminated clothing and wash it before reuse.

Further information can be found in the safety data sheet.

13.2. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

For information about the packaging materials refer to PACKAGING MATERIALS.

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
Free T3

Test Preparation

Prepare reagents and samples as described.
Establish the distribution and identification plan for all specimens and controls.
Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Standard 0	Standard 1	Standard 2	Standard 3	Standard 4	Standard 5	Sample
Standard 0	50 µl	-	-	-	-	-	-
Standard 1	-	50 µl	-	-	-	-	-
Standard 2	-	-	50 µl	-	-	-	-
Standard 3	-	-	-	50 µl	-	-	-
Standard 4	-	-	-	-	50 µl	-	-
Standard 5	-	-	-	-	-	50 µl	-
Sample	-	-	-	-	-	-	50 µl
Conjugate	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl
Cover wells with foil supplied in the kit Swirl the plate gently for 20 – 30 sec Incubate for 1 h room temperature (22...28°C) Wash each well three times with 300 µl diluted Wash Solution (if you use automated equipment, wash the wells at least 5 times).							
TMB Substrate	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl
Incubate for exactly 15 min at room temperature in the dark							
Stop Solution	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl	100 µl
Mix gently for 15 – 20 sec Photometric measurement at 450 nm, 620-630 nm							

Gold Standard Diagnostics Frankfurt GmbH

Waldstrasse 23 A6

63128 Dietzenbach, Germany

Tel.: +49 6074 23698-0

Fax: +49 6074 23698-900

E-Mail: info.frankfurt@eu.goldstandarddiagnostics.com

Website: clinical.goldstandarddiagnostics.com

Date of issue: 2024-10-11
DNOV051_FT3_IFU_rev02_fromLot_6023BN1