Ferritin

Only for in-vitro diagnostic use

Instructions for use / Notice d'utilisation / Instrucciones de uso

English	2
Français	9
Español	16
Bibliography / Bibliographie / Bibliografía	23
Abbreviations / Abréviations / Abreviaciones	26
Packaging materials / Matériels d'emballage / Materiales de embalaje	26
Symbols Key / Explication des symboles / Símbolos	27
Summary of Test Procedure / Résumé de la procédure de test / Resumen de la técnica	28

REF

DNOV100 (96 Determinations)

ENGLISH

1. INTRODUCTION

Ferritin is a protein mainly present in the cytoplasm but also found in serum and plasma as an iron-carrier protein. It contains, stores and releases iron in a controlled fashion. One molecule of ferritin is capable of binding between 4000 and 5000 atoms of iron, making ferritin the major iron storage protein for the body^{1,2}.

Iron is an essential micronutrient as it is required for an adequate erythropoietic function, oxidative mechanism and cellular immune response^{1,2,3}. Under physiologic conditions there is a balance between iron absorption, iron transport and iron storage in the human body. When this balance is interrupted, two main conditions may present: iron deficiency and iron overload^{1,2,4}. In humans, ferritin acts as a buffer against iron deficiency and overload.

Iron deficiency is the most common and widespread nutritional disorder in the world, affecting 2 billion people and in particular children, women of childbearing age with heavy menstrual flow and miscarriages, and premenopausal women. The WHO estimates that 30 % of nonpregnant and 42 % of pregnant women are affected by iron deficiency and anemia⁵⁻¹⁰.

Since the concentration of ferritin is directly proportional to the total iron stores in the body, serum ferritin concentrations are considered as a common diagnostic tool in the diagnosis of diseases affecting iron metabolism¹.

In addition to ferritin levels additional markers assessed as part of the routine diagnosis of iron deficiency are transferrin saturation and the assessment of haemoglobin levels. All these parameters should also be considered together with other parameters such as the familial history of the patients and the evaluation of clinical symptoms.

2. INTENDED USE

The Ferritin ELISA is a manual *in vitro* diagnostic device intended for the quantitative determination of ferritin in human serum or plasma from an adult population. Results are to be used in conjunction with other clinical and laboratory data as an aid in the diagnosis of diseases affecting iron metabolism.

3. PRINCIPLE OF THE ASSAY

The Ferritin assay is based on simultaneous binding of human Ferritin to two monoclonal antibodies; one is immobilized on the microplate; the other is soluble and conjugated with horseradish peroxidase (HRP). Microtiter strip wells are precoated with anti-Ferritin antibodies. Ferritin in samples and standards binds to the immobilised antibodies on the surface of the microtiter wells and the second, soluble anti-Ferritin antibody conjugated with HRP binds to the immobile antibody-Ferritin-complex during the first incubation. Afterwards a bound/free separation is performed by solid-phase washing. The immune complex is visualized by adding Tetramethylbenzidine (TMB) substrate, develops a blue colour that changes into yellow when the Stop Solution (H₂SO₄) is added. The intensity of this product is proportional to the amount of Ferritin in samples and standards. Sulphuric acid is added to stop the reaction. Absorption at 450 nm is read using an ELISA microwell plate reader. Ferritin concentration in the sample is calculated through a calibration curve.

4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 breakapart 8-well snap-off strips coated with anti-Ferritin antibody; in resealable aluminium foil.
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, < 5 % (w/v); avoid any skin contact.
- **Conjugate:** 1 bottle containing 21 mL of horseradish peroxidase (HRP) labelled anti-Ferritin antibodies; < 0.06 % (v/v) CMIT/MIT (3:1).
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (H₂O₂-TMB 0.26 g/L) (avoid any skin contact); < 0.0015 % (v/v) CMIT/MIT (3:1).
- **Wash Solution 10x conc.:** 1 bottle containing 50 mL Phosphate buffer 0.2 M, pH 7.4; < 0.06 % (v/v) CMIT/MIT (3:1).
- **Control:** 1 bottle containing 1 mL of a lot-specific, ready to use control solution. The concentration is stated on the label; < 0.06 % (v/v) CMIT/MIT (3:1).
- **Standards:** 6 bottles, 3 mL of Standard 0, 1 mL each of all other standards. The standards are ready to use and are provided at the following concentrations; < 0.06 % (v/v) CMIT/MIT (3:1).

Standard 0: 0 ng/mL

Standard 1: 5 ng/mL

Standard 2: 20 ng/mL

Standard 3: 100 ng/mL

Standard 4: 400 ng/mL

Standard 5: 800 ng/mL

For hazard and precautionary statements see 16.1

4.2. Materials supplied

- 1 Cover foil
- 1 Instructions for use (IFU)

4.3. Materials and Equipment needed

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450 nm, 620-630 nm
- Manual or automatic equipment for rinsing wells
- Precision Pipetting Devices
- Pipettes to deliver volumes between 10 and 1000 µL
- Roller shaker
- Distilled water
- Disposable tubes
- Timer

5. STABILITY AND STORAGE

The reagents are stable until the expiration date indicated on the label when stored at 2...8 °C. Once opened, the kit is stable for 5 months (not beyond the expiration date indicated on the label) if stored at 2...8 °C.

6. REAGENT PREPARATION

It is very important to bring all reagents, samples and standards to room temperature (22...28 °C) for at least 30 minutes before starting the test run! At the end of the assay, store immediately the reagents at 2...8 °C; avoid long exposure to room temperature.

6.1. Coated snap-off Strips

The ready to use break apart snap-off strips are coated with anti-Ferritin antibodies. Store at 2...8 °C. Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C; stability until expiry date. Once opened, the strips are stable for 5 months (not beyond the expiration date indicated on the label).

6.2. Conjugate

The conjugate is ready to use. Mix gently, for 5 minutes, on a roller shaker.

6.3. Standards

The standards are ready to use. After first opening standards are stable for another 5 months at 2...8° C.

6.4. TMB Substrate Solution

The bottle contains 15 mL of a tetramethylbenzidine/hydrogen peroxide system. The reagent is ready to use and has to be stored at 2...8 °C in the dark. *The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.*

6.5. Stop Solution

The bottle contains 15 mL of a ready to use sulphuric acid solution (H314).

6.6. Wash Solution

Dilute the content of each vial of the "10X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2...8°C.

In concentrated wash solution it is possible to observe the presence of crystals; in this case mix at room temperature until the complete dissolution of crystals; for greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care to transfer completely the crystals, then mix until crystals are completely dissolved.

6.7. Control

The bottle contains 1 mL of a lot-specific, ready to use control solution.

7. SPECIMEN COLLECTION AND PREPARATION

The assay should be performed using serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin or potassium EDTA) samples. Specimen can be stored at 2...8 °C for a short time (max three days). For longer storage the specimen should be aliquoted and stored deep-frozen (-20 °C). If samples are stored frozen, mix thawed samples well before testing. *Avoid repeated freezing and thawing (max 3 cycles).*

Store the sample at -20°C if the determination is not performed on the same day of the sample collection. Before using, mix gently, for 5 minutes, with a roller shaker.

7.1. Precaution

- Avoid the exposure of TMB substrate to direct sunlight, metal or oxidants. Do not freeze the solution.
- Maximum precision is required for dilution and dispensation of the reagents.

8. ASSAY PROCEDURE

8.1. Test Preparation

Please read the test instructions for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instructions for use as described. Prior to commencing the assay, the distribution and identification plan for all specimens and standards should be carefully established. Select the required number of microtiter strips or wells and insert them into the holder. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If it lasts more than ten minutes, please follow the same order during the dispensation. If more than one plate is used, it is recommended to repeat the dose response curve. Please allocate at least:

1 well	(e.g. A1)	for the substrate blank
2 wells	(e.g. B1+C1)	for standard 0
2 wells	(e.g. D1+E1)	for standard 1
2 wells	(eg. F1+G1)	for standard 2
2 wells	(eg. H1+A2)	for standard 3
2 wells	(eg. B2+C2)	for standard 4
2 wells	(eg. D2+E2)	for standard 5
2 wells	(e.g. F2+G2)	for control

It is recommended to determine standards and patient samples in duplicate.

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each standard and each patient sample.

1. Dispense 10 µL standards, control and samples into their respective wells. Add 200 µL conjugate to each well. Leave well A1 for substrate blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour at room temperature (22...28 °C)**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL diluted Wash Solution. Avoid overflows from the reaction wells.

Important note: During each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: In case you use automatic equipment, wash the wells at least 5 times.

Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.

5. Dispense 100 µL TMB Substrate Solution into all wells.
6. **Incubate for exactly 15 min at room temperature (22...28 °C) in the dark.**
7. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution. Shake the microplate gently.
Any blue colour developed during the incubation turns into yellow.
8. Measure the absorbance of the specimen at 450 nm against a reference wavelength of 620-630 nm or against blank within 5 minutes.

8.2. Measurement

Adjust the ELISA Microwell Plate Reader **to zero** using the **substrate blank in well A1**.

If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard and sample.

Where applicable calculate the **mean absorbance values** of all duplicates.

9. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A cubic spline or 4-parameter logistic (4PL) curve fit, **including Standard 0 is recommended.**

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Standard 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

In order for the assay results to be considered valid the kit standards and control must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid and samples must be retested.

10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit control provided in the kit should be tested as unknown and is intended to assist in assessing the validity of results obtained with each assay plate.

We recommend the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the Certificate of Analysis (CoA) should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories¹².

11. MEASURING RANGE

The assay measuring range (AMR) is 2.84 – 800 ng/mL.

Any value that reads below 2.84 ng/mL should be reported as “< 2.84 ng/mL”. Any value that reads above 800 ng/mL should be reported as “> 800 ng/mL”.

12. METROLOGY AND TRACEABILITY

The Ferritin ELISA has been standardized against internal reference standards (serum matrix) which have been value assigned to another commercially available test method.

13. EXPECTED VALUES

The following ranges were determined using the Ferritin ELISA and are provided for information only. The 95 % reference interval for apparently healthy adults were calculated by a non-parametric method following guidance from CLSI C28-A “Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory”.

		n°	Median (ng/mL)	Reference interval (ng/mL)
Women	Pre-menopausal	42	33.5	8.4– 156.9
	Post-menopausal	42	44.8	16.8 – 180.8
Men		84	122.5	21.7 – 276.4

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

14. PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

14.1. Detection Capability

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A, “Protocols for Determination of Limits of Detection and Limits of Quantitation” using 6 blanks and 6 low level samples.

Sensitivity	Concentration
Limit of Blank (LoB)	0.35 ng/mL
Limit of Detection (LoD)	2.00 ng/mL
Limit of Quantitation (LoQ)	2.84 ng/mL

14.2. Trueness

Trueness has been demonstrated through method comparison of the Ferritin ELISA to a commercially available LC-MS/MS using native donor samples – refer to section 14.5.

14.3. Precision

Precision of the Ferritin ELISA was determined by performing a complex precision study.

Repeatability (within laboratory)

A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators. Data from one representative lot is shown below:

Sample	n	Mean Conc. (ng/mL)	Within run (Repeatability)	
			SD	CV%
1	75	9.21	1.23	13.3%
2	75	68.30	2.47	3.6%
3	75	178.38	9.44	5.3%
4	75	380.97	26.54	7.0%
5	75	467.72	25.58	5.5%
6	75	619.17	51.33	8.3%

Reproducibility

A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators.

Data from one representative lot is shown below:

Sample	n	Mean Conc. (ng/mL)	Within Laboratory (Reproducibility)	
			SD	CV%
1	150	8.96	1.21	13.5%
2	150	67.85	4.31	6.4%
3	150	177.92	14.14	7.9%
4	150	377.41	32.10	8.5
5	150	460.17	46.83	10.2%
6	150	614.18	75.43	12.3%

14.4. Linearity

Linearity was evaluated based on CLSI EP-06, "Evaluation of the Linearity of Quantitative Measurement Procedures". For ferritin concentration by Ferritin ELISA, the measurement procedure shows linearity for the interval from 2.5 to 849 ng/mL within the allowable deviation of linearity (ADL) of $\pm 15\%$.

14.5. Method comparison

The Ferritin ELISA was compared against a commercially available quantitative automated ferritin assay, following CLSI EP-9A, "Method Comparison and Bias Estimation Using Patient Samples". A total of 100 samples, selected to represent a wide range of ferritin concentrations [15.65 to 782 ng/mL], were assayed by each method. Passing-Bablok regression analysis was performed on the comparative data:

n	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
100	0.96 [0.92 to 1.01]	0.24 [-5.46 to 3.07]	0.96

14.6. Interference

The following substances do not interfere with a bias of $> \pm 15\%$ in the Ferritin ELISA when the concentrations are below the stated threshold presented in the following table.

Potentially Interfering Substance	Threshold Concentration
Bilirubin, conjugated	15 mg/dL
Bilirubin, unconjugated	15 mg/dL
Hemoglobin	1000 mg/dL
Total Protein	10 g/dL
Triglycerides	3000 mg/dL

14.7. Serum-plasma study

The Ferritin ELISA matrix comparison study was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI EP9-A3 guidelines. A total of 20 samples covering the assay range were evaluated. Passing-Bablok regression analysis was performed on the comparative data:

Sample Type	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
SST	1.02 [0,99 to 1,25]	-1.02 [-4.98 – 0.28]	1.00
Lithium Heparin	1.04 [0.95 – 1.09]	-1.61 [-3.33 – 0.42]	1.00
Sodium Heparin	1.07 [0.98 – 1.21]	-1.92 [-4.51 – -0.22]	1.00
EDTA	0.97 [0.88 – 1.06]	-1.02 [-3.24 – 0.66]	1.00

14.8. Hook Effect

The Ferritin ELISA, a competitive enzyme immunoassay, shows no Hook Effect up to 10,000 ng/mL.

15. LIMITATIONS OF THE PROCEDURE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with in vitro immunoassays¹³. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

16. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use by professional persons. Not for internal or external use in humans or animals.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- Some reagents (standards, control, conjugate and wash solution) contain small amounts of CMIT/MIT as preservative. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- The ELISA is only designed for qualified personnel following the standards of good laboratory practice (GLP).

16.1. Safety note for reagents containing hazardous substances

Reagents may contain CMIT/MIT (3:1) (refer to 4.1)

Therefore, the following hazard and precautionary statements apply.

Warning 	H317	May cause an allergic skin reaction.
	P261	Avoid breathing gas/mist/vapours/spray.
	P280	Wear protective gloves/protective clothing/eye protection/face protection.
	P302+P352	IF ON SKIN: Wash with plenty of water.
	P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
	P362+P364	Take off contaminated clothing and wash it before reuse.

Further information can be found in the safety data sheet

16.2. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

For information about the packaging materials refer to PACKAGING MATERIALS.

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
BIBLIOGRAPHY / BIBLIOGRAPHIE / BIBLIOGRAFÍA

1. Muñoz M, García-Erce JA, Remacha ÁF. Disorders of iron metabolism. Part II:iron deficiency and iron overload. *J Clin Pathol*. 2011 Apr;64(4):287-96.
2. Moyer TP, Highsmith WE, Smyrk TC, Gross JB Jr. Hereditary hemochromatosis: laboratory evaluation. *Clin Chim Acta*. 2011 Aug 17;412(17-18):1485-92.
3. Crownover BK, Covey CJ. Hereditary hemochromatosis. *Am Fam Physician*. 2013 Feb 1;87(3):183-90.
4. Daru J, Colman K, Stanworth SJ, De La Salle B, Wood EM, Pasricha SR. Serum ferritin as an indicator of iron status: what do we need to know? *Am J Clin Nutr*. 2017 Dec;106(Suppl 6):1634S-1639S.
5. Peyrin-Biroulet L, Williet N, Cacoub P. Guidelines on the diagnosis and treatment of iron deficiency across indications: a systematic review. *Am J Clin Nutr*. 2015 Dec;102(6):1585-94.
7. De Franceschi L, Iolascon A, Taher A, Cappellini MD. Clinical management of iron deficiency anemia in adults: Systemic review on advances in diagnosis and treatment. *Eur J Intern Med*. 2017 Jul;42:16-23.
8. Cappellini MD, Musallam KM, Taher AT. Iron deficiency anaemia revisited. *J Intern Med*. 2020 Feb;287(2):153-170.
9. Friedman AJ, Shander A, Martin SR, Calabrese RK, Ashton ME, Lew I, Seid MH, Goodnough LT. Iron deficiency anemia in women: a practical guide to detection, diagnosis, and treatment. *Obstet Gynecol Surv*. 2015 May;70(5):342-53.
10. Lopez A, Cacoub P, Macdougall IC, Peyrin-Biroulet L. Iron deficiency anaemia. *Lancet*. 2016 Feb 27;387(10021):907-16.
11. Daru J, Allotey J, Peña-Rosas JP, Khan KS. Serum ferritin thresholds for the diagnosis of iron deficiency in pregnancy: a systematic review. *Transfus Med*. 2017 Jun;27(3):167-174.
12. Basic QC Practices On-line Course; <http://www.Westgard.com>.
13. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27–33

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
Ferritin

Test Preparation

Prepare reagents and samples as described.
 Establish the distribution and identification plan for all specimens and controls.
 Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate blank	Standard 0 - 5	Control	Sample
Standard 0 - 5	-	10 µL	-	-
Control	-	-	10 µL	-
Sample	-	-	-	10 µL
Conjugate	-	200 µL	200 µL	200 µL
Cover wells with foil supplied in the kit Incubate for 1 h at room temperature Wash each well three times with 300 µL diluted wash solution In case you use automatic equipment, wash the wells at least 5 times.				
TMB Substrate	100 µL	100 µL	100 µL	100 µL
Incubate for exactly 15 min at room temperature in the dark				
Stop Solution	100 µL	100 µL	100 µL	100 µL
Shake the microplate gently. Photometric measurement at 450 nm, 620-630 nm				

Gold Standard Diagnostics Frankfurt GmbH

Waldstrasse 23 A6

63128 Dietzenbach, Germany

Tel.: +49 6074 23698-0

Fax: +49 6074 23698-900

E-Mail: info.frankfurt@eu.goldstandarddiagnostics.com

Website: clinical.goldstandarddiagnostics.com