Cortisol Saliva

CE

Only for in-vitro diagnostic use

Instructions for use / Instrucciones de uso

English	2
Español	8
Bibliography / Bibliografía	14
Abbreviations / Abreviaciones	14
Packaging materials / Materiales de embalaje	14
Symbols Key / Símbolos	15
Summary of Test Procedure / Resumen de la técnica	16

REF

DSNOV20 (96 Determinations)

ENGLISH

1. INTRODUCTION

Cortisol is a steroid hormone released from the adrenal cortex in response to a hormone called ACTH (produced by the pituitary gland), it is involved in the response to stress; it increases blood pressure, blood sugar levels, may cause infertility in women, and suppresses the immune system.

Cortisol acts through specific intracellular receptors and has effects in numerous physiologic systems, including immune function, glucose-counter regulation, vascular tone, substrate utilization and bone metabolism. Cortisol is excreted primarily in urine in an unbound (free) form.

The majority of cortisol in saliva is not-bound and enters the saliva via intracellular mechanisms. Salivary cortisol levels are unaffected by salivary flow rate or salivary enzymes.

It is a high correlation between serum and saliva cortisol levels

These normal endogenous functions are the basis for the physiological consequences of chronic stress - prolonged cortisol secretion causes muscle wastage, hyperglycaemia, and suppresses immune/ inflammatory responses. The same consequences arise from long-term use of glucocorticoid drugs.

2. INTENDED USE

Competitive immunoenzymatic colorimetric method for quantitative determination of Cortisol in saliva from an adult population.

3. PRINCIPLE OF THE ASSAY

Microtiter strip wells are precoated with anti-Cortisol antibodies (solid-phase). Cortisol in the sample competes with added horseradish peroxidase labelled Cortisol (enzyme-labelled antigen) for antibody binding. After incubation a bound/free separation is performed by solid-phase washing. The immune complex formed by enzyme-labelled antigen is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product. The intensity of this product is **inversely** proportional to the amount of Cortisol in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorption at 450 nm is read using an ELISA microwell plate reader.

Cortisol concentration in the sample is calculated through a calibration curve.

4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 breakapart 8-well snap-off strips coated with anti-Cortisol; in resealable aluminium foil.
- **Stop Solution:** 1 bottle containing 15 ml sulphuric acid, 0.15 mol/l (avoid any skin contact).
- **Conjugate conc.:** 1 bottle containing 1 ml of concentrated horseradish peroxidase labelled Cortisol; contains < 0.06 % (v/v) CMIT/MIT (3:1) and 0.1 % BSA.
- **Incubation Buffer:** 1 bottle containing 30 ml phosphate buffer (50 mM, pH 7.4, BSA 1 g/l); contains < 0.06 % (v/v) CMIT/MIT (3:1)
- **TMB Substrate Solution:** 1 bottle containing 15 ml 3,3',5,5'-tetramethylbenzidine (H₂O₂-TMB 0.26 g/l) (avoid any skin contact); contains ≤ 0.0015 % (v/v) CMIT/MIT (3:1).
- **Wash solution 10x conc:** one bottle containing 50 ml (phosphate buffer 0,2M, ph 7.4; contains < 0.06 % (v/v) CMIT/MIT (3:1)
- **Control A:** 1 bottle containing 1 ml of a lot-specific control solution. The concentration is indicated on the label of the bottle; contains < 0.06 % (v/v) CMIT/MIT (3:1) and 0.1 % BSA.
- **Control B:** 1 bottle containing 1 ml of a lot-specific control solution. The concentration is indicated on the label of the bottle; contains < 0.06 % (v/v) CMIT/MIT (3:1) and 0.1 % BSA.
- **Standards:** 7 bottles, 1 ml each; contains < 0.06 % (v/v) CMIT/MIT (3:1) and 0.1 % BSA.
 - Standard 0: 0 ng/ml
 - Standard 1: 0.5 ng/ml
 - Standard 2: 1 ng/ml
 - Standard 3: 5 ng/ml
 - Standard 4: 10 ng/ml
 - Standard 5: 20 ng/ml
 - Standard 6: 100 ng/ml

For SI units: ng/ml x 2.76 = nmol/l

For hazard and precautionary statements see 14.1

4.2. Materials supplied

- 1 Cover foil
- 1 Instructions for use (IFU)

4.3. Materials and Equipment needed

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450 nm, 620-630 nm
- Manual or automatic equipment for rinsing wells
- Pipettes
- Distilled water

- Glass centrifuge tubes
- Plastic straw
- Timer

5. STABILITY AND STORAGE

The reagents are stable up to the expiry date stated on the label when stored at 2...8 °C in the dark

6. PRECAUTION

- Avoid the exposure of TMB substrate solution to direct sun light, metals or oxidants. Do not freeze the solution.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- This method allows the determination of Cortisol from 0.5 – 100 ng/ml.
- Treatment of the patient with corticosteroids, natural or synthetic steroids can impair Cortisol determination.

7. REAGENT PREPARATION

It is very important to bring all reagents, samples and standards to room temperature (22...28 °C) before starting the test run!

7.1. Microtiterplate

The ready to use break apart snap-off strips are coated with anti-Cortisol antibodies. Store at 2...8 °C. Open the bag only when it is at room temperature. *Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C; stability until expiry date.*

7.2. Conjugate

Dilute the Cortisol conjugate immediately before use. Add 10 µl conjugate to 1.0 ml Incubation Buffer. Mix gently. Stable for 3 hours at room temperature (22-28° C).

At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.

7.3. Standards

Mix gently for at least 5 minutes with rotating mixer. The standards are ready to use and have the concentration mentioned in 4.1. After first opening stable for another 6 months at 2-8°C.

7.4. TMB Substrate Solution

The bottle contains 15 ml of a tetramethylbenzidine/hydrogen peroxide system. The reagent is ready to use and has to be stored at 2...8°C in the dark. *The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.*

7.5. Stop Solution

The bottle contains 15 ml 0.15 M sulphuric acid solution . This ready to use solution has to be stored at 2...8°C.

7.6. Wash Solution

Dilute the concentrated solution with distilled water to a final volume of 500 ml prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2...8°C. In the concentrated solution it is possible to observe the presence of crystals, in this case mix at room temperature until complete dissolution of crystals. For greater accuracy dilute the whole bottle of concentrated wash solution to 500 ml and take care that all crystals are transferred by washing the bottle, then mix until crystals are completely dissolved.

7.7. Controls

Ready to use. The bottles contain 1 ml of a lot-specific control solution. The concentration is indicated on the label.

8. SPECIMEN COLLECTION AND PREPARATION

The determination of Cortisol can be performed in saliva. It is recommended to collect saliva samples with a centrifuge glass tube and a plastic straw or with the "Salivette" (Sarstedt, Ref. 51.1534.500).

Other equipment of sample collection commercially available has not been tested.

8.1. Method and Limitation

Collect saliva samples at the times indicated.

If no specific instructions have been given oral fluid (saliva) samples may be collected at any time. The following should be noted:

- a) If saliva collection is carried out in the morning ensure that this is carried out prior to brushing teeth
- b) During the day allow 1 hour after any food or drink before collecting saliva samples
- c) It is very important that a good clear sample is received – i.e. no contamination with food, lipstick, blood (bleeding gums) or other extraneous materials.

8.2. Processing Instructions

1. Let the saliva flow down through the straw into the centrifuge glass tube
2. Centrifuge at 3000 rpm for 15 min.
3. Store at -20°C for at least 1 hour
4. Defrost samples
5. Centrifuge again for 15 min. at 3000 rpm
6. The saliva sample is now ready to be tested.
7. Store the sample at 2...8 °C for one week or at -20°C for longer time.

8.3. Saliva Processing Instructions with Salivette Sarstedt

1. Remove the swab from the suspended insert of the Salivette.
2. Gently chewing the swab for 1 minute produces a sufficient quantity of Saliva.
3. Replace the swab into the Salivette and firmly close the tube using the stopper.
4. Centrifuge the Salivette for 2 minutes at 1000g (rcf) for saliva generation.
5. Remove the insert complete with the swab from the centrifuge vessel and discard. The clear Saliva is now
6. ready for analysis (at least 1 ml of saliva should be recovered with this method)

9. ASSAY PROCEDURE

9.1. Test Preparation

Please read the test instructions for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instructions for use as described. Prior to commencing the assay, the distribution and identification plan for all specimens and standards should be carefully established. Select the required number of microtiter strips or wells and insert them into the holder. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than one plate is used, it is recommended to repeat the dose response curve. Please allocate at least:

1 well	(e.g. A1)	for the substrate blank
2 wells	(e.g. B1+C1)	for standard 0
2 wells	(e.g. D1+E1)	for standard 1
2 wells	(e.g. F1+G1)	for standard 2
2 wells	(e.g. H1+A2)	for standard 3
2 wells	(e.g. B2+C2)	for standard 4
2 wells	(e.g. D2+E2)	for standard 5
2 wells	(e.g. F2+G2)	for standard 6
2 wells	(e.g. H2+A3)	for control A
2 wells	(e.g. B3+C3)	for control B

It is recommended to determine standards, controls and patient samples in duplicate.

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each standard, each control and each patient sample.

1. Dispense 25 µl standards, controls and samples into their respective wells.
2. Add 200 µl diluted Cortisol-HRP Conjugate to each well. Leave well A1 for substrate blank.
3. Cover wells with the foil supplied in the kit.
4. **Incubate for 1 hour at 37 °C (± 0.5 °C).**
5. When incubation has been completed, remove the foil, aspirate the content of the wells. Wash the wells 3 times with 300 µl of diluted wash solution. Avoid overflows from the reaction wells. During each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: In case you use an automatic washer, it is advised to do 6 washing steps.

Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.

6. Dispense 100 µl TMB Substrate Solution into all wells.
7. **Incubate for exactly 15 min at room temperature (22...28°C) in the dark.**
8. Dispense 100 µl Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution. Shake the microplate gently.
Any blue colour developed during the incubation turns into yellow.
9. Measure the absorbance of the specimen at 450 nm against reference wavelength of 620-630 nm or against Blank within 5 minutes.

9.2. Measurement

Adjust the ELISA Microwell Plate Reader **to zero** using the **substrate blank in well A1**.

If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard and sample.

Where applicable calculate the **mean absorbance values** of all duplicates.

10. RESULTS

10.1. Calculation of results

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) curve fit, **including Standard 0 is required**. Other curve fitting algorithms are not recommended.

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Standard 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

In order for the assay results to be considered valid the kit standards and controls must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid and samples must be retested.

Conversion of units

To convert results to SI units:

$$\text{nmol/L} = \text{ng/mL} \times 2.76$$

To convert results to mass units:

$$\text{ng/mL} = \text{nmol/L} \times 0.36$$

10.2. Reference Value

The following values can be used as preliminary guideline until each laboratory established its own normal ranges.

A.M.	3.0 – 10.0 ng/ml
P.M.	0.6 – 2.5 ng/ml

Please pay attention to the fact that the determination of a range of expected values for a “normal” population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

11. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit control provided in the kit should be tested as unknown and are intended to assist in assessing the validity of results obtained with each assay plate.

We recommend the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the Certificate of Analysis (CoA) should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories⁸.

12. SPECIFIC PERFORMANCE CHARACTERISTICS

12.1. Precision

Intra Assay Variation

Within run variation was determined by replicate determination (20x) of three different saliva control in one assay. The within assay variability is $\leq 10\%$.

Inter Assay Variation

Between run variation was determined by replicate measurements (10x) of three different saliva control in different lots. The between assay variability is $\leq 8.3\%$.

12.2. Specificity

The cross reaction of the antibody calculated at 50% according to Abraham is:

Cortisol	100 %
Prednisolone	46.2 %
11-Deoxycortisol	4 %
Cortisone	3.69 %
Prednisone	3.10 %
11 α OH Progesterone	1 %
Progesterone	< 0.1 %
Aldosterone	< 0.1 %
Pregnenolone	< 0.1 %
17 beta Estradiol	< 0.1 %
Estrone 3-sulfate	< 0.1 %
Estriol	< 0.1 %
Testosterone	< 0.1 %
Spirolactone	< 0.1 %
DHEA	< 0.1 %
DHEA-S	< 0.1 %
Androstenedione	< 0.1 %
Androsterone	< 0.1 %
DHT	< 0.1 %
Danazol	< 0.1 %
Cholesterol	< 0.1 %
Dexamethasone	< 0.1 %

12.3. Sensitivity

The lowest detectable concentration of Cortisol that can be distinguished from standard 0 is 0.12 ng/ml at the 95% confidence limit.

12.4. Accuracy

The recovery of 6.25 – 12.5 – 25 – 50 ng/ml of Cortisol added to a saliva sample gave an average value (\pm SD) of 95.42% \pm 9.11% with reference to the original concentrations.

12.5. Correlation

The Cortisol Saliva ELISA kit was compared to a previous version of the Cortisol saliva ELISA kit. 35 saliva samples were analysed.

The linear regression curve was calculated:

$$Y = 0.94 \cdot X - 0.11$$

$$r^2 = 0.799$$

13. LIMITATIONS OF THE PROCEDURE

Sample(s), which are contaminated microbiologically, should not be used in the assay. Highly lipemic or haemolysed specimen(s) should similarly not be used. It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than one plate is used, it is recommended to repeat the dose response curve. Addition of the substrate solution initiates a kinetic reaction, which is terminated by the addition of the stop solution. Therefore, the addition of the substrate and the stopping solution should be added in the same sequence to eliminate any time deviation during reaction. Plate readers measure vertically. Do not touch the bottom of the wells. Failure to remove adhering solution adequately in the aspiration or decantation wash step(s) may result in poor replication and spurious results.

14. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the test kits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with in vitro immunoassays⁹. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive. Nevertheless, all materials should still be regarded and handled as potentially infectious.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.

- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Some reagents (standards, controls, incubation buffer, conjugate and wash solution) contain small amounts of CMIT/MIT as preservative. Avoid contact with skin or mucosa.
- The ELISA is only designed for qualified personnel who are familiar with good laboratory practice.

14.1. Safety note for reagents containing hazardous substances

Reagents may contain CMIT/MIT (3:1) (refer to 4.1)

Therefore, the following hazard and precautionary statements apply.

	Warning	H317	May cause an allergic skin reaction.
		P261	Avoid breathing gas/mist/vapours/spray.
		P280	Wear protective gloves/protective clothing/eye protection/face protection.
		P302+P352	IF ON SKIN: Wash with plenty of water.
		P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
		P362+P364	Take off contaminated clothing and wash it before reuse.

Further information can be found in the safety data sheet.

14.2. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

For information about the packaging materials refer to PACKAGING MATERIALS.

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Cortisol Saliva

Test Preparation

Prepare reagents and samples as described.
Establish the distribution and identification plan for all specimens and controls
Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate blank	Standard 0 - 6	Controls	Sample
Standard 0 - 6	-	25 µl	-	-
Controls	-	-	25 µl	-
Sample	-	-	-	25 µl
Conjugate	-	200 µl	200 µl	200 µl
Cover wells with foil supplied in the kit Incubate for 1 h at 37°C Wash each well three times with 300 µl diluted wash solution In case you use an automatic washer, it is advised to do 6 washing steps.				
TMB Substrate	100 µL	100 µL	100 µL	100 µL
Incubate for exactly 15 min at room temperature (22...28°C) in the dark				
Stop Solution	100 µL	100 µL	100 µL	100 µL
Shake the microplate gently Photometric measurement at 450 nm against a reference wavelength of 620-630 nm or against blank within 5 minutes.				

Gold Standard Diagnostics Frankfurt GmbH

Waldstrasse 23 A6

63128 Dietzenbach, Germany

Tel.: +49 6074 23698-0

Fax: +49 6074 23698-900

E-Mail: info.frankfurt@eu.goldstandarddiagnostics.com

Website: clinical.goldstandarddiagnostics.com

Date of issue: 2024-09-13
DSNOV20_CORT Saliva_IFU_rev02_fromLot_6162BN1