# Androstenedione Saliva

CE

Enzyme immunoassay for the quantitative determination of Androstendione in human saliva

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### 1. INTRODUCTION

---

Androstenedione (also known as  $\Delta^4$ -androstenedione) is a steroid hormone produced in the adrenal glands and the gonads as an intermediate step in the biochemical pathway that produces the androgen testosterone and the estrogens estrone and estradiol. It is the common precursor of male and female sex hormones. Some androstenedione is also secreted into the plasma, and may be converted in peripheral tissues to testosterone and estrogens.

Androstenedione has relatively weak androgenic activity, estimated at ~ 20% of testosterone. Secretion and production rates also exceed those of testosterone in women in whom significant extra-adrenal conversion of androstenedione to testosterone occurs. In premenopausal women the adrenal glands and ovaries each produce about half of the total androstenedione (about 3 mg/day). After menopause androstenedione production is about halved, primarily due to the reduction of steroid secreted by the ovary. Nevertheless, androstenedione is the principal steroid produced by the postmenopausal ovary.

The high serum-saliva correlation for androstenedione suggests that individual differences in serum androstenedione levels may be accurately estimated using saliva as a non-invasive alternative specimen.

### 2. INTENDED USE

---

Competitive immunoenzymatic colorimetric method for quantitative determination of Androstenedione in human saliva.

### 3. PRINCIPLE OF THE ASSAY

---

Microtiterplates are precoated with anti-Androstenedione antibodies (solid phase). Androstenedione in the sample competes with added horseradish peroxidase labelled Androstenedione (enzyme-labelled antigen) for antibody binding. After incubation a bound/free separation is performed by solid-phase washing. The immune complex formed by enzyme-labelled antigen is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product. The intensity of this product is **inversely** proportional to the amount of Androstenedione in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorption at 450 nm is read using an ELISA Microtiterplate reader. Androstenedione concentration in the sample is calculated through a calibration curve.

### 4. MATERIALS

---

#### 4.1. Reagents supplied

- **Microtiterplate:** 12 breakapart 8-well snap-off strips coated with anti-Androstenedione antibodies; in resealable aluminium foil.
  - **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, < 5 % (w/v); (avoid any skin contact).
  - **Conjugate conc.:** 1 bottle containing 1 mL concentrated horseradish peroxidase labelled Androstenedione.
  - **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine ( $H_2O_2$ -TMB 0.26 g/L) (avoid any skin contact).
  - **Wash solution 50x conc:** 1 bottle containing 20 mL (NaCl 45 g/L; Tween-20 55 g/L)
  - **Incubation buffer:** 1 bottle containing 30 mL phosphate buffer pH 7.5, BSA 1 g/L, stabilizer.
  - **Control A:** 1 bottle containing 1 mL of a lot-specific control solution. The concentration is indicated on the label.
  - **Control B:** 1 bottle containing 1 mL of a lot-specific control solution. The concentration is indicated on the label.
  - **Standards:** 5 bottles, 1 mL each
    - Standard 0: 0 pg/mL
    - Standard 1: 20 pg/mL
    - Standard 2: 100 pg/mL
    - Standard 3: 400 pg/mL
    - Standard 4: 1000 pg/mL
- For SI units: pg/mL x 3.487= pmol/L

#### 4.2. Materials supplied

- 1 Cover foil
- 1 Instruction for use (IFU)
- 1 Distribution and identification plan

#### 4.3. Materials and Equipment needed

- 37 °C Incubator
- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450 nm, 620 – 630 nm
- Centrifuge (3000 rpm) with rotor for glass tubes
- Manual or automatic equipment for rinsing Microtiterplates
- Pipettes to deliver volumes between 10 and 1000  $\mu$ L
- Roller shaker
- Distilled water
- Plastic straw
- Glass tubes for centrifugation
- Timer

#### 4.4. Precautions

- Some reagents contain Proclin 300® as preservative. Avoid contact with skin or mucosa.
- The TMB substrate contains an irritant, which may be harmful if inhaled, ingested, or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- This method allows the determination of Androstenedione from 5 – 1000 pg/mL.
- Treatment of the patient with cortisone, natural or synthetic steroids can impair Androstenedione determination.
- Avoid the exposure of TMB substrate solution to direct sunlight, metals, or oxidants. Do not freeze the solution.

#### 5. STABILITY AND STORAGE

---

The originally closed reagents are stable up to the expiry date stated on the label when stored at 2...8 °C in the dark.

#### 6. REAGENT PREPARATION

---

*It is very important to bring all reagents, samples, standards, and controls to room temperature (22...28 °C) before starting the test run! At the end of the assay, store immediately the reagents at 2...8 °C; avoid long exposure to room temperature.*

##### 6.1. Microtiterplate

The ready to use break apart snap-off strips are coated with anti-Androstenedione antibodies. Store at 2...8 °C. Open the bag only when it is at room temperature. *Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C; stability until expiry date.*

##### 6.2. Conjugate

Prepare immediately before use. Add 10 µL of conjugate to 1 mL of incubation buffer. Mix gently. The diluted conjugate is stable for 3 h at 22...28° C.

##### 6.3. Standards

Before use mix 5 min with a rotating mixer. The standards are ready to use and have the following concentration of Androstenedione:

|             |            |
|-------------|------------|
| Standard 0: | 0 pg/mL    |
| Standard 1: | 20 pg/mL   |
| Standard 2: | 100 pg/mL  |
| Standard 3: | 400 pg/mL  |
| Standard 4: | 1000 pg/mL |

For samples with Androstenedione concentration greater than 1000 pg/mL dilute the sample (1:2) with standard 0.

Once opened the standards are stable for another 6 months at 2...8 °C.

##### 6.4. TMB Substrate Solution

The bottle contains 15 mL of a tetramethylbenzidine/hydrogen peroxide system. The reagent is ready to use and has to be stored at 2...8 °C in the dark. *The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.*

##### 6.5. Stop Solution

The bottle contains 15 mL sulphuric acid solution (H314). This ready to use solution has to be stored at 2...8 °C.

##### 6.6. Wash Solution

Dilute the contents of each vial of the buffered wash solution concentrate (50x) with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio. The diluted wash solution is stable for 30 days at 2...8 °C.

##### 6.7. Controls

Ready to use. The bottles contain 1 mL of a lot-specific control solution. The concentration is indicated on the label.

#### 7. SPECIMEN COLLECTION AND PREPARATION

---

The determination of Androstenedione can be performed in saliva. It is recommended to collect saliva samples with a centrifuge glass tube and a plastic straw. Other commercially available sample collector devices have not been tested.

##### 7.1. Method and Limitations

Collect saliva samples at the times indicated.

If no specific instructions have been given, oral fluid (saliva) samples may be collected at any time. The following should be noted:

- a) If saliva collection is carried out in the morning, ensure that this is carried out prior to brushing teeth
- b) During the day allow 1 hour after any food or drink, oral intake of pharmaceutical drugs or tooth cleaning before collecting saliva samples.
- c) It is very important that a good clear sample is received – i.e., no contamination with food, lipstick, blood, chewing gum or other extraneous materials.

## 7.2. Saliva Processing Instructions

1. Let the saliva flow down through the straw into the centrifuge glass tube
2. Centrifuge at 3000 rpm for 15 min.
3. Store at -20 °C for at least 1 hour
4. Defrost samples
5. Centrifuge again for 15 min. at 3000 rpm
6. The saliva sample is now ready to be tested.
7. Store the sample at 2...8 °C for one week or at -20 °C for longer time.

## 8. ASSAY PROCEDURE

---

### 8.1. Test Preparation

Please read the test protocol carefully **before** performing the assay. Result reliability depends on strict adherence to the test protocol as described. Prior to commencing the assay, the distribution and identification plan for all specimens and standards should be carefully established on the result sheet supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than one plate is used, it is recommended to repeat the dose response curve. Please allocate at least:

- |         |              |                         |
|---------|--------------|-------------------------|
| 1 well  | (e.g. A1)    | for the substrate blank |
| 2 wells | (e.g. B1+C1) | for standard 0          |
| 2 wells | (e.g. D1+E1) | for standard 1          |
| 2 wells | (e.g. F1+G1) | for standard 2          |
| 2 wells | (e.g. H1+A2) | for standard 3          |
| 2 wells | (e.g. B2+C2) | for standard 4          |
| 2 wells | (eg. D2+E2)  | for control A           |
| 2 wells | (eg. F2+G2)  | for control B           |

*It is necessary to determine standards, controls and patient samples in duplicate.*

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each standard and each patient sample.

1. Dispense 50 µL standards, controls and samples into their respective wells.
2. Add 150 µL diluted Conjugate to each well. Leave well A1 for substrate blank.
3. Cover wells with the foil supplied in the kit.
4. **Incubate for 1 hour at 37 °C.**
5. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well 3 times with 300 µL of diluted wash solution. Avoid overflows from the reaction wells.

**Important note:** During each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

**Automatic washer:** In case you use automatic equipment, wash the wells at least 5 times.

*Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.*

6. Dispense 100 µL TMB Substrate Solution into all wells.
7. **Incubate for exactly 15 min at room temperature (22...28 °C) in the dark.**
8. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution. Shake the Microtiterplate gently.  
*Any blue colour developed during the incubation turns into yellow.*
9. Measure the absorbance of the specimen at 450 nm against a reference wavelength of 620-630 nm or against blank within 5 minutes.

### 8.2. Measurement

Adjust the ELISA Microtiterplate Reader **to zero** using the **substrate blank in well A1**.

*If - due to technical reasons - the ELISA Microtiterplate reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!*

**Measure the absorbance** of all wells at **450 nm** against a reference wavelength of 620-630 nm or against the blank within 5 min. Record the absorbance values for each standard and patient sample in the distribution and identification plan.

Where applicable calculate the **mean absorbance values** of all duplicates.

## 9. RESULTS

---

### 9.1. Calculation of Results

Calculate the mean absorbance for each point of the standard curve and each sample. Plot the mean value of absorbance of the standards against concentration. Draw the best-fit curve through the plotted points. (e.g., Four Parameter Logistic).

Interpolate the values of the samples on the standard curve to obtain the corresponding values of the concentrations expressed in pg/mL.

## 9.2. Reference values

As the values of salivary Androstenedione have a circadian pattern, we suggest collecting the samples at the same hour (8 a.m.). Each laboratory must establish its own normal ranges based on patient population.

The following value should be considered as a guideline:

|       |                |                 |
|-------|----------------|-----------------|
| Women | normal         | 20 – 160 pg/mL  |
|       | P.C.O.-Hirsute | 120 – 300 pg/mL |
| Men   |                | 20 – 150 pg/mL  |

Please pay attention to the fact that the determination of a range of expected values for a “normal” population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore, each laboratory should consider the range given by the manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

## 10. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of Androstenedione for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the standard curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

## 11. SPECIFIC PERFORMANCE CHARACTERISTICS

### 11.1. Precision

#### Intra Assay Variation

Within run variation was determined by replicate determination (16x) of two different saliva controls in one assay. The within assay variability is  $\leq 8.5\%$ .

#### Inter Assay Variation

Between run variation was determined by replicate measurements (10x) of two different saliva control with different lots of kit. The inter assay variability is  $\leq 11\%$ .

### 11.2. Cross Reactivity

The cross reaction of the antibody calculated at 50% according to Abraham is:

|                                |        |
|--------------------------------|--------|
| Androstenedione                | 100%   |
| Testosterone                   | 1.2%   |
| Epitetosterone                 | 0.2%   |
| 5 $\alpha$ Dihydrotestosterone | 0.1%   |
| Dehydroepiandrosterone (DHEA)  | 0.1%   |
| Progesterone                   | 0.001% |
| Estrone                        | 0.001% |
| Cortisol                       | 0.001% |

### 11.3. Analytic Sensitivity

The lowest detectable concentration of Androstenedione that can be distinguished from the standard 0 is 5 pg/mL at the 95 % confidence limit.

### 11.4. Accuracy

The recovery of 50 – 200 – 500 pg/mL of Androstenedione added to sample gave an average value ( $\pm$ SD) of  $102.60\% \pm 13.23\%$  with reference to the original concentrations.

### 11.5. Correlation with serum sample

The NovaTec Androstenedione Saliva ELISA kit was compared to another commercially available Androstenedione saliva assay. 38 saliva samples were analysed according in both test systems. The linear regression curve was calculated:

$$y = 0.46x + 5.51$$

$$r^2 = 0.983$$

y = Androstenedione Saliva NovaTec ELISA kit

x = reference Salivary Androstenedione ELISA kit

## 12. PRECAUTIONS AND WARNINGS

---

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the test kits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use by professional persons. Not for internal or external use in humans or animals.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive. Nevertheless, all materials should still be regarded and handled as potentially infectious.
- Do not interchange reagents or strips of different production lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- Sample(s), which are contaminated microbiologically, should not be used in the assay. Highly lipemic or haemolysed specimen(s) should similarly not be used.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If it lasts longer than ten minutes, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate.
- Addition of the substrate solution initiates a kinetic reaction, which is terminated by the addition of the stop solution. Therefore, the addition of the substrate and the stopping solution should be added in the same sequence to eliminate any time deviation during reaction.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- The NovaTec ELISA is only designed for qualified personnel who are familiar with good laboratory practice.

### 12.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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## Androstenedione Saliva

### Test Preparation

Prepare reagents and samples as described.  
Establish the distribution and identification plan for all specimens and controls on the result sheet supplied in the kit.  
Select the required number of microtiter strips or wells and insert them into the holder.

### Assay Procedure

|                                                                                                                                                                                                                                  | Substrate blank | Standard 0 - 4 | Controls | Sample |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|----------------|----------|--------|
| Standard 0 - 4                                                                                                                                                                                                                   | -               | 50 µL          | -        | -      |
| Controls                                                                                                                                                                                                                         | -               | -              | 50 µ L   | -      |
| Sample                                                                                                                                                                                                                           | -               | -              | -        | 50 µL  |
| Conjugate                                                                                                                                                                                                                        | -               | 150 µL         | 150 µL   | 150 µL |
| Cover wells with foil supplied in the kit<br><b>Incubate for 1 h at 37°C</b><br>Wash each well three times with 300 µl diluted wash solution<br><b>In case you use an automatic washer, it is advised to do 5 washing steps.</b> |                 |                |          |        |
| TMB Substrate                                                                                                                                                                                                                    | 100 µL          | 100 µL         | 100 µL   | 100 µL |
| <b>Incubate for exactly 15 min at room temperature (22...28°C) in the dark</b>                                                                                                                                                   |                 |                |          |        |
| Stop Solution                                                                                                                                                                                                                    | 100 µL          | 100 µL         | 100 µL   | 100 µL |
| Shake the microplate gently<br>Photometric measurement at 450 nm against a reference wavelength of 620-630 nm or against blank within 5 minutes.                                                                                 |                 |                |          |        |

### NovaTec Immundiagnostica GmbH, part of Gold Standard Diagnostics Europe

Waldstraße 23 A6

63128 Dietzenbach, Germany

Tel.: +49 6074 23698-0

Fax: +49 6074 23698-900

E-Mail: [info@goldstandarddiagnostics.eu](mailto:info@goldstandarddiagnostics.eu)

Website: [www.NovaTec-ID.com](http://www.NovaTec-ID.com)

DSNOV27-Androstenedione S-en,it,es\_20221206\_gr.\_Lot\_5094A