Serazym[®] Astrovirus

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Serazym[®] Astrovirus

Enzyme immunoassay for the qualitative detection of the capsid polyprotein of human-pathogenic astroviruses in stool samples of human origin

REF	E-045		96
REF	E-045-A2		2x96
IVD	In-vitro-diagnostic medical device		CE

Seramun Diagnostica GmbH • Sprenhagener Str. 1 • 15754 Heidesee • Germany •
 T +49 33767 791-10 • info@seramun.com • www.seramun.com

 In-vitro diagnostic medical device	 Unique device identifier	 Manufacturer
 Country of manufacture and date of manufacture	 Article number	 Serial number
 Keep away from sunlight	 Humidity limitation	 Batch code
 Consult instructions for use	 Temperature limit	 Do not reuse
 Sufficient for <i>n</i> tests	 Biohazard	 Use-by date
		 Attention

Intended Use

Serazym® Astrovirus is an IVD test for the qualitative determination of the capsid polyprotein of human-pathogenic astroviruses in stool samples of human origin by a laboratory professional user.

It is intended to aid in the diagnosis of astrovirus-associated gastroenteritis in samples of patients with symptoms of a gastroenteritis.

The test must not be used with specimen materials such as biopsy specimens, culture suspensions and materials other than stool samples of human origin, for diagnosis, monitoring, prognosis, prediction, screening, as companion diagnostic, in the near-patient setting and by lay persons.

Principle of the Test

Serazym® Astrovirus is an enzyme immunoassay based on polyclonal and monoclonal antibodies against astrovirus-specific antigens. Diluted, untreated stool samples as well as negative and positive control are dispensed simultaneously with peroxidase (HRP)-labeled monoclonal anti-astrovirus antibodies, into the wells of the microtiter plate coated with polyclonal anti-astrovirus antibodies. After incubation, unbound components are removed by a washing step, then HRP converts the colorless substrate solution to a blue reaction product in the following enzymatic reaction step. After incubation the reaction is stopped by addition of the stop solution, resulting in a color change from blue to yellow. The optical density (OD) of the reaction product measured at 450 nm measurement and ≥ 620 nm reference filter, respectively, is directly proportional to the concentration of specifically bound astrovirus antigens.

Test Components (Delivery Scope)

			For 96 wells	For 2x 96 wells
1	WELLS	Microtiter plate coated with polyclonal anti-astrovirus antibodies (sheep)	12 single breakable 8-well strips, light blue color marking, vacuum- sealed with desiccant	2x 12 single breakable 8-well strips, light blue color marking, vacuum-sealed with desiccant
2	WASHBUF (10x)	Wash buffer (10x) Seramun® Wash buffer A TRIS-based buffer	100 mL concentrate for 1000 mL solution, colorless, white cap	2x 100 mL concentrate for 2 x 1000 mL solution, colorless, white cap
3	DIL	Sample diluent Seramun® Sample diluent A Phosphate-based buffer	100 mL, ready to use, colored yellow black cap	2x 100 mL, ready to use, colored yellow black cap
4	CONTROL +	Positive control Astrovirus positive sample	1.5 mL, ready to use, colored blue, red cap	3.0 mL, ready to use, colored blue, red cap
5	CONTROL -	Negative control TRIS-based buffer	1.5 mL, ready to use, colored blue, green cap	3.0 mL, ready to use, colored blue, green cap
6	CONJ HRP	HRP conjugate HRP-labeled, monoclonal anti-astrovirus antibodies (mouse)	12 mL, ready to use, colored green, brown cap	24 mL, ready to use, colored green, brown cap

7	SUBSTR	Substrate SeramunBlau® automat fast < 0.1 % 3,3',5,5'- tetramethylbenzidine; < 0.05 % hydrogen peroxide	15 mL, ready to use, colorless, blue cap	28 mL, ready to use, colorless, blue cap
8	STOP	Stop solution SeramunBlau® stop 0.25 M sulphuric acid	15 mL, ready to use, colorless, yellow cap	28 mL, ready to use, colorless, yellow cap
9	COVER	Covering film	1 piece	-
10		Certificate of Analysis	1 piece	1 piece
11		Instructions for Use	1 piece	1 piece

Additional Materials and Aids Required for the Test Procedure

Adjustable single-channel micropipette • 8-channel micropipette or multichannel micropipette with pipette tips • 8-channel wash comb with vacuum pump and waste bottle or microplate washer • microplate reader with 450 nm measuring filter and ≥ 620 nm reference filter • deionized water • measuring cylinder • tubes for sample preparation

Important Information

This device is for *in-vitro* diagnostic use only. Follow the instructions carefully. The kit may be performed by health professionals only.

Do not use reagents from damaged packages or bottles. The shelf life specified must be observed. Do not mix components with reagents from other manufacturers.

Mixing of test kit components of different lots is only allowed for wash buffer (10x), sample diluent, substrate and stop solution.

Wash buffer (10x), sample diluent, substrate and stop solution are universally applicable for Serazym® stool ELISA Adenovirus (E-017), Astrovirus (E-045), Campylobacter (E-093), Clostridium difficile GDH (E-107), Clostridium difficile Toxin A+B (E-040), Cryptosporidium parvum (E-039), Entamoeba histolytica (E-018), Giardia (E-106), H. pylori 2nd Gen. (E-114), Norovirus (E-061) and Rotavirus (E-020).

All serious incidents occurring in relation with Serazym® Astrovirus must be reported to the manufacturer and the competent authority of the EU member state in which user and/or patient are located.

Information on Assay Procedure

All reagents should be stored at 2...8 °C. Bring all test components to room temperature before use. Reagents that appear contaminated should not be used.

Each well of a microtiter plate can only be used once. Each sample and control have to be pipetted with a new pipette tip. Positive and negative control are ready to use.

For larger sample series, pipetting reagents from liquid reservoirs using a multi-channel micropipette is recommended to avoid time delays and contaminations. Follow the pipetting scheme and time schedules of the protocol.

The aspiration and washing steps can be performed manually or with the help of a microplate washer or waterjet pump. Wash solution should be allowed a minimum reaction time of 5 s in the wells per wash cycle. Remove wash buffer residues by thoroughly aspirating or tapping out the wells!

Protect substrate from light!

Safety Instructions

Reagents must not be swallowed. Contact with skin or mucous membranes should be avoided. Handle all components and patient samples as if potentially hazardous and infectious. Additional information may be taken from the Safety Data Sheet.

Product contains the following hazardous component/-s:

Test component	Hazard labeling and supplementary information on ingredients
WELLS	Contains material of animal origin.
WASHBUF (10x)	EUH208: Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. EUH210: Safety data sheet available on request. Preservatives: < 0.0015 % reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1); < 0.1 % 5-bromo-5-nitro-1,3-dioxane
DIL	Contains material of animal origin. Preservatives: < 0.1 % sodium azide
CONTROL +	Contains material of microbial and animal origin. Preservatives: < 0.1 % sodium azide
CONTROL -	Contains material of animal origin. Preservatives: < 0.01 % sodium azide
CONJ HRP	EUH210: Safety data sheet available on request. Contains material of animal origin. Preservative: < 0.1 % 5-bromo-5-nitro-1,3-dioxane
SUBSTR	Hazard component: 2-pyrrolidone Signal word: Danger  H319: Causes serious eye irritation. H360: May damage fertility or the unborn child. P201: Obtain special instructions before use. P280: Wear protective gloves/protective clothing/eye protection/face protection. P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P308+P313: IF exposed or concerned: Get medical advice/attention. Restricted to professional users. Preservatives: < 0.00015 % reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1)
STOP	Hazard component: sulphuric acid 2.5 % Signal word: Warning  H290: May be corrosive to metals.

Limitations of the Procedure

The qualitative enzyme immunological detection of astrovirus antigens in stool samples does not allow a correlation between the measured OD and the severity of an infection. Also, it is not allowed to correlate absorbances of samples with the absorbance of the positive control. Cross contamination of reagents and samples may result in false results. Incorrect dilutions, insufficiently homogenized samples, and particles not sedimented by centrifugation may cause false positive as well as false negative test results. Sample taking within the acute phase of an infection is recommended as the number of excreted particles in this phase is expected to reach its maximum. A negative test result obtained with Serazym® Astrovirus does not exclude an astrovirus infection. False negative tests may result from improper timing of sample collection or inhomogeneous antigen distribution in the sample. The overall interpretation of the ELISA test result should consider the full clinical picture.

Sample Treatment

Sample Collection

Collect stool sample in suitable sampling container.

Sample Shelf Life and Storage

Stool samples should be stored at 2...8 °C immediately after collection and examined within 72 h or stored frozen at -20 °C. Repeated (> 3x) freezing and thawing of samples should be avoided due to the risk of incorrect results. Stool samples that have already been diluted in Seramun® Sample diluent A according to the instructions for use can be stored at 2...8 °C for up to 48 h and subsequently analyzed by ELISA.

Sample Preparation

Mix untreated stool samples well and dilute 1 : 6 with sample buffer.

Example: Pipette 500 µL sample buffer into a reaction tube. For solid or semi-solid stool samples transfer 100 mg (approx. 2 – 3 mm diameter) with a disposable stick, for liquid stool samples transfer 100 µL into the sample buffer and mix thoroughly. If necessary, sediment suspended particles by centrifugation in a microcentrifuge for 1 min at maximum speed.

Reagent Treatment

Reagent Shelf Life and Storage

The complete test kit with sealed reagent bottles and microtitration strips can be stored at 2...8 °C until the printed expiration date. All opened test kit components are stable for up to 2 months when stored properly at 2...8 °C. The diluted wash buffer can be stored at 2...8 °C for up to 1 month.

Reagent Preparation

Microtiter plate with breakable 8-well strips is vacuum sealed with desiccant. Allow packaging to reach room temperature before opening. Protect unused wells from moisture and store refrigerated with desiccant in the original bag carefully resealed. Dilute wash buffer (10x) 1 : 10 with deionized water.

Example: 10 mL wash buffer (10x) + 90 mL deionized water.

Assay Procedure

1. Allow test reagents and required number of wells to reach room temperature (RT). Shake reagents gently before use. Avoid foaming.
2. Pipette 2 drops (or 75 μL) **CONJ HRP** HRP conjugate per well.
3. Add 75 μL **CONTROL +** Positive control
75 μL **CONTROL -** Negative control
50 μL diluted stool specimen each and mix gently.
4. Cover the plate and incubate for 60 min at RT.
5. Decant, then wash each well 5x with 300 μL diluted wash buffer. Tap dry onto absorbent paper if necessary.
6. Add 2 drops (or 75 μL) **SUBSTR** substrate per well.
7. Incubate for 10 min at RT **protected from light**.
8. Add 2 drops (or 75 μL) **STOP** stop solution per well, mix gently.
9. Read OD at 450 nm measuring filter and ≥ 620 nm reference filter with a microplate reader within 30 min following reaction stop.

Evaluation of Results

Qualitative Evaluation

Cut-off determination: OD negative control + 0.10

Samples showing OD values equal with or higher than the cut-off are considered positive, samples with OD values below cut-off are considered negative for astrovirus antigens.

The test run is valid if:

- mean OD of the negative control is ≤ 0.15 (manual processing)
 ≤ 0.30 (automatic processing)
- mean OD of the positive control is ≥ 1.00

If the above-mentioned quality criteria are not met, test should be repeated strictly following the test procedure (incubation times and temperatures, sample and wash buffer dilution, wash steps, etc.). In case of repeated failure of the quality criteria contact the manufacturer.

Interpretation of Results

Positive	\geq cut-off
Negative	$<$ cut-off

It is recommended that each laboratory establishes its own normal and pathological reference ranges.

Performance Characteristics

Precision

To determine precision, 4 stool samples were measured multiple times. For the determination of the intra-assay coefficient of variation (CV) samples were measured in an 8-fold determination in one test run.

The determination of the inter-assay coefficient of variation was done by an 8-fold determination in 6 different test runs.

Sample	Intra-assay coefficient of variation		Inter-assay coefficient of variation	
	\bar{x} OD	CV (%)	\bar{x} OD	CV (%)
1	1.667	8.9	1.853	3.8
2	0.994	6.4	1.019	5.8
3	0.443	6.1	0.583	11.9
4	0.185	9.8	0.350	9.7

Detection Limit

The lower detection limit of Serazym® Astrovirus was determined by titration of purified astrovirus antigen at 6 ng/mL.

Sensitivity and Specificity

Sensitivity and specificity of Serazym® Astrovirus have been determined in a retrospective study with 98 stool specimens in comparison to a commercially available ELISA.

n = 98	ELISA positive	ELISA negative
Serazym® ELISA positive	49	0
Serazym® ELISA negative	2	47

Sensitivity: 96.1 %

Specificity: 100 %

Cross Reactivity

Stool samples positive for one of the following pathogens have been tested with the Serazym® Astrovirus and showed no cross reactivity:

Adenovirus, *Campylobacter coli*, *Campylobacter jejuni*, *Clostridioides difficile*, *Giardia lamblia*, norovirus, rotavirus, *Salmonella enteritidis*.

Negative stool suspensions were spiked with the following microorganisms with a bacterial count of $\geq 10^8$ colony-forming units per mL in sample buffer and tested negative in the Serazym® Astrovirus (450 nm measurement and ≥ 620 nm reference filter < cut-off):

<i>Aeromonas hydrophila</i>	(ATCC 7966)	<i>Klebsiella pneumoniae</i>	(ATCC 13883)
<i>Bacillus cereus</i>	(ATCC 11778)	<i>Peptostreptococcus anaerobius</i>	(ATCC 27337)
<i>Bacillus subtilis</i>	(ATCC 6633)	<i>Proteus vulgaris</i>	(ATCC 8427)
<i>Bacteroides fragilis</i>	(ATCC 25285)	<i>Pseudomonas aeruginosa</i>	(ATCC 10145)
<i>Candida albicans</i>	(ATCC 10231)	<i>Salmonella enterica</i> Serovar <i>enteritidis</i>	(ATCC 13076)
<i>Campylobacter coli</i>	(ATCC 33559)	<i>Salmonella enterica</i> Serovar <i>typhimurium</i>	(ATCC 14028)
<i>Campylobacter jejuni</i>	(ATCC 33291)	<i>Shigella flexneri</i>	(ATCC 12022)
<i>Citrobacter freundii</i>	(ATCC 8090)	<i>Shigella sonnei</i>	(ATCC 25931)
<i>Clostridium sordellii</i>	(ATCC 9714)	<i>Staphylococcus aureus</i>	(ATCC 25923)
<i>Enterobacter aerogenes</i>	(ATCC 13048)	<i>Staphylococcus epidermidis</i>	(ATCC 12228)
<i>Enterobacter cloacae</i>	(ATCC 13047)	<i>Vibrio parahaemolyticus</i>	(ATCC 17802)
<i>Enterococcus faecalis</i>	(ATCC 29212)	<i>Vibrio cholerae</i>	Clinical isolate
<i>Escherichia coli</i>	(ATCC 25922)	<i>Yersinia enterocolitica</i> Serotyp O3, O9	Clinical isolate

Interference

None of the following substances in the indicated concentrations added to astrovirus positive and negative stool samples did show a significant impact on the test result:

Barium sulfate (5 %), Buscopan® (2 mg/mL), Cyclamate (5 %), Diclofenac (2 mg/mL), human haemoglobin (5 mg/mL), Hylak® N (5 %)*, Iberogast® (5 %), Immodium® akut duo (0.2/12.5 mg/mL), Loperamide (0.2 mg/mL), Metronidazole (2 mg/mL), Mucin (5 mg/mL), Nexium® (2 mg/mL), palmitic acid (20 %), Pentofuryl® (2 mg/mL), Pepto-Bismol (1 mg/mL), Perenterol (2.5 mg/mL), Rennie® (8 mg/mL), Simagei® (2 mg/mL), stearic acid (20 %), Vancomycin (2 mg/mL).

*The addition of 5 % (v/v) Hylak® N to astrovirus positive stool suspensions may decrease OD values.


Application

Automatic Processing

Performing Serazym® Astrovirus on fully automated microplate processors (e.g., DS2®, DSX®; Dynex Technologies) may cause elevated absorbance values in comparison to the manual procedure caused by differences in the wash procedures and technical specifications of the equipment. In these cases, a maximum value of OD = 0.3 is permissible for the negative control. It is recommended to program a wash protocol with 10 s soak time per strip and wash step. A final wash step with deionized water and a soak time of 10 s is recommended after each wash cycle. If necessary, the number of wash steps may be increased to 7x or 8x.

Correlation: manual – automatic processing

A panel of 96 stool specimens was processed manually and automatically in parallel (DS2®, Dynex Technologies). The correlation coefficient was calculated at $r = 0.993$.

Change History

Version	Section	Modifications
2022-08_v01	Entire document	Updating of Intended Use Conversion of subsections Insertion of Safety Instructions Updating of sample dilution Removal of Pipetting Scheme
2023-09_v02	Entire document	Update of Important Information; Safety Instructions Editorial changes

References

1. Cukor, G and Blacklow, NR (1984): "Human Viral Gastroenteritis". *Microbiological Reviews* 48(2): 157-179.
2. Gaggero, A, O’Ryan, M et al. (1998): "Prevalence of Astrovirus Infection among Chilean Children with Acute Gastroenteritis". *Journal of Clinical Microbiology* 36(12): 3691-3693.
3. Koukou, G, Niendorf, S, Hornei, B et al (2019): "Human astrovirus infection associated with encephalitis in an immunocompetent child: a case report". *J Med Case Reports* 13(341). <https://doi.org/10.1186/s13256-019-2302-6>
4. Palombo, EA and Bishop, RF (1996): "Annual Incidence, Serotype Distribution and Genetic Diversity of Human Astrovirus Isolates from Hospitalized Children in Melbourne, Australia". *Journal of Clinical Microbiology* 34(7): 1750-1753.
5. Rohwedder, A (2000): "Virale Gastroenteritiden, Erreger und Diagnostik". *Mikrobiologie* 10: 121-126.