RF-Absorbent

To be used with Infectious disease IgM ELISAs.

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IBL@IBL-International.com
www.IBL-International.com

Pretreatment with RF-Absorbent

	In order to avoid interferences of specific IgG and rheumatoid factors, patient sera should be treated with RF-Absorbent (REF KIRF561). Do not treat standards and controls with RF-Absorbent!
---	--

For sample pretreatment with RF-Absorbent exist two variants:

A) Dilution of RF-Absorbent in Buffer

The volumes stated below are for one run with 4 strips (32 determinations).

Preparation of Components

Dilution	Relation	Remarks	Storage	Stability
1 mL RF-Absorbent + 20 mL Diluent Buffer	1:21	Incubate ≥ 1 min (< 20 min) at RT (18-25°C)	2-8°C	8 w

Sample dilution

Sample	to be diluted with	Relation	Remarks
Serum / Plasma	Diluent Buffer (+ RF-Absorbent)	1:101	z.B. 5 μ L + 500 μ L DILBUF

Samples with RF-Absorbent: Do not incubate > 20 min to avoid adsorption of specific antibodies. Pretreated samples may be turbid.

B) Direct addition of RF-Absorbent to the diluted sample

Sample dilution

Sample	to be diluted with	Relation	Remarks
Serum / Plasma	Diluent Buffer	1:101	z.B. 5 μ L + 500 μ L DILBUF

Addition of RF-Absorbent

Sample treatment	Remarks
Pipette 20 μ L of RF-Absorbent to 400 μ L of 1:101 diluted sample.	Mix well. Incubate ≥ 1 min (< 20 min) at RT (18-25°C).

Samples with RF-Absorbent: Do not incubate > 20 min to avoid adsorption of specific antibodies. Pretreated samples may be turbid.

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COMPLAINTS: Complaints may be submitted initially written or vocal. Subsequently they need to be filed including the test performance and results in writing in case of analytical reasons.

WARRANTY: The product is warranted to be free from material defects within the specific shelf life and to comply with product specifications delivered with the product. The product must be used according to the Intended use, all instructions given in the instructions for use and within the product specific shelf life. Any modification of the test procedure or exchange or mixing of components of different lots could negatively affect the results. These cases invalidate any claim for replacement.

LIMITATION OF LIABILITY: IN ALL CIRCUMSTANCES THE EXTENT OF MANUFACTURER'S LIABILITY IS LIMITED TO THE PURCHASE PRICE OF THE KIT(S) IN QUESTION. IN NO EVENT SHALL MANUFACTURER BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING DAMAGES FOR LOST PROFITS, LOST SALES, INJURY TO PERSON OR PROPERTY OR ANY OTHER INCIDENTAL OR CONSEQUENTIAL LOSS.

	IBL International GmbH	Tel.: + 49 (0) 40 532891 -0 Fax: -11
	Flughafenstr. 52A, 22335 Hamburg, Germany	E-MAIL: IBL@IBL-International.com
		WEB: http://www.IBL-International.com