DHEA-S

Enzyme immunoassay for the quantitative determination of DHEA-S in human serum or plasma

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ENGLISH

1. INTRODUCTION

Dehydroepiandrosterone sulfate (DHEA-S), is a natural steroid hormone found atop of the kidneys in the human body. DHEA-S derived from enzymatic conversion of DHEA in adrenal and extradrenal tissues. DHEA-S is also produced in the gonads, adipose tissue and the brain. It is the most abundant hormone in the human body, and it is precursor of all sex steroids

As the most DHEA-S is produced by the zona reticularis of the adrenal, it is argued that there is a role in the immune and stress response. DHEA-S may have more biological roles. Its production in the brain suggests that it also has a role as neurosteroid.

Measurement of serum DHEA-S is a useful marker of adrenal androgen synthesis. Abnormally low levels have been reported in hypoadrenalism, while elevated levels occur in several conditions, e.g., virilising adrenal adenoma and carcinoma, 21-hydroxylase and 3 β -hydroxysteroid dehydrogenase deficiencies and in some cases of female hirsutism. Women with polycystic ovary syndrome tend to have normal or mildly elevated levels of DHEA-S. As very little DHEA-S is produced by the gonads, measurement of DHEA-S levels may aid in the localization of androgen source in virilising conditions.

DHEA-S shows no diurnal variation.

2. INTENDED USE

Competitive immunoenzymatic colorimetric method for quantitative determination of DHEA-S in human serum or plasma (EDTA).

3. PRINCIPLE OF THE ASSAY

Microtiter strip wells are precoated with anti-DHEA-S antibodies (solid-phase). DHEA-S in the sample competes with added horseradish peroxidase labelled DHEA-S (enzyme-labelled antigen) for antibody binding. After incubation a bound/free separation is performed by solid-phase washing. The immune complex formed by enzyme-labelled antigen is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product. The intensity of this product is **inversely** proportional to the amount of DHEA-S in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorption at 450 nm is read using an ELISA microwell plate reader.

4. MATERIALS

4.1. Reagents supplied

- **Anti-DHEA-S Coated Wells:** 12 breakapart 8-well snap-off strips coated with anti-DHEA-S antibody; in resealable aluminium foil.
- **Stop Solution:** 1 bottle containing 15 ml sulphuric acid, 0.15 mol/l (avoid any skin contact).
- **DHEA-S-HRP Conjugate:** 1 bottle containing 12 ml horseradish peroxidase labelled DHEA-S.
- **TMB Substrate Solution:** 1 bottle containing 15 ml 3, 3', 5, 5'-tetramethylbenzidine (H₂O₂-TMB 0.26 g/l) (avoid any skin contact).
- **Diluent:** 1 bottle containing 100 ml HEPES 187 mM pH 7.5; BSA 0.5 g/l.
- **Control A:** 1 bottle containing 1 ml of a lot-specific control solution. The concentration is indicated on the label of the bottle.
- **Control B:** 1 bottle containing 1 ml of a lot-specific control solution. The concentration is indicated on the label of the bottle.
- **Wash solution 10 x conc.:** 1 bottle containing 50 ml of a 10 x concentrated solution of phosphate buffer 0.2 M, pH 7.4.
- **DHEA-S Standards:** 6 bottles, 1 ml each
 - Standard 0: 0 μ g/ml
 - Standard 1: 0.1 μ g /ml
 - Standard 2: 0.4 μ g /ml
 - Standard 3: 1.0 μ g /ml
 - Standard 4: 4.0 μ g /ml
 - Standard 5: 10.0 μ g /ml

4.2. Materials supplied

- 1 Cover foil
- 1 Instructions for use (IFU)

4.3. Materials and Equipment needed

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450 nm, 620- 630 nm
- Manual or automatic equipment for rinsing wells
- Pipettes to deliver volumes between 10 and 1000 μ l
- Vortex tube mixer
- Distilled water
- Disposable tubes
- Timer

5. STABILITY AND STORAGE

The reagents are stable up to the expiry date stated on the label when stored at +2...+8° C in the dark.

6. REAGENT PREPARATION

It is very important to bring all reagents, samples and standards to room temperature (+22...+28° C) for at least 30 minutes before starting the test run! At the end of the assay, store immediately the reagents at 2-8° C; avoid long exposure to room temperature.

6.1. Coated snap-off Strips

The ready to use break apart snap-off strips are coated with anti-DHEA-S IgG antibodies. Store at +2...+8° C. Open the bag only when it is at room temperature. *Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at +2...+8° C; stability until expiry date.*

6.2. DHEA-S-HRP Conjugate

The DHEA-S-HRP conjugate is ready to use.

6.3. Diluent

The diluent is ready to use.

6.4. Wash Solution

Dilute the concentrated solution with distilled water to a final volume of 500 ml prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted solution is stable for 30 days at 2...8° C. In the concentrated solution it is possible to observe the presence of crystals, in this case mix at room temperature until complete dissolution of crystals. For greater accuracy dilute the whole bottle of concentrated wash solution to 500 ml and take care that all crystals are transferred by washing the bottle, then mix until crystals are completely dissolved.

6.5. DHEA-S Standards

The standards are ready to use and are stable until the expiry date printed on the label. Opened standards are stable for six months at +2...+8° C.

6.6. TMB Substrate Solution

The bottle contains 15 ml of a tetramethylbenzidine/hydrogen peroxide system. *The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.*

6.7. Stop Solution

The bottle contains 15 ml 0.15 M sulphuric acid solution.

6.8. Controls

The bottles contain 1 ml of a lot-specific control solution. The concentration is indicated on the label. The controls are ready to use.

7. SAMPLE COLLECTION AND PREPARATION

The determination of DHEA-S can be performed in human plasma as well as in serum of patients. Store the sample at -20° C if the determination is not performed on the same day of the sample collection. *Avoid repeated freezing and thawing.*

7.1. Sample Dilution

Immediately before use dilute each sample 1:50 with the diluent (i.e. add 20 µl of sample to 980 µl of diluent). Mix well.

7.2. Precaution

- Avoid the exposure of TMB substrate to direct sunlight, metal or oxidants. Do not freeze the solution.
- Maximum precision is required for dispensation of the reagents.
- This method allows the determination of DHEA-S from 0.1 to 10 µg/ml
- Treatment of the patient with cortisone, natural or synthetic steroids can impair DHEA-S determination.

8. ASSAY PROCEDURE

8.1. Test Preparation

Please read the instructions for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instructions for use as described. Prior to commencing the assay, the distribution and identification plan for all specimens and standards should be carefully established. Select the required number of microtiter strips or wells and insert them into the holder. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than one plate is used, it is recommended to repeat the dose response curve. Please allocate at least:

1 well	(e.g. A1)	for the substrate blank
2 wells	(e.g. B1+C1)	for standard 0
2 wells	(e.g. D1+E1)	for standard 1
2 wells	(e.g. F1+G1)	for standard 2
2 wells	(e.g. H1+A2)	for standard 3
2 wells	(e.g. B2+C2)	for standard 4
2 wells	(e.g. D2+E2)	for standard 5
2 wells	(e.g. F2+G2)	for control A
2 wells	(e.g. H2+A3)	for control B

It is recommended to determine standards, controls and patient samples in duplicate.

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each standard and each patient sample.

1. Dispense 30 µl standards, controls and diluted samples into their respective wells.
2. Add 100 µl DHEA-S-HRP Conjugate to each well. Leave well A1 for substrate blank.
3. Cover wells with the foil supplied in the kit.
4. **Incubate for 1 hour at 37° C.**
5. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µl of diluted wash solution. Avoid overflows from the reaction wells.

Important note: During each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: In case you use automatic equipment, wash the wells at least 5 times.

Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.

6. Dispense 100 µl TMB Substrate Solution into all wells.
7. **Incubate for exactly 15 min at room temperature (+22 to +28 °C) in the dark.**
8. Dispense 100 µl Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution. Shake the microplate gently.
Any blue colour developed during the incubation turns into yellow.
9. Measure the absorbance of the specimen at 450 nm against a reference wavelength of 620-630 nm or against blank within 5 minutes.

8.2. Measurement

Adjust the ELISA Microwell Plate Reader **to zero** using the **substrate blank in well A1**.

If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard and patient sample.

9. RESULTS

9.1. Calculation of Results

Calculate the mean of the absorbencies for each point of the standard curve and of each sample.

Plot the values of absorbance of the standards against concentration.

Draw the best-fit curve through the plotted points (esp. Four Parameter Logistics). Interpolate the values of the samples on the standard curve to obtain the corresponding values of the concentration expressed in µg/ml.

9.2. Reference Values

The serum and plasma DHEA-S reference values are:

	Woman [µg/ml]	Man [µg/ml]
Newborns	0.9 – 1.8	0.9 – 1.8
Before puberty	0.25 – 1.0	0.25 – 1.0
Adults	0.9 – 3.6	0.9 – 3.6
After menopause	< 0.25 – 1.0	--
Pregnancy	0.25 – 1.8	--

Please pay attention to the fact that the determination of a range of expected values for a “normal” population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore, each laboratory should consider the range given by the manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

10. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of DHEA-S for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the standard curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

11. SPECIFIC PERFORMANCE CHARACTERISTICS

11.1. Precision

Intra Assay Variation

Within run variation was determined by replicate determination (16x) of three different serum samples in one assay. The within assay variability is $\leq 7.9\%$.

Inter Assay Variation

Between run variations was determined by replicate measurements (20x) of three different control sera in different lots. The between assay variability is $\leq 10.4\%$.

11.2. Specificity

The cross reaction of the antibody calculated at 50% according to Abraham:

Cross-reagent	Cross reactivity (%)
DHEA-S	100 %
DHEA	100 %
Androstenedione	59 %
Androsterone	16 %
DHT	1.0 %
Testosterone	0.63 %
Estrone	0.3 %
Progesterone	0.27 %
Pregnenolone	0.18 %
17-OH Progesterone	0.15 %
11-deoxycortisol	0.08 %
Cortisone	0.013 %
Cortisol	< 0.01 %
17 beta- Estradiol	< 0.01 %
Corticosterone	< 0.01 %
17a-Estradiol	< 0.01 %
Cholesterol	< 0.001 %
Estriol	< 0.001 %
Estradiol Sulphate	< 0.001 %
Aldosterone	< 0.001 %
Estradiol-3-Sulphate-17-Glucuronide	< 0.001 %

11.3. Specificity: interfering substances

Interference by Bilirubin, Hemoglobin and Triglycerides has been investigated on the DHEA-S ELISA kit:

Substance	Assayed Conc.	Interference
Bilirubin	0.2 mg/ml	No
Hemoglobin	2 mg/ml	No
Triglycerides	6 mg/ml	No

No interference has been observed with the substances under investigation; following good laboratory practices, it is anyway advised to avoid using highly lipemic or haemolysed samples.

11.4. Specificity: plasma and SST tube

Interference in plasma and SST (serum separation tube) samples has been evaluated. Serum obtained from the same patient has been used as reference.

Sample	Interference
SST (serum separation tube)	No
EDTA plasma	No
Lithium heparin plasma	No
Sodium-heparin plasma	No

No interference has been observed.

11.5. Sensitivity

The lowest detectable concentration of DHEA-S that can be distinguished from standard 0 is 0.04 µg/ml at the 95% confidence limit.

11.6. Accuracy

The recovery of 0.6 - 1.25 - 2.5 and 5 µg/ml of DHEA-S added to sample gave an average value (±SD) of 102.87% ± 8.63 % with reference to the original concentrations.

The dilution test performed on three sera diluted 2 - 4 and 8 times gave an average value (±SD) of 100.15% ± 9.02 %.

11.7. Correlation with RIA

The DHEA-S ELISA was compared to another commercially available DHEA-S assay. 42 serum samples have been analysed.

The linear regression curve was calculated:

$$y = 0.93x + 0.28$$

$$r^2 = 0.961$$

12. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use by professional persons. Not for internal or external use in Humans or Animals.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive. Nevertheless, all materials should still be regarded and handled as potentially infectious.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- Microbiologically contaminated samples should not be used in the assay. Highly lipemic or haemolysed specimens should also not be used.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Some reagents contain small amounts of Proclin 300® as preservative. Avoid contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.

12.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

For information about the packaging materials refer to PACKAGING MATERIALS.

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
DHEA-S

Test Preparation

Prepare reagents and samples as described.
 Establish the distribution and identification plan for all specimens and controls.
 Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate blank	Standard 0 - 5	Controls	Sample
Standard 0 - 5	-	30 µl	-	-
Controls	-	-	30 µl	-
Sample	-	-	-	30 µl
Conjugate	-	100 µl	100 µl	100 µl
Cover wells with foil supplied in the kit Incubate for 1 hour at 37° C Wash each well three times with 300 µl of diluted wash solution In case you use automatic equipment, wash the wells at least 5 times.				
TMB Substrate	100 µl	100 µl	100 µl	100 µl
Incubate for exactly 15 min at room temperature (22...28° C) in the dark				
Stop Solution	100 µl	100 µl	100 µl	100 µl
Shake the microplate gently Photometric measurement at 450 nm, 620-630 nm				

Gold Standard Diagnostics Frankfurt GmbH

Waldstrasse 23 A6

63128 Dietzenbach, Germany

Tel.: +49 6074 23698-0

Fax: +49 6074 23698-900

E-Mail: info.frankfurt@eu.goldstandarddiagnostics.com

Website: clinical.goldstandarddiagnostics.com

DNOV005_DHEA-S_IFU_rev01_fromLot_5990BN