Free Testosterone

Only for in-vitro diagnostic use

Instructions for use / Notice d'utilisation / Istruzioni per l'uso/ Instrucciones de uso

| | |
|---|----|
| English | 2 |
| Français | 9 |
| Italiano | 17 |
| Español | 25 |
| Bibliography/Bibliographie/Bibliografia/Bibliografía | 34 |
| Packaging materials / Matériels d'emballage / Materiali d'imballaggio / Materiales de embalaje | 34 |
| Symbols Key / Explication des symboles / Legenda / Símbolos | 35 |
| Summary of Test Procedure / Résumé de la procédure de test / Schema della procedura / Resumen de la técnica | 36 |

| |
|-----|
| REF |
|-----|

DNOV009 (96 Determinations)

ENGLISH

1. INTRODUCTION

Testosterone is found in circulation predominantly linked to carrier proteins, the most common of which being sex-hormone binding globulin (SHBG). Testosterone plays a key role in the development of primary and secondary sexual characteristics in males and is involved in the production of female sexual hormones.

Only 1-2 % of testosterone in circulation is not bound to any protein and is biologically active – this is referred to as 'free testosterone' (FT). Bioavailable testosterone refers to the sum of FT and the testosterone bound to serum albumin, since it is bound with low affinity and readily able to dissociate to become available for its biological function.

In males elevated levels of testosterone are associated with several conditions such as, early (precocious) puberty, congenital adrenal hyperplasia (CAH), androgen insensitivity syndrome (AIS), steroid use and testicular or adrenal tumours. Whereas the major causes of suppressed levels include Klinefelter's syndrome, testicular damage, pituitary disorders etc. In females of all ages, elevated testosterone levels can be associated with a variety of virilising conditions including adrenal tumours and polycystic ovarian syndrome (PCOS).

These clinical conditions are associated with either a lack or excess of testosterone in circulation (hypoandrogenism or hyperandrogenism). Diagnosis of these disorders involve the quantification of total testosterone (TT) in association with other clinical evidence and laboratory data. However, clinical manifestations of androgen disorders are often associated with normal levels of TT. In such cases, additional information may be gained by the assessment of the biologically active, FT level. Several androgen disorders can be caused by alteration of SHBG production which affects the levels of FT available in serum.

Measurement of FT can be considered useful in the diagnosis of several conditions including androgen deficiency in men and androgen excess in women¹. Assessment of free testosterone levels may prove beneficial² and may avoid an incorrect diagnosis of hypogonadism in cases when low concentrations of total testosterone are determined and alterations of SHBG are suspected.

There is an observed and well documented circadian variation of testosterone levels in men with circulating concentrations being higher in the morning and declining throughout the day³. Testosterone levels also decline in ageing males (andropause) and is associated with loss of muscle and bone mass, leading to osteoporosis, loss of libido, erectile dysfunction, depression and impaired cognitive function⁴.

2. INTENDED USE

Free Testosterone ELISA is a manual *in vitro* diagnostic device intended for the quantitative determination of free testosterone in human serum or plasma. Results are to be used in conjunction with other clinical and laboratory data as an aid in the diagnosis and monitoring of disorders involving the male sex hormones (androgens).

3. PRINCIPLE OF THE ASSAY

The Free Testosterone ELISA is a competitive enzyme immunometric assay (ELISA) where free testosterone (antigen) in the sample competes with the antigenic testosterone conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti-testosterone coated on the Microtiterplate (solid phase).

After the incubation, the bound/free separation is performed by a simple solid phase washing. Then, the enzyme HRP in the bound fraction reacts with the Substrate (H_2O_2) and the TMB Substrate and develops a blue colour that changes into yellow when the Stop Solution (H_2SO_4) is added. The colour intensity is inversely proportional to the free testosterone concentration of in the sample.

Free testosterone concentration in the sample is calculated through a calibration curve.

4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 breakapart 8-well snap-off strips coated with anti-Testosterone, in resealable aluminium foil.
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.15 M (avoid any skin contact).
- **Conjugate:** 1 bottle containing 15 mL of horseradish peroxidase labelled Testosterone.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (H_2O_2 -TMB 0.26 g/L), (avoid any skin contact).
- **Wash solution 10x conc.:** 1 bottle containing 50 mL of a 10x concentrated solution of phosphate buffer 0.2 M, pH 7.4.
- **Control A:** 1 bottle containing 1 mL of a lot-specific, ready to use control solution. The concentration is stated on the label.
- **Control B:** 1 bottle containing 1 mL of a lot-specific, ready to use control solution. The concentration is stated on the label.

- **Standards:** 6 bottles, 1 mL each
 - Standard 0: 0.0 pg/mL
 - Standard 1: 0.2 pg/mL
 - Standard 2: 1.0 pg/mL
 - Standard 3: 4.0 pg/mL
 - Standard 4: 20.0 pg/mL
 - Standard 5: 100.0 pg/mL

4.2. Materials supplied

- 1 Cover foil
- 1 Instructions for use (IFU)

4.3. Materials and Equipment needed

- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450 nm, 620-630 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing Microtiterplate
- Pipettes
- Roller shaker
- Distilled water
- Disposable tubes
- Timer

5. STABILITY AND STORAGE

The reagents are stable up to the expiry date stated on the label when stored at 2...8 °C in the dark. Once opened, the kit is stable at 2...8 °C for 6 months.

6. REAGENT PREPARATION

It is very important to bring all reagents, samples and standards to room temperature (22...28 °C) for at least 30 minutes before starting the test run! At the end of the assay store the reagents immediately at 2...8 °C; avoid long exposure to room temperature.

6.1. Microtiterplate

The ready to use break apart snap-off strips are coated with anti-Testosterone antibodies. Store at 2...8 °C. Open the bag only when it is at room temperature. *Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C; stability until expiry date.*

6.2. Conjugate

Testosterone-HRP Conjugate is a ready to use solution.

6.3. Standards

Before use mix 5 min. with a rotating mixer. The standards are ready to use. After first use the standards are still stable for another 6 months if stored at 2...8 °C.

6.4. Controls

Before use mix 5 min. with a rotating mixer. The bottles contain 1 mL of a lot-specific control solution, ready to use. The concentration is indicated on the label.

6.5. TMB Substrate Solution

The bottle contains 15 mL of a tetramethylbenzidine/hydrogen peroxide system. The reagent is ready to use and has to be stored at 2...8 °C in the dark. *The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.*

6.6. Stop Solution

The bottle contains 15 mL 0.15 M sulphuric acid solution. This ready to use solution has to be stored at 2...8 °C.

6.7. Wash Solution

Dilute the content of the concentrated Wash Solution with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2...8 °C. In concentrated wash solution it is possible to observe the presence of crystals, in this case mix at room temperature until complete dissolution of crystals. For greater accuracy dilute the whole bottle of concentrated wash solution to 500 mL on taking care also to transfer crystals with washing of the bottle, then mix until crystals are completely dissolved.

7. SPECIMEN COLLECTION AND PREPARATION

The assay should be performed using serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin or potassium EDTA) samples.

| Sample Storage | Duration |
|-----------------------|-----------------|
| 2...8 °C | 24 hours |
| Freeze/thaw cycles | 1 cycle |

8. ASSAY PROCEDURE

8.1. Test Preparation

Please read the test instructions for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instructions for use as described. Prior to commencing the assay, the distribution and identification plan for all specimens and standards should be carefully established. Select the required number of microtiter strips or wells and insert them into the holder. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If it lasts more than ten minutes, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve. Please allocate at least:

| | | |
|---------|--------------|-------------------------|
| 1 well | (e.g. A1) | for the substrate blank |
| 2 wells | (e.g. B1+C1) | for standard 0 |
| 2 wells | (e.g. D1+E1) | for standard 1 |
| 2 wells | (e.g. F1+G1) | for standard 2 |
| 2 wells | (e.g. H1+A2) | for standard 3 |
| 2 wells | (e.g. B2+C2) | for standard 4 |
| 2 wells | (e.g. D2+E2) | for standard 5 |
| 2 wells | (e.g. F2+G2) | for control A |
| 2 wells | (e.g. H2+A3) | for control B |

It is recommended to determine standards, controls and patient samples in duplicate.

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each standard and each patient sample.

Adjust the incubator to 37 °C.

1. Dispense 20 µL standards, controls and samples into their respective wells. Add 100 µL conjugate to each well. Leave well A1 for substrate blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour at 37 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL diluted wash solution. Avoid overflows from the reaction wells. During each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: In case you use automatic equipment, wash the wells at least 5 times.

Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.

5. Dispense 100 µL TMB Substrate Solution into all wells.
6. **Incubate for exactly 15 min at room temperature (22...28 °C) in the dark.**
7. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution. Shake the Microtiterplate gently.
Any blue colour developed during the incubation turns into yellow.
8. Measure the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against blank within 5 minutes.

8.2. Measurement

Adjust the ELISA Microtiterplate Reader to **zero** using the **substrate blank in well A1**.

If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** against a reference wavelength of 620-630 nm or against the blank within 5 min. Record the absorbance values for each standard and patient sample.

Where applicable calculate the **mean absorbance values** of all duplicates.

9. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) curve fit, **including Standard 0 is required.**

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Standard 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

Conversion of units

To convert results to SI units:

$$\text{pmol/L} = \text{pg/mL} \times 3.47$$

To convert results to mass units:

$$\text{pg/mL} = \text{pmol/L} \times 0.289$$

10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit controls provided in the kit should be tested as unknowns and are intended to assist in assessing the validity of results obtained with each assay plate.

We recommend the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the Certificate of Analysis (CoA) should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories⁵.

11. MEASURING RANGE

The assay measuring range (AMR) is 0.34 - 70.0 pg/mL (1.18 - 242.9 pmol/L).

Any value that reads below 0.34 pg/mL (1.18 pmol/L) should be reported as "< 0.34 pg/mL (< 1.18 pmol/L)". Any value that reads above 70.0 pg/mL (242.9 pmol/L) should be reported as "> 70.0 pg/mL (242.9 pmol/L)".

12. METROLOGY AND TRACEABILITY

The Free Testosterone ELISA has been standardised against internal reference standards (serum matrix) which have been value assigned to another commercially available test method.

13. REFERENCE VALUES

The following ranges were determined using the Free Testosterone ELISA and are provided for information only.

| | No. of subjects | Median (pg/mL) | Reference Interval (pg/mL) |
|-----------------|-----------------|----------------|----------------------------|
| Males | | | |
| 21 – 49 years | 120 | 14.13 | 5.01 – 27.78 |
| > 50 years | 120 | 12.75 | 4.11 – 21.85 |
| Females | | | |
| pre-menopausal | 120 | 0.55 | < LOQ – 1.70 |
| post-menopausal | 120 | 0.75 | < LOQ – 2.34 |

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

14. SPECIFIC PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

14.1. Analytical Sensitivity

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A, "Protocols for Determination of Limits of Detection and Limits of Quantitation" using 6 blanks and 6 low level samples.

| Sensitivity | Concentration |
|-----------------------------|---------------|
| Limit of Blank (LoB) | 0.10 pg/mL |
| Limit of Detection (LoD) | 0.20 pg/mL |
| Limit of Quantitation (LoQ) | 0.34 pg/mL |

14.2. Trueness

Trueness has been demonstrated through method comparison of the Free Testosterone ELISA to a commercially available assay using native donor samples – refer to section 14.5.

14.3. Precision

Precision of the Free Testosterone ELISA was determined by performing a complex precision study.

Repeatability: A total of 6 serum samples were assayed in 5 replicates, once per day for 5 days by 3 operators. Data from one representative lot is shown below:

| Sample | n | Mean Conc. (pg/mL) | Within run (Repeatability) | |
|--------|----|--------------------|----------------------------|------|
| | | | SD | CV% |
| 1 | 75 | 0.68 | 0.08 | 11.3 |
| 2 | 75 | 1.62 | 0.10 | 5.9 |
| 3 | 75 | 5.25 | 0.19 | 3.7 |
| 4 | 75 | 10.43 | 0.48 | 4.6 |
| 5 | 75 | 34.99 | 1.17 | 3.3 |
| 6 | 75 | 68.65 | 4.62 | 6.7 |

Reproducibility: A total of 6 serum samples were assayed in 5 replicates, once per day for 5 days by 3 operators.

Results for the combined data from two lots is shown below:

| Sample | n | Mean Conc. (pg/mL) | Within Laboratory (Reproducibility) | |
|--------|-----|--------------------|-------------------------------------|------|
| | | | SD | CV% |
| 1 | 150 | 0.75 | 0.14 | 18.0 |
| 2 | 150 | 1.74 | 0.24 | 14.1 |
| 3 | 150 | 5.48 | 0.52 | 9.4 |
| 4 | 150 | 10.72 | 0.94 | 8.7 |
| 5 | 150 | 37.41 | 4.83 | 12.9 |
| 6 | 150 | 73.81 | 10.83 | 14.7 |

14.4. Linearity

Linearity was evaluated based on CLSI EP-06, "Evaluation of the Linearity of Quantitative Measurement Procedures". For free testosterone concentration by Free Testosterone ELISA, the measurement procedure shows linearity for the interval from 0.27 to 89.83 pg/mL within the allowable deviation of linearity (ADL) of $\pm 15\%$.

14.5. Method comparison

The Free Testosterone ELISA was compared against a commercially available quantitative manual ELISA test, following CLSI EP-9A, "Method Comparison and Bias Estimation Using Patient Samples". A total of 95 samples, selected to represent a wide range of free testosterone concentrations, was assayed by each method. Passing-Bablok regression analysis was performed on the comparative data:

| n | Slope [95% CI] | Intercept (pg/mL) [95% CI] | Correlation coefficient (r) |
|----|------------------------|-------------------------------|--------------------------------|
| 95 | 0.94 [0.88 to 0.99] | -0.02 [-0.14 to 0.06] | 0.93 |

14.6. Analytical Specificity

The specificity was assessed with the following cross-reactants.

| Cross-reactant | Concentration tested (unit) | Mean %Cross reactivity |
|------------------------------|-----------------------------|------------------------|
| 11-keto-testosterone | 100 ng/mL | 0.1% |
| 11-β-hydroxy-testosterone | 10 ng/mL | 0.2% |
| 17αOH-Progesterone | 500 ng/mL | 0.0% |
| Aldosterone | 3000 ng/mL | 0.0% |
| Androstenedione | 100 ng/mL | 0.0% |
| Cortisol | 1000 ng/mL | 0.0% |
| Cortisone | 1000 ng/mL | 0.0% |
| Danazol | 1000 ng/mL | 0.0% |
| Dexamethasone | 2000 ng/mL | 0.0% |
| DHEA | 1000 ng/mL | 0.0% |
| DHEA-S | 10 000 ng/mL | 0.0% |
| 5α-dihydrotestosterone (DHT) | 500 ng/mL | 0.0% |
| Estradiol | 1000 ng/mL | 0.0% |
| Estriol | 100 ng/mL | 0.0% |
| Estrone | 1000 ng/mL | 0.0% |
| Ethisterone | 100 ng/mL | 0.0% |
| Norgestrel | 100 ng/mL | 0.0% |
| Prednisone | 1000 ng/mL | 0.0% |
| Pregnenolone | 5000 ng/mL | 0.0% |
| Progesterone | 1000 ng/mL | 0.0% |
| Testosterone propionate | 1000 ng/mL | 0.0% |

The following substances do not interfere with a bias of $\pm 15\%$ in the Free Testosterone ELISA when the concentrations are below the stated threshold presented in the following table.

| Potentially Interfering Reagent | Threshold Concentration |
|---------------------------------|-------------------------|
| Bilirubin, conjugated | 15 mg/dL |
| Bilirubin, unconjugated | 15 mg/dL |
| Haemoglobin | 200 mg/dL |
| Total Protein | 7 g/dL |
| Triglyceride | 500 mg/dL |

14.7. Serum-plasma study

The Free Testosterone ELISA matrix comparison study was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI EP9-A3 guidelines. A total of 21 samples (18 native, 3 spiked) to cover the range were evaluated. Linear regression analysis was performed on the comparative data:

| Sample type | Slope [95% CI] | Intercept (pg/mL) [95% CI] | Correlation coefficient (r) |
|-----------------|------------------------|-------------------------------|--------------------------------|
| SST | 0.97 [0.92 to 1.02] | -0.02 [-1.01 to 0.98] | 0.99 |
| Lithium Heparin | 0.99 [0.97 to 1.02] | -0.01 [-0.40 to 1.02] | 1.00 |
| Sodium Heparin | 0.93 [0.87 to 0.99] | -0.01 [-1.12 to 1.10] | 0.99 |
| EDTA | 1.03 [0.97 to 1.08] | 0.27 [0.74 to 1.28] | 0.99 |

15. LIMITATIONS OF THE PROCEDURE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays⁶. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

16. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive. Nevertheless, all materials should still be regarded and handled as potentially infectious.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- Use appropriate personal protective equipment while working with the reagents provided.
- Some reagents contain small amounts of Proclin 300[®] as preservative. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.
- The clinical significance of Free Testosterone determination can be invalidated if the patient was treated with cortisone or natural or synthetic steroids.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- The ELISA is only designed for qualified personnel following the standards of good laboratory practice (GLP).

16.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

For information about the packaging materials refer to PACKAGING MATERIALS.

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
(96 Determinations)


BIBLIOGRAPHY/BIBLIOGRAPHIE/BIBLIOGRAFIA/BIBLIOGRAFIA

1. Shea JL, Wong PY, Chen Y. Free testosterone: clinical utility and important analytical aspects of measurement. *Adv Clin Chem*. 2014;63:59-84.
2. Diver MJ. Analytical and physiological factors affecting the interpretation of serum testosterone concentration in men. *Ann Clin Biochem*. 2006 Jan;43(Pt 1):3-12.
3. Brambilla DJ, Matsumoto AM, Araujo AB and McKinlay JB. The Effect of Diurnal Variation on Clinical Measurement of Serum Testosterone and Other Sex Hormone Levels in Men. *J Clin Endocrinol Metab*. 2009 Mar; 94(3): 907–913.
4. Rajfer J. Decreased Testosterone in the Aging Male. *Rev Urol*. 2003;5(suppl 1):S1–S2.
5. Basic QC Practices On-line Course; <http://www.Westgard.com>.
6. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27–33

PACKAGING MATERIALS / MATERIELS D'EMBALLAGE / MATERIALI D'IMBALLAGGIO / MATERIALES DE EMBALAJE

| | | | |
|--|--|--|---|
|  PAP 21 |  PAP 21 |  PAP 22 | <div style="border: 1px solid black; padding: 2px; text-align: center;">MTP</div>  ALU / LDPE 90 |
|--|--|--|---|

SYMBOLS KEY / EXPLICATION DES SYMBOLES / LEGENDA / SIMBOLOS

| | |
|---|--|
|  | Manufactured by / Fabriqué par / Prodotto da / Fabricado por |
| IVD | In Vitro Diagnostic Medical Device / Dispositif médical de diagnostic <i>in vitro</i> / Diagnostico <i>in vitro</i> / Producto para diagnóstico In vitro |
| LOT | Lot Number / Numéro de lot / Lotto/Numero de lote |
|  | Expiration Date / Date de péremption / Scadenza / Fecha de caducidad |
|  | Storage Temperature / Température de conservation / Temperatura di conservazione / Temperatura de almacenamiento |
|  | Keep away from sunlight / Protéger de rayonnement solaire / Tenere lontano dalla luce del sole / Mantener alejado de la luz solar |
| CE | CE Mark / Marquage CE / Marchio CE / Marca CE |
| REF | Catalogue Number / Référence du catalogue / Numero di codice / Número de Catálogo |
|  | Consult Instructions for Use / Consulter la notice d'utilisation / Consultare le istruzioni per l'uso / Consulte las Instrucciones de Uso |
| MTP | Microtiterplate / Plaque de Microtitrage / Piastra di Microtitolazione / Placa de Microtitulación |
| CONJ | Conjugate / Conjugué / Coniugato / Conjugado |
| CONTROL A | Control A/ Contrôle A/ Controllo A / Contro A |
| CONTROL B | Control B / Contrôle B / Controllo B / Contro B |
| CAL | Calibrator resp. Standard / Calibrateur resp Étalon / Calibratore ossia Standard / Calibrador o Estándar |
| SOLN STOP | Stop solution / Solution stop / Soluzione bloccante / Solución de parada |
| SUB TMB | TMB Substrate solution / Substrat TMB / soluzione substrato TMB / solución substrato TMB |
| WASH BUF 10x | Washing solution 10x concentrated / Solution de lavage concentré 10 x / soluzione di lavaggio concentrazione x10 / solución de lavado concentrado x10 |
|  | Contains sufficient for "n" tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para "n" tests |

**SUMMARY OF TEST PROCEDURE / RÉSUMÉ DE LA PROCEDURE DE TEST /
SCHEMA DELLA PROCEDURA / RESUMEN DE LA TÉCNICA**

SCHEME OF THE ASSAY

Free Testosterone

Test Preparation

Prepare reagents and samples as described.
Establish the distribution and identification plan for all specimens and controls.
Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

| | Substrate blank | Standard 0 - 5 | Controls | Sample |
|--|-----------------|----------------|----------|--------|
| Standard 0 - 5 | - | 20 µL | - | - |
| Controls | - | - | 20 µL | - |
| Sample | - | - | - | 20 µL |
| Conjugate | - | 100 µL | 100 µL | 100 µL |
| Cover wells with foil supplied in the kit Incubate for 1 h at 37 °C Wash each well three times with 300 µL diluted wash solution (if you use automated equipment wash the wells at least 5 times) | | | | |
| TMB Substrate | 100 µL | 100 µL | 100 µL | 100 µL |
| Incubate for exactly 15 min at room temperature in the dark | | | | |
| Stop Solution | 100 µL | 100 µL | 100 µL | 100 µL |
| Shake the microplate gently. Photometric measurement at 450 nm against a reference wavelength of 620-630 nm or against blank within 5 minutes. | | | | |

Gold Standard Diagnostics Frankfurt GmbH

Waldstrasse 23 A6

63128 Dietzenbach, Germany

Tel.: +49 6074 23698-0

Fax: +49 6074 23698-900

E-Mail: info.frankfurt@eu.goldstandarddiagnostics.com

Website: clinical.goldstandarddiagnostics.com

Date of issue: 2024-03-15

DNOV009_Free_TESTO_IFU_rev01_fromLot_6179AN